CE

# AKTIVLINE® ECO

Aluminium-Version  
Aluminum version

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### Revision status:

AKTIVLINE-ECO\_Gebrauchsanweisung\_interco\_  
de\_01

Please read these instructions for use carefully.

- ▶ Observe all safety instructions.
- ▶ If you have any questions, please contact us, interco Reha, or your specialist dealer.

Congratulations on your purchase of a quality product by interco Reha!

We have described how to handle the individual components in these instructions for use in order to make it easier for you to use the seat system. Please be aware that these guidelines also provide descriptions of components that may not be included with your product. It is important to adapt the product to meet the needs of each individual patient. The list of accessories and their corresponding descriptions are also a little more detailed in order to meet these many demands.

Please read the sections concerning your particular model carefully and keep these instructions for use in a safe place so that you can refer to them later. Please follow all the safety instructions and symbols to prevent any type of injury. If you should have any questions about your product, please ask us or your local qualified medical supplies dealer.

Please observe the following in particular:

- ▶ Parents, accompanying persons, nurses and family members who operate this product should be instructed in the use of the product by authorized special personnel using these instructions for use.
- ▶ This product was customized to the needs of the user. Subsequent changes to the product may only be made by trained specialist personnel. We recommend regular checks of the product and the product adaptation to ensure optimal long-term operation.

- ▶ This product may only be combined with the combination options listed later on. interco will not accept any liability for combinations with medical devices and/or accessories from other manufacturers.
- ▶ Technical modifications may only be made with the consent of the manufacturer.
- ▶ The manufacturer will not accept liability for modifications, changes to parts of the frame of the product and unapproved repairs to this product.
- ▶ The manufacturer reserves the right to technical changes to the model described in these instructions for use.

## 1.2. Product description

We are very pleased that you have chosen the AKTIVLINE® ECO seat system, a high-quality custom product made by interco Group.

The AKTIVLINE® ECO is the latest member of the AKTIVLINE® product family. With the dynamic arrangement of its seat shell on the undercarriage, the AKTIVLINE® ECO has the capacity to divert the uncontrolled muscle spasms suffered by spastic patients, for example, because the base of its seat system is able to open upwards. This results in a significant reduction in the stresses experienced by the patient. Other positive impacts on the patient were able to be observed in connection with this new concept of care. In addition, you will find that the AKTIVLINE® ECO model is considerably sturdier than conventional fixed/rigid care seat shells. The integrated dynamic range within the system also reduces wear on the materials. It is our sincere hope that your AKTIVLINE® ECO brings you great pleasure and makes your daily life easier!

### 1.3. Product overview

### 1.4. Delivery

The product is delivered pre-assembled in a crate or on a pallet, in some cases with detached rear wheels. If necessary, attachments such as head rests, arm rests and/or therapy tray tables are enclosed separately. For correct assembly and adjustment of the attachments please refer to Chapter 7 "Assembly/adjustment/operation". Please read all notes and instructions carefully before assembling.

As a rule, scope of delivery consists of:

- ▶ An undercarriage, roadworthy or for indoor use
- ▶ A seat system consisting of a seat and back unit, connected via physiological pivot points
- ▶ Additional accessories (customized to the needs of the user) depending on the order

- ▶ instructions for use with warranty card

#### Important!

Assembly and adjustments may only be carried out by qualified and trained specialist personnel.

### 1.5. Environmental conditions

The AKTIVLINE® ECO medical device must always be stored in a dry place to prevent possible corrosion. Please keep to ambient temperatures of 10 °C to +40 °C. Sand or other dirt particles should not be allowed to penetrate the product. This could damage the wheel bearings or mechanics. Use in salty water is not permissible.

## 2. Medical requirements and indications

### 2.1. Indications

Definition of indications according to the medical aids catalog of the statutory health insurers dated 26.11.2003: seat aids serve to compensate pronounced incorrect sitting postures and/or sitting posture instabilities. They are to enable permanent sitting in a physiological posture free from complaints. The seat system can be the basis for all further treatment/rehabilitation measures.

Seat shells allow patients suffering from considerably reduced or lacking stability of the torso or with pronounced deformities of the torso to sit correctly and without stress. In children with lacking stability of the torso, timely and appropriate provision of a seat shell can counteract the formation of poor posture and deformities.

Seat shells stabilize body posture in the therapeutically desired position through the full body design of these aids. This way they facilitate the patient's perception and contact to his/her surroundings, improve vital functions such as breathing and cardiovascular circulation, reduce spasms and allow active use of arms and hands for autonomous mobility and activities.

## 2.2. Scope of indications

*Extract from the National Association of Statutory Health Insurance Funds*

In case of medical necessity and to expand the area of application compared to the basic model and to adapt to specific clinical pictures and disabilities.

In any case, apart from the indications and contraindications listed here, the individual symptoms must be taken into account, accompanied by a regularly trained medical device advisor.

- ▶ ICD-10-G80 (Infantile Cerebral Palsy)
- ▶ ICD-10-G81.1 (spastic hemiparesis and hemiplegia)
- ▶ ICD-10-G82.1-G82.4 (spastic paraparesis and paraplegia / spastic tetraparesis and tetraplegia)
- ▶ ICD-10-P11 (Other Birth Injuries of the Central Nervous System)
- ▶ ICD-10-Q00-Q07 (congenital malformations of the nervous system)
- ▶ ICD-10-Q65-Q79 (congenital malformations and deformities of the musculoskeletal system)
- ▶ ICD-10-G10 (Huntington's Disease)
- ▶ ICD-10-G95.83 (spinal spasticity of the striated muscles)

- ▶ ICD-10-G24 (dystonia)
- ▶ ICD-10-R25.2, R25.8 (cramps and spasms of the muscles / other and unspecified abnormal involuntary movements)
- ▶ ICD-10-R29.2, R29.8 (Abnormal reflexes / other and unspecified symptoms affecting the nervous system and the musculoskeletal system)
- ▶ ICD-10-G25.5 (other chorea)

### Contra-indications and exclusions:

- ▶ ICD-10-Q78 (glass bone)
- ▶ ICD-10-H81 (vestibular / perception disorder)
- ▶ ICD-10-R42 (dizziness and staggering)
- ▶ ICD-10-F98.4 (rocking movements)
- ▶ ICD-10-F41 (anxiety disorders)
- ▶ ICD-10-R63.4, R63.5 (abnormal weight fluctuation)
- ▶ ICD-10-M25.5 (joint pain)
- ▶ ICD-10-M25.65 (joint stiffness in the pelvis)

## 2.3. Intended use

The dynamic seating system AKTIVLINE® ECO is used for impaired sitting, for severe functional and/or structural damage to the skeletal muscles and for spastic posture. It enables the controlled change between different postures via physiological pivot points, promotes breathing, circulation, digestion and the mobility of the body joints. It has a preventive effect on the development of pressure ulcer and other injuries. The system accompanies the patient in a controlled manner during the stretching movement and brings the patient back to the starting position without losing position. According to §33 German SGB V paragraph 1, the AKTIVLINE seating system serves to compensate the disability as far as possible.

## 3. Rating plate

The rating plate is located on the rear connecting tube of the base frame. The rating plate provides the following information:

- A** Type designation
- B** Patent number and country of manufacture
- C** Manufacturer's data/address
- D** Serial number (please always quote in the event of complaints or other questions)
- E** Observe warnings/safety instructions
- F** Read instructions for use before use
- G** Information on material strength
- H** Maximum weight / maximum load

## 4. Safety instructions

You are responsible for the safety of the patient. Your patient's safety could be at risk if you do not follow these guidelines. Nevertheless, it is not possible to cover all potential circumstances and unforeseeable situations. Common sense, caution, and discretion are factors that this product cannot provide; they must be

exercised by the person operating the product. Such persons must understand all of the instructions provided. You should explain how to operate the product (and the associated accessories) to everyone who will be using it.

Please contact us or your local specialist dealer if the instructions are unclear and a more detailed explanation is required or if you have any other questions.

### Meaning of the safety instructions

#### Warning!

Warning against possible severe damage due to accidents or injuries.

#### Caution!

Warning against possible damage due to accidents or injuries.

#### Note!

Warning for possible technical damage.

#### Caution!

*Risk of injury due to operating errors or lacking knowledge of the instructions for use*

- ▶ Please read the instructions for use carefully first.
- ▶ Familiarize yourself with the handling and functions of the product and practice appropriate handling procedures with the product prior to initial use. The safety of the user could be at risk if these guidelines are not followed.
- ▶ The maximum load for the AKTIVLINE® ECO is 70 kg.
- ▶ Do not use the foot rest when getting into or out of the wheelchair.

#### Warning!

*Unattended parking*

- ▶ We would like to expressly point out that seat shell occupants must be secured properly at all times by means of suitable restraint systems (straps).
- ▶ Never leave the user unattended, not even when the user has been securely strapped in and the brakes have been set.
- ▶ In order to avoid the risk of injury to other persons, never leave the wheelchair with its occupant in the care of untrained individuals.

#### Warning!

- ▶ Incorrect handling of packaging material
- ▶ Risk of suffocation! Keep the packaging materials out of reach of children.

#### Caution!

*Incorrect center of gravity / shift in center of gravity*

- ▶ Test the effects of shifts in the center of gravity on the behavior of the wheelchair on inclines, declines and lateral slopes, for example, or while negotiating obstacles only with the reliable support of an assistant. It is definitely advisable for inexperienced wheelchair users to employ an anti-tipping device.

- ▶ Fundamentally, the seat system should only be placed on horizontal, flat surfaces. If it should become necessary to set it down on a slope, then please make sure that the seat has been placed in an upright position. There is a risk of tipping backwards if the occupant is in the supine position on an incline.
- ▶ Ensure that the seat should only be reclined if the anti-tipping system is deployed (see chapter "Reclining the seat via two gas-pressure springs").
- ▶ Make sure that the patient does not lean too far out of the seat shell when reaching for objects (behind, to the side or in front of the wheelchair) because there is a risk of tipping or rolling over due to the shift in the center of gravity.
- ▶ The anti-tipping system is a device intended to prevent the patient from inadvertently tipping backwards. Under no circumstances is the anti-tipping system to be used in place of rollers, for example to transport a person in the wheelchair with the rear wheels removed.
- ▶ Do not hang heavy shopping bags or similar on the push bar or the levers! Risk of tipping.

#### Warning!

- ▶ Negotiating and using stairs
- ▶ Stairs are only permitted to be negotiated with the assistance of accompanying persons. Other facilities for this purpose, such as ramps or elevators, are to be used wherever available. If such facilities are not available, then the assistants must carry the wheelchair over smaller obstacle.
- ▶ Carrying techniques are explained in detail in chapter 10, it is essential to follow these.
- ▶ Make sure that the accompanying persons grip the seat shell wheelchair on firmly mounted or welded components only.
- ▶ **Do not** lift with the push bar, foot rest, arm rests or similar parts, this can result in a definite risk to the user.
- ▶ Never use escalators with a user in the wheelchair.

#### Caution!

*Regular check of brakes and wheels*

- ▶ Make sure that the clamp levers and quick-releases on the height-adjustable push bars are firmly tightened.
- ▶ Use the brakes to secure the seat shell undercarriage on uneven surfaces or when transferring the patient (to a motor vehicle, for example).
- ▶ The braking action of the knee lever brake, as well as the general handling of the wheelchair, depend on the air pressure in the tires. The seat shell wheelchair is considerably easier and safer to maneuver with properly inflated rear wheels and equal air pressure in both tires on one axle. Therefore, check the wheels for proper air pressure before using the wheelchair. The correct air pressure value is printed on the tire and/or listed in the air pressure table in chapter 9. Knee lever brakes are only effective with suitable air pressure and proper adjustment (max. 10 mm clearance, subject to technical modifications without notice).
- ▶ Please make sure that the Bowden cables for operating the drum brakes and the trigger for

the gas-pressure springs are not protruding into the spokes.

- ▶ Make sure that the wheels have adequate tread depth.

### Warning!

*Self-initiated modification of the settings*

- ▶ Adhere to the settings made by the specialist dealer. Only the settings described in chapter 7 may be adjusted by yourself.
- ▶ If you experience difficulties with the settings, please contact your specialist dealer or us.
- ▶ Discuss all changes to the settings properly with the specialist personnel or therapist.
- ▶ It is not permissible for any parties other than the manufacturer to make additions or modifications to your AKTIVLINE® ECO seat system and any cases that should become known can be subject to legal action. Likewise, the warranty obligations of interco Group shall cease to apply in such cases and any liability claims directed at the manufacturer shall no longer be valid.

### Caution!

*Public road traffic*

- ▶ Please note that you must observe traffic regulations when traveling on public streets.
- ▶ In dark conditions wear bright clothing whenever possible in order to be more visible. We also recommend installing an active light source.

### Caution!

*Risk of injury*

- ▶ In order to prevent hand injuries do not reach between the rear wheel and the knee lever brake when operating the wheelchair.
- ▶ Be aware of your fingers during operation and all adjustment and assembly work (risk of injury).

### Caution!

*When using the product*

- ▶ Make sure that screws are re-tightened after all adjustments.
- ▶ Always make sure that the seat shell is securely connected to the undercarriage.
- ▶ Always consult your local specialist dealer or the manufacturer in case of questions associated with the AKTIVLINE® ECO.
- ▶ Carry out regular visual inspections of all connections on the AKTIVLINE® ECO. Consult your local specialist dealer in case of doubt.
- ▶ Keep in mind that the AKTIVLINE® ECO is a customized care system that has been designed and approved for use by a single person only. The name of this person is included on the rating plate of the AKTIVLINE® ECO. The AKTIVLINE® ECO is intended for use by this person only and is not permitted to be used by anyone else.

### Caution!

- ▶ General cautionary measures
- ▶ The AKTIVLINE® ECO is only approved for its intended use. Avoid actions, such as driving the wheelchair against an obstacle (stair, curb) without using the brake or down from elevated surfaces.
- ▶ Use caution when handling fire (e.g. a burning cigarette)! Seat and back coverings could catch fire.
- ▶ Make sure to keep children from playing with the wheelchair or any parts thereof.
- ▶ There is a risk of injury due to trapping and crushing fingers.
- ▶ Children must not play next to, underneath or in front of the wheelchair when it is in motion and/or the seat or back unit is in use.

### Important!

- ▶ Prevent interventions by children and/or other persons!

## 5. Caution symbols

**Note:** For adjusting the angle, open control lever at least 3 turns!

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*For adjusting the angle, open control lever at least 3 turns!*

**Note:** Assure appropriate tilting protection when adjusting seat angle!

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

**Note:** Do not carry!  
**Location:** push bar right

**Note:** Climbing stairs with patient is prohibited!  
**Location:** push bar left

**Note:** Please do not use the upper frame for securing wheelchair for transport!  
**Location:** Top frame, undercarriage left and right

**Note:** It is not permitted to use this wheelchair as a car seat!  
**Location:** bottom frame, undercarriage left and right

**Note:** All screw connections that were loosened by changing settings must be re-tightened!  
**Location:** Upper frame left and right

## 6. Equipment/Options

**6.1. Standard equipment with interco undercarriage**  
A custom configured AKTIVLINE® ECO consists of a seat and back unit which is combined using physiological pivot points, as well as an aluminum undercarriage.

The seat and back unit are fitted with good anatomical basic padding as well as adduction guides and lateral side supports. The various accessories are selected and added according to individual requirements.

AKTIVLINE® ECO dynamic seat-backrest-system made of aluminium as a standard version with integrated physiological pivots in hip and knee area including:

- ▶ pelvic retraction unit and system fixation at 90° hip angle
- ▶ anatomical molded seat with abduction pommel with lateral hip support, as well as Velcro side pads for seat width adjustment
- ▶ anatomical molded backrest with (height and width) adjustable side supports for upper body, with lumbar pad (under cover), with extra backpolder for seat depth adjustment (2 cm)
- ▶ fabric cover: DOTS BLACK  
frame color: PEARLDARKGREY RAL 9023
- ▶ 5 cm seat depth regulation with continuously adjustable seat plate (3 cm) and option of adjusting the backrest unit (2 cm) and adjustment of the pivot points
- ▶ ready for fitting armrests
- ▶ ready for fitting head support and fixation aids
- ▶ anti-tip and support roller
- ▶ seating system mount on slide rail with retraction unit and slide rail in the back

An extensive range of accessories is available and is selected with care according to the needs of the patient and accompanying person by the specialist personnel.

## 6.2. All equipment and options

### 6.2.1. Lower extremities

<b>Foot rest system</b>	One-piece foot rest
	Split foot rest
<b>Foot rest cushion</b>	
<b>Calf support</b>	One-piece
	Split
<b>Foot rest mechanism locking option</b>	
<b>Foot pedals</b>	
<b>Heel stops</b>	
<b>Heel stop cushion</b>	
<b>Foot fixation</b>	Foot gaiters
	Foot gaiters with forefoot restraint

### 6.2.2. Pelvis

<b>Hip belt</b>	2-point design with neoprene cushion and press buckle
	4-point design with neoprene cushion and press buckle
<b>PELVI-LOCK 3DS pelvic retraction belt</b>	
<b>Abduction pommel</b>	Adjustable, removable
<b>Side cushions for the seat, attachable with Velcro</b>	

### 6.2.3. Upper body

<b>Chest harness made of neoprene</b>	Standard version
	With zipper
	Vario design
	Removable design
<b>Strap guide plates</b>	
<b>Lateral side supports</b>	Foldable with swing-away unit
	Height and width-adjustable
<b>Back extension, attachable by slipping on</b>	

### 6.2.4. Upper extremities

<b>Arm rests</b>	Height and angle-adjustable, ALULINE
	Height-adjustable, made of steel
<b>Therapy tray table</b>	Standard version made of wood
	Semi-circular design made of Neopolen
	Acrylic glass design
	Tabletop pad, attachable with Velcro
	Cover for therapy tray table upholstery
	Swing-away unit Slip-on mechanism

### 6.2.5. Head

<b>Mounting for head rest bracket</b>	
<b>Head rest bracket</b>	Rigid design with one fastening element
	Rigid design with two fastening elements
	Laterally adjustable design with one fastening element
	ALULINE, with depth adjustment
<b>Head rest cushion</b>	Shell-shaped Standard shape
<b>Fabric head rest cover</b>	
<b>Back extension incl. cover</b>	

## 6.2.6. Undercarriages

### AMF force application point system

#### Anti-tip

Solid tires	Rear wheels
	Front wheels

#### Wheel protectors

#### Spoke protectors

### Knee lever brake

Supporting rollers Standard

## 6.2.7. Custom items

So

Seat depth regulation in 5 cm increments via infinitely adjustable seat plate

# 7. Assembly/adjustment/operation

## 7.1. Standard equipment

### Parking brakes/drum brakes

AKTIVLINE® ECO is equipped with two independently operating parking drum brakes. These can be activated on both rear wheels. The brake levers are equipped with a three-level ratchet mechanism, which provides for varying levels of braking action.

#### Activating the parking brake

To activate the parking brakes, proceed as follows:

- ▶ Pull both brake handles on the push bar (Fig. 1).  
Make sure that the brake handle does not return to the starting position (locking mechanism will click audibly into place); otherwise the parking brake has not been activated.
- ▶ Check the braking action by attempting to push the AKTIVLINE® ECO.

#### Important!

If the braking action is no longer sufficient, then have the brakes readjusted by your medical supplies dealer.

- ▶ You can also use the parking brake for better braking control on descents. In such cases, engage the parking brake, as needed, into the first or second ratchet position only.


#### Releasing the parking brake

To release the parking brake, proceed as follows: deactivate the locking mechanism by activating the red ratchet lever and allowing the brake lever to glide back into the starting position. You may need to discharge the ratchet lever by pulling slightly on the brake lever before the brake can be released.

## Height- and angle-adjustable push bar

### Settings and adjustment options for the push bar

The push bar (Fig. 2) has various adjustment options. This makes it possible to achieve optimal customization to the person, who is pushing.


The following adjustment options are possible:

- ▶ Adjusting the angle of the top piece of the push bar
- ▶ Adjusting the height of the push bar unit
- ▶ Adjusting the inclination angle of the push bar unit relative to the vertical

#### Adjusting the angle of the top piece of the push bar

To adjust the angle of the push bar, proceed as follows:

- ▶ Activate the parking brakes.
- ▶ Loosen the screw connections (Fig. 3) on the right and left ratchet joints on the push bar by turn-

ing the screw nuts counterclockwise until the top piece of the push bar can be moved.

- ▶ Move the push bar into the desired position.
- ▶ Hold the push bar in the desired position with one hand and secure it in place by tightening the corresponding screw nut with the other hand.
- ▶ Secure the screw connection on the other side as well by tightening the screw nut.
- ▶ When tightening the ratchet joints make sure that the teeth engage properly.
- ▶ If necessary, release the parking brakes.

#### *Adjusting the height of the push bar unit*

Adjustment of the height of the push bar is always necessary when adapting the undercarriage to suit the person pushing the chair or if the inclination angle of the seat shell has been changed.

To adjust the height of the push bar, proceed as follows:

- ▶ First, loosen the quick-release fasteners on the fork tubes by opening the two quick-release levers outwards.
- ▶ Then adjust the push bar to the desired height. To do this, pull or press the push bar to the desired height.
- ▶ Secure the push bar in position by tightening the quick-release connectors. To do so, tilt the lever on the quick-release fasteners back in.
- ▶ If the clamping pressure is not sufficient (push bar not secured in place), then turn the screw across from the quick-release fastener clamp lever (with lever in open position) counter-clockwise, until sufficient clamping pressure is achieved by closing the clamp lever. Repeat this process, as necessary.

#### *Adjusting the inclination angle of the back unit (including the push bar) relative to the vertical*

In addition to the adjustment options already described above, it is also possible to swivel the entire push bar back by up to 30°. This can become necessary, for example, if a seat shell is being used with a seat angle of more than 90°. This enables the gap between the back of the seat shell and the push bar to be as uniformly small as possible. This adjustment option is not intended for individual adjustment to the person pushing the chair. To adjust the inclination angle of the push bar, proceed as follows:

- ▶ Loosen the two clamp levers located in front of the pivot point for the push bar folding mechanism by turning them counterclockwise.
- ▶ Apply (axial) pressure to the clamp lever and press the coach bolts out of the guide on the ratcheting element until it is possible to freely adjust the inclination adjustment angle.
- ▶ Move the push bar into the desired position.

#### **Important!**


Apply pressure to the rounded head in order to press the square head of the screw into the ratcheting element at the required inclination angle. 

- ▶ Check that the carriage bolt is correctly seated in the ratcheting element.
- ▶ Secure the unit in this position by turning both clamp levers clockwise to tighten them. It is also possible to maintain the inclination angle of the

back unit, including the push bar unit, when the push bar is folded down (e.g. during transport).

#### **Trapezoid-shaped adapter receptacle and adapter**

The seat shell is mounted to the undercarriage by means of an adapter system. This system consists of a trapezoidal adapter, which is mounted to the underside of the shell, and a trapezoid-shaped adapter receptacle which is connected to the undercarriage.


The adapter is securely locked by a lever located on the rear of the adapter holder (Fig. 4). Always make sure that it is firmly locked.

#### **Tilt-in-space mechanism via two gas-pressure springs**

SIMPLY and SIMPLY LIGHT are equipped with a gas-pressure spring adjustment for the inclination/tilting of the seat shell.

To adjust the tilt of the seat, proceed as follows:

- ▶ Activate the parking brakes.
- ▶ Pull both gas-pressure spring release levers mounted to the right and left of the push bar frame (Fig. 5) simultaneously.

- ▶ Adjust the inclination angle of the seat shell by pulling the push bar up or pressing it down while pulling the release lever for both gas-pressure springs.
- ▶ Let go of the release lever to lock the seat shell into the adjusted inclination position.
- ▶ Deactivate the parking brake, as appropriate.

### Important!

Never pull the release lever without a firm grip on the push bar if a patient is seated in the wheelchair. Prior to inclination of the seat shell, ensure that the support roller always face in direction of the driving surface. In a further undercarriage variant, the tilt of the seat shell is adjusted with a telescopic unit. First the mounting of the unit is loosened and the telescopic tube moved either upwards or downwards, depending on how the tilt is to be adjusted. Then the unit is refastened.

## Dynamic sitting

### Locking of the AKTIVLINE® ECO mechanism in the pelvic region incl. pelvis adjustment unit

The entire motion cycle of the shell is controlled via the physiological pivotal points of the hip and knee. The seat and back unit are mounted on sliding bars and allow opening of the hip angle to correspond to the patient's motion sequences. The knee-angle and foot rest units permit simultaneous stretching tendencies for the patient. A readjusting spring system both in the seat and the leg system returns the patient automatically into the physiological starting position.

#### Locking of the pelvis mechanism

Activate the pelvis mechanism (mobility in the hip region) by pulling the snap locks on the side of the AKTIVLINE® ECO below the seat shell (Fig. 6) from the groove and rotating them through 90°. However, sometimes, for example during meals or transport, it is desirable to keep the patient in a passive condition.

To deactivate the pelvis mechanism, rotate the snap locks through 90° and return them to their starting position until they engage and the pelvis mechanism is firmly locked.

### Important!

In case of a fault in the mobility of the seat shell, always check the pelvis mechanism first and then the foot rest mechanism!

## 7.2. Options


### 7.2.1. Lower extremities

#### Foot rest system

The foot rest system of the AKTIVLINE® ECO must be customized to the lower leg length of the patient. During adjustment it must be ensured that the patient's feet are comfortably positioned on the foot rest. The foot rest system of the AKTIVLINE® ECO is available in both individual as well as single-piece design.

In order to adjust the foot rest system to the lower leg length of the patient, release the screw on the clamp block of the foot rest and shift the foot rest up or down, as necessary, on the centered tubing of the foot rest system; this also applies to the foot rest system with individual foot supports. Then tighten the screw after the adjustment is complete.

#### Calf support

Calf supports are already fitted ex-works.

#### Foot rest retraction unit

The foot rest retraction unit is located at the foot rest suspension and mounted under the seating unit. The stretching movement of the legs is restricted by a gas-pressure spring or by an expander, i.e. the gas-pressure spring or the expander restrains stretching of the legs through resistance. If you wish to release the foot rest retraction unit to allow free movement for the patient's legs, then proceed as follows:

- ▶ The mounting of the foot rest retraction unit is fitted to the foot rest tube.
- ▶ Release the quick pin by pressing the red button and keeping it depressed. During pressing, pull the pin from the opening.
- ▶ This now hangs loosely attached to a suspension.
- ▶ If you wish to refit the foot rest retraction unit, hold the gas-pressure spring into the eyelet, which is attached to the middle tube of the foot rest, with one hand.
- ▶ Then press the red button of the quick pin and guide this into the opening/mounting intended for this purpose.
- ▶ Pull slightly to check whether the quick pin has engaged and whether stretching of the legs is again restrained by the gas-pressure spring.

### Important!

Make sure that you do not stand or knee directly in front of the foot rest during assembly. Sudden stretching movements by the patient could result in your being kicked or injured unexpectedly.

#### Foot rest mechanism locking option

Activate the foot rest mechanism (mobility in the leg area) by turning the clamp lever on the outside of the

seat shell (Fig. 7) counterclockwise until the foot rest can be moved. Make sure not to turn the clamp lever more than necessary, otherwise there is a risk of the lever popping out of the thread. To deactivate the foot rest mechanism, turn the clamp lever clockwise until the foot rest is firmly locked again.

**Important!**

In case of a fault in the mobility of the seat shell, always check the pelvis mechanism first and then the foot rest mechanism!

*Foot pedals*

The foot pedals are bolted in place onto the foot rest, these are already fitted ex-works and need not be adjusted further. However, as slipping of the legs is unavoidable for foot pedals, fixation of the feet is generally recommended.

*Heel stops*

Heel stops are already permanently fitted ex-works and no longer require any adjustment. Heel stops prevent slipping of the feet backwards.

*Foot fixation*

Foot restraints (Fig. 8) support the dynamic range of the AKTIVLINE® ECO foot rest system. They should be gently placed around the ankles to achieve a snug fit without exerting any pressure. Tensioning of the straps is only effective if they bend the foot to 90°. This requires leg and foot guidance which enables a right angle.

Fig. 8

The strap runs across the ankle joint and both ends of the strap are screwed to the foot rest. The images illustrate how the foot restraints are to be applied.

**7.2.2. Pelvis**

*Hip belt*

The AKTIVLINE® ECO comes equipped with a pelvic restraint (Fig. 9), which is of decisive importance for the seat system. Only proper adjustment of the pelvic restraint enables the patient to move the AKTIVLINE® ECO correctly, thus achieving the desired effects.

Various versions of the pelvic restraint are customized to meet the needs of the patient. There is a choice of 2-point or 4-point hip belts. As a rule, both versions are bolted in place onto the right and left of the shell or attached to the bottom shell edge with a snap buckle and clamping plate for length adjustment.

Fig. 9

Using the snap buckle and clamping plate, the length can be adjusted individually to the size of the patient, to do this, proceed as follows:

- ▶ Open the snap buckle on one side, adjust the belt to the desired length and reclose the snap buckle. Now proceed in the same sequence on the opposite side.
- ▶ Or: open the clamping plate with an Allen key, adjust the belt to the desired length, and then securely screw the clamping plate back again. Then repeat this procedure on the opposite side.

The hip belt itself is fitted with one, or, if necessary, two tensioning straps. These enable pressure to be applied to the pelvis, as desired or tolerated. To do this, first close the hip belt and then pull the flap, or both flaps simultaneously, to adjust the desired pressure. Please observe the following instructions when applying the pelvic restraint:

When applying the pelvic restraint, the patient must be seated as far back as possible in the seat shell, i.e. the patient's pelvis must be in contact with the back of the seat shell. In order to ensure this, it may be helpful to recline the seat shell using the undercarriage tilt-in-space mechanism until the patient slides comfortably into the desired position.

**Important!**

The patient must always be secured and under observation at all times.

*Abduction pommel*

An integrated abduction pommel is located under the cover and cannot be adjusted. This is always required in the case of a foot rest with individual foot supports. An adjustable or removable abduction pommel can be selected for a single-piece foot rest, proceed as follows: To loosen the abduction pommel, turn the clamp lever - mounted at the front of the abduction pommel - counterclockwise, this loosens the abduction pommel and it can simply be pulled out. Then reinsert the abduction pommel in the intended receptacle and turn the clamp lever clockwise until tight. Then check whether the abduction pommel has its desired position again.

*Side cushions for the seat, attachable with Velcro*

Side cushions with Velcro attachment serve to adapt to growth. These can be removed by a strong pull of the cushion. Zippers allow removal of the covers for washing. Please follow the exact washing instructions provided for washing the relevant materials.

### 7.2.3. Upper body

#### Chest harness

The chest harness is used to restrain the patient's entire upper body. The chest harness is tailored individually to the needs of the patient. To avoid pressure areas, the fasteners are always sewn to the harness ends. This makes slipping of the fasteners to the side impossible.

The lower end of the harness runs across the two spines (the two upper, anterior iliac spines) to the side of the pelvis. The middle part of the chest harness is guided across the torso at the front and ends at the height of the sternum. The two upper harness ends pass along the side of the neck over the shoulders.

The chest harness is an optimal correction aid to achieve symmetrical alignment of the pelvis and the torso. At the same time it continues straightening the pectoral girdle and the neck. By placing the chest harness symmetrically on the spines (straps pull to straighten the pelvis), the pull upwards helps to straighten the shoulders and to symmetrically stretch the spine, which supports the erection the head.

#### Important!

Always adjust the chest harness, working from the hard back of the shell outwards and using the tighteners located there.

#### Restraining the patient in the seat shell

**Precondition:** place the patient in the seat shell and lay the chest harness over the back edge of the shell.

To restrain the patient, proceed as follows:

- ▶ Position the patient correctly in the seating shell with the pelvis on the inner back of the shell.
- ▶ Lay the chest harness facing forwards. The head is then positioned between the two shoulder straps.
- ▶ The lower end of the chest harness is equipped with two tighteners for fastening the straps.
- ▶ Open the flap on one of the black strap tighteners (if this version applies).
- ▶ Now push the belt webbing through the seat cushion and the back cushion into the slot on the corresponding side between the moving flap and the fixed part.
- ▶ Pull the belt webbing as far as it will go.
- ▶ Hold the belt webbing in this position.
- ▶ Then close the flap on the tightener until it lies flat against the fixed part.
- ▶ Or use the plug lock, if available, to close the strap.
- ▶ To open the lock, pull the flap on the tightener upwards and pull out the belt webbing. Or open the plug lock with two fingers.
- ▶ Repeat this procedure on the other side.
- ▶ Repeat these steps in reverse order to release restraint by the chest harness.

Some chest harnesses are fitted with tensioning straps, proceed here as follows.

- ▶ Close the clip buckles or plug locks on both sides.
- ▶ Now pull on the tensioning straps on the right and left to adjust the desired pressure.
- ▶ Always start first in the pelvic region and then in the region of the shoulders.

Through correct guidance of the slots in the seat, the strap is pulled outwards/downwards and thus stabilizes the pelvis in central position. The pull direction of the straps should be an angled extension of the harness upwards/to the back (erection of shoulders).

To adjust the length of the straps of the chest harness, proceed as follows:

- ▶ Adjust the lower ends of the chest harness first.
- ▶ Open the lower tighteners on the hard shell back or under the seat. To do so, simply pull hard on the free end of the belt webbing and the tightener will open automatically.
- ▶ Pull the belt webbing to the desired length on one side.
- ▶ Hold the belt webbing in this position.
- ▶ Then close the flap on the tightener until it lies flat against the fixed part of the tightener.
- ▶ Repeat this procedure to adjust the other side.
- ▶ Repeat these steps to adjust the upper ends of the chest harness.
- ▶ Adjust these via the upper tighteners on the hard back of the seat shell. The procedure is identical as for the lower ends of the chest harness.

#### Important!

Ensure correct strap arrangement. Straps arranged too closely around the neck lead to a risk of strangulation. If the straps are arranged too loosely around the shoulders, the straps could slip off the patient's shoulders; the patient is then no longer secured and could fall out of the chair, resulting in possible injury.

#### Important!

The patient must always be secured and under observation at all times.

#### Strap guide plates

Strap guide plates (Fig. 10) are used to optimally adjust the chest harness for each patient. These allow individual adjustment of the strap guide. The strap guide plates are already mounted ex-works to the upper, inner edges of the back.

#### Important!

If strap guide plates are present, always use these for correct strap guidance. To do this, pull the straps on the right and left through the provided slot in the center of the plates and then lock the strap in the clip buckle. Make sure that the straps are not twisted or knotted.

The correct height of the straps is to be adjusted as follows:

- ▶ Loosen the two screws on the strap guide plates with an Allen key.
- ▶ Push the plate upwards or downwards into the best position.
- ▶ Tighten all screws again with an Allen key.
- ▶ Check strap arrangement on the patient.
- ▶ Readjust, as necessary.

**Important!**

Ensure correct strap arrangement. Straps arranged too closely around the neck lead to a risk of strangulation. If the straps are arranged too loosely around the shoulders, the straps could slip off the patient's shoulders; the patient is then no longer secured and could fall out of the chair, resulting in possible injury.

*Reclining supports with Velcro fastening*

Reclining supports are shaped to the patient, starting from the top. These pass from the edge of the shell over the collar bone. The pressure of reclining supports passes from the upper ribcage parallel to the sternum further outwards, upwards and to the back. This pressure helps the patient when erecting the pectoral belt. Reclining supports are screwed firmly to the back of the shell and are rotated sideways using the lock, this allows for easy transfer of the user.

*Lateral side supports*

Lateral side supports support the torso symmetrically, or displaced the spine accordingly, for example, in case of an S-scoliosis. Various versions of lateral side supports allow individual adaptation to the patient's needs. Adjustable lateral side supports can be fitted directly on site. The mounting for the adjustable lateral side supports is on the hard edge of the rear back part.

Fig. 11

It is best to adjust the lateral side supports with the patient seated in the chair. In order to achieve an ideal level of support, adjust one side first and then the other. To do so, proceed as follows:

- ▶ Loosen the two screws on the rear plates of the side supports (Fig. 11).
- ▶ Now adjust the lateral side supports to the patient's measurements. The lateral side supports can be adjusted both for height as well as width.
- ▶ Tighten all screws again.
- ▶ Proceed in the same manner on the opposite side.
- ▶ Check the fit of the lateral side supports and make any necessary changes.

*Lumbar spine support with Velcro attachment*

Support of the lumbar spine leads to better erection of the torso. The lumbar spine support with Velcro attachment is located underneath the cover.

Lumbar support allows lordotic bending of the spine in the lumbar region. This type of support is particularly suitable for patients who tend towards a kyphotic back posture and require support in torso erection.

To adjust the lumbar spine support correctly, proceed as follows:

- ▶ Lumbar spine support under the cover: open and remove the cover of the back part.
- ▶ Detach the lumbar spine support from the Velcro fleece connection.
- ▶ Place the lumbar spinal support in the correct desired position and attach again with Velcro.
- ▶ Check the fit of the lumbar spine support and make any necessary changes.
- ▶ Reclose the cover of the back cushion.

*Back extension, attachable by slipping on*

If the back length no longer provides adequate due to the patient's growth, then slip-on back extension is used. This is slipped onto the rigid back and is located between the back part and the head rest, thus extending the back height by the length of the back extension. Proceed as follows:

- ▶ Take the head rest from its mounting (see section "Head rest").
- ▶ Then insert the back extension with the hook facing downwards behind the cover of the back part.
- ▶ Now insert the head rest back into the head rest bracket through the back extension into the actual back unit.
- ▶ Then adjust the head rest correctly and secure all clamp levers.

**7.2.4. Upper extremities**

*Arm rests*

The AKTIVLINE® ECO seat system can be fitted with different arm rests of the ALULINE range. Arm rests give the patient sufficient support and stability and support erection of the torso. Mount the arm rests by releasing the clamp levers of the arm rest holders (Fig. 12/1) on the right and left sides of the seating shell, insert the arm rests and tighten the locking levers again. To adjust the height, change the insertion depth on the oval tube universal holder.


Fig. 12

To adjust the angle of the arm rest, proceed as follows:

- ▶ Open the clamp on one side of the arm rest bracket by turning the clamp lever (Fig. 12/2) counterclockwise.
- ▶ The angle of the arm rest to the seat shell can then be adjusted.
- ▶ Hold the arm rest in the desired position with one hand and secure it in place by tightening the clamp lever clockwise with the other hand.
- ▶ This will maintain the arm rest in position.
- ▶ Repeat this procedure on the other side.

#### Therapy tray table

Various shapes and materials were taken into consideration in the selection of the therapy tray table. All variants used can be fastened differently, nonetheless the AKTIVLINE® ECO has already been fitted with a mounting ex-works and this only needs to be adjusted to meet the patient's needs.

#### Variant with slip-on mechanism

The holders on the therapy tray table, two tubes on the right and left of the table tray, as well as the receptacles for these tubes which are mounted under the arm rests, are already pre-fitted on delivery. To apply the therapy tray table, proceed as follows:

- ▶ Insert the therapy tray table and the two tubes fitted to the therapy tray table resp. into the intended receptacles under the arm rests.
- ▶ Make sure that the therapy tray table is not pushed too close to the patient. Leave the patient sufficient space in the stomach region to avoid sitting trapped during stretching.
- ▶ You can adjust the height of the table via the arm rests.
- ▶ Mount the table in the correct position by turning the wing screws on the mountings under the arm rests clockwise.
- ▶ Repeat these steps in reverse order to detach the therapy tray table. Loosen the wing screws by turning counterclockwise and pull the therapy tray table from the holder away from the patient.

#### Variant with swing-away unit

To attach the therapy tray table, proceed as follows:

- ▶ There is a receiving tube under the table. Insert this in the intended holder under the arm rest.
- ▶ Make sure that the therapy tray table is not pushed too close to the patient. Leave the patient sufficient space in the stomach region to avoid sitting trapped during stretching.
- ▶ You can adjust the height of the table via the arm rests.
- ▶ Mount the table in the correct position by turning the wing screw on the mounting under the arm rest clockwise.
- ▶ Repeat these steps in reverse order to detach the therapy tray table. Loosen the wing screw by turning counterclockwise and pull the therapy tray table from the holder away from the patient.
- ▶ Once you have adjusted the therapy tray table into the correct position, you no longer need to remove it from the holder for transfers and the like.
- ▶ In the variant with swing-away unit, simply fold the table tray to the side. Always pay attention to the patient. Should the radius for swiveling not be sufficient, loosen the wing screw and pull the table tray slight forwards to increase the radius.

#### Variant with dynamic co-motion, with Velcro fleece connection and strap:

In the variant with dynamic co-motion, a Velcro fleece connection is fitted to the underside of the therapy tray table and the arm rests.

- ▶ Simply place the therapy tray table in its correct position on the arm rests and press gently.
- ▶ Now close the strap which is fitted to the sides of the therapy tray table by pressing the strap ends of the therapy tray table onto the strap ends of the receptacle mounted to the back.
- ▶ Adjust the correct height of the therapy tray table via the arm rests.
- ▶ Repeat these steps in reverse order to detach the therapy tray table, release the straps and simply pull the table tray from the arm rests in your direction.


## 7.2.5. Head

#### Head rest systems

Once the patient has been restrained in the correct position in the AKTIVLINE® ECO with the hip belt, adjust the head rest. The holder for the head rest is built into the back of the shell on the AKTIVLINE® ECO and was already fitted before delivery.

The AKTIVLINE® ECO seat system can be fitted with different head rest systems of the ALULINE range (standard or shell-shaped). Both can be adjusted in height and depth. Mount the head rest by loosening the securing lever at the head rest holder (Fig. 13/1) located in the upper part of the back rest, insert the head rest bracket and retighten the lever. To adjust the height, change the insertion depth on the oval tube universal holder.

To adjust the depth, proceed as follows:

- ▶ Loosen the clamp lever (Fig. 13/2) at the upper oval tubing by turning the clamp lever counterclockwise.
- ▶ Shift the tube of the headrest bracket until the desired position has been reached.
- ▶ Secure it into position by turning the clamp lever clockwise until it is no longer possible to slide the head rest bracket tube.

#### Standard-shaped head rest

This head rest is fitted with a core aluminum plate. The head rest can be bent individually to the needs of the patient. In case of pronounced asymmetry, the usual U-shape may sometimes require stronger bending on one side to counteract the sideways inclination of the head. The head rest is already supplied curved

ex-works, if you wish to alter the setting, please contact your local specialist dealer.

Ensure that the side elements are not bent too much to avoid exerting pressure on the head. Leave the patient enough space of movement to allow him/her to perceive his/her surroundings adequately and ensure that the patient's ears are not trapped or squashed.

#### Shell-shaped head rest

The shell-shaped head rest virtually positions the patient's head automatically in the center, due to the shaped side and centering in the recess, the head automatically "slips" into this position.

#### Warning!

The head rest is always to be adjusted such, that the gap between the head rest and the back has a maximum breadth of two fingers. There is a risk of strangulation if the patient's head slips between the head rest and the back!

#### Caution!

The gap between the head rest (bottom edge) and the upper edge of the back rest must be less than 60 mm (for children) or 120 mm (for adults) or greater than 250 mm according to DIN EN 12182.

#### Warning!

The combination of the AKTIVLINE® ECO with third party head rests (head rest systems not manufactured by interco) is not permitted. There is a risk of injury!

### 7.2.6. Undercarriage

*AKTIVLINE® ECO seat system, separable from the undercarriage*

#### Attaching the seat shell (if option exists)

The seat shell is attached to the trapezoid-shaped adapter receptacle on the undercarriage using the trapezoidal adapter. To do this, perform the following work steps and check secure mounting on the trapezoid-shaped adapter receptacle after placing the seat shell. First release the seat shell fixation by pulling out the snap lock and turning it by 90 degrees. Position the back end of the trapezoidal adapter to the front side of the trapezoid-shaped adapter receptacle so that the lock is pressed downwards.

Now push the set shell to the back until this engages audibly. Relock the seat shell fixation by rotating the snap lock through 90 degrees until it engages. Finally, check the seat shell and the seat shell locking system to ensure they are secure.

Now fix the security pin which connects the back unit of the seat shell to the undercarriage.

#### Important!

Always check whether the seat shell is firmly fitted to the trapezoid-shaped

adapter receptacle. Check this by pulling the seat shell and simultaneously pressing on the undercarriage. The seat shell may not slip out.

#### Removing the seat shell

First, pull out the two snap locks and turn through 90 degrees. Remove the security pin which connects the back unit of the seat shell to the undercarriage. You will now be able to move the seat shell system forward and detach it from the adapter.

Before reinstalling the shell, remember to turn the locking pins 90° again in order to reset them in the 'closed' position. Be careful that the security pin connecting the back unit with the undercarriage is again moved into place and that it is correctly set and adjusted.

#### AMF force application point system

See chapter on "Transport"

#### Anti-tip

The anti-tip is already factory-mounted. It is not possible to perform any readjustments.

#### Wheel protectors

The wheel protectors cover the upper part of the wheels. They consist of a plastic mudguard which is screwed to an aluminum mounting plate. For fitting to the chassis, the free end of the mounting plate is inserted into the upper slot of the wheel mounting.

To mount the wheel protector, proceed as follows:

- ▶ Take the wheel protector and place the groove-shaped part of the plate into the upper groove of the axle block.
- ▶ Insert the wheel protector to the stop. You can feel a slight click when the plate engages in the guide.
- ▶ Repeat this procedure on the other side.
- ▶ For disassembly, pull the wheel protector off in upward direction.

#### Spoke protectors

The spoke protectors are already factory-mounted. It is not possible to perform any readjustments. If it is necessary to perform repairs to the spoke protectors, please contact your medical supplies dealer.

#### Knee lever brake (air brake)

Knee lever brakes (Fig. 14) are used in addition to parking brakes (drum brakes). These are factory-mounted and adjusted with your order. If the braking action is diminished by changes in distance (between wheel and the knee-lever brake), you can readjust this by shifting the sliding nuts on the C-rail.

### Anti-tipping wheels

The rear end of the frame is equipped with an anti-tipping wheel (Fig. 15).

#### Important!

In normal driving position, the anti-tipping wheels must always point in the direction of the driving surface!

It may be necessary to swivel the anti-tipping wheels upwards when negotiating high curbs and for transport (reduces the size of the wheelchair for transport).

To do so, proceed as follows:

- ▶ Activate the parking brakes.
- ▶ Position yourself behind the undercarriage.
- ▶ Pull on one of the anti-tipping brackets first and slide it partly out of the adapter bracket groove.
- ▶ Then move the anti-tipping bracket upwards until it engages in the upper mounting groove.
- ▶ Now allow the anti-tipping bracket to slide back into the groove of the adapter bracket unit. Repeat the same steps with the second anti-tipping wheel.
- ▶ Release the parking brake.
- ▶ Follow the same procedure for returning into the original position.

#### Note!

Always check to ensure that the two anti-tipping wheels are always facing in the same direction. This is the only way to ensure effective protection against tipping!

### 7.2.7. Custom items

*Seat depth regulation in 5 cm increments via infinitely adjustable seat plate*

Seat depth regulation should be adjusted exclusively by your local specialist dealer or by us!

### 7.3. Optional accessories

#### 7.3.1. Lower extremities

*Rubber pads*

If a foot rest retraction unit has not been chosen, then rubber pads are used. These are already fitted ex-works and cannot be adjusted. Two rubber pads are located below the seating unit, these prevent the foot rest from hitting against the material in the rest phase after stretching.

#### 7.3.2. Undercarriages

*Solid tires*

These tires are solid and already pre-fitted on delivery. Adjustments cannot be made.

### 7.3.3. Custom items

*Hook for tube-fed nutrition*

The hook for tube-fed nutrition is located at the rear cross-section of the push bar. This can be shifted easily to the right or left.

#### Important!

Use the hook exclusively for attaching tube-fed nutrition! Not intended for rucksacks, carrier bags, etc. – risk of tipping.

## 8. Combination options

- 8.1. Combination options with pushing aids
- 8.2. Combination options with stair climbing aids
- 8.3. Combination options with electric wheelchairs

The combination of AKTIVLINE® ECO seating care systems with other products not manufactured by interco is not permitted. interco Group has made arrangements with some other manufacturers in terms of combining various additional options, for example, pushing aids or stair-climbing aids.

## 9. Technical data

### Important!

The maximum load (workload) is 80 kg for the AKTIVLINE® ECO seat system.

	Size 1 SB 40	Size 2 SB 40	Size 3 SB 40	Size 4 SB 45
Frame width undercarriage	40 cm	40 cm	40 cm	45 cm
Total weight incl. headrest and armrests	ca. 38 kg	ca. 40 kg	ca. 42 kg	ca. 44 kg
Seat width	20 - 25 cm	25 - 30 cm	30 - 35 cm	35 - 40 cm
Seat depth	25 - 30 cm	30 - 35 cm	35 - 40 cm	40 - 45 cm
Back height	40 cm	45 cm	50 cm	55 cm
Total width incl. armrests	68 cm	68 cm	68 cm	73 cm
Totals height incl. headrest	113 cm	118 cm	123 cm	128 cm
Height from ground to seat unit	60 cm	60 cm	60 cm	60 cm
Total length incl. stabilizer rollers	75 cm	80 cm	85 cm	90 cm
Max. load	80 kg	80 kg	80 kg	80 kg

Tire type	Air pressure max. [bar]
Front wheels, Ø 200 x 50 mm	2.5
Rear wheels 12 inch	2.5
Rear wheels 16 inch	2.5
Rear wheels 22 inch	4.0
Rear wheels 24 inch	4.0

### Note!

Please note that in choosing the appropriate tire pressure you will need to reach a compromise between riding comfort for the patient (low air pressure) and pushing comfort for accompanying persons (high air pressure). It may be necessary to experiment here. The factory-assembled tires are delivered with an air pressure of 2.5 bar.

## 10. Transport/transport lock/climbing stairs

---

### 10.1. General information on transport

Only the tested and approved wheelchair and occupant restraint system may be used for transporting the AKTIVLINE® ECO seat system in motor vehicles for the transportation of the disabled. In addition, a 4-point restraint system is required for transport in motor vehicles.

The AKTIVLINE® ECO seat system thus complies with the requirements of DIN 75078. For information on how to correctly secure the force application point system in the respective vehicle, please see the separate manufacturer's instructions and the following chapter 10.2 "Transport with an AMF force application point system".

#### Warning!

Never secure the transport lock to the frame or the back of the AKTIVLINE® ECO seat system!

### 10.2. Transport with an AMF force application point system

The AKTIVLINE® ECO seat system was subjected to a crash test as standard version with maximum fittings in accordance with ISO 7176-19.

The wheelchair is marked on four points which indicate the mounting points for the tension belts for transportation with motor vehicles. Only use the provided mounting points for strapping with tension belts. When fitted correctly, the head rest of the AKTIVLINE® ECO seat system is very robust, but cannot replace head rests fitted to motor vehicles.

The AKTIVLINE® ECO seat system has been testing without an additional drive. If an additional electrical drive is to be added, then this accessory is to be disassembled if the AKTIVLINE® ECO seat system is to be transported in a motor vehicle.

### 10.3. Climbing stairs

Stairs or high curbs are only to be negotiated with the assistance of two accompanying persons. Assistants are only permitted to grip the front frame tube and the back brace of the frame, without lifting or carrying the AKTIVLINE® ECO seat system. The anti-tipping device must be retracted and secured. When ascending, the assistant behind the AKTIVLINE® ECO seat system pulls the wheelchair up the step. The assistant at the front stabilizes the position by pressing the AKTIVLINE® ECO seat system evenly into the steps on both sides. When descending, the assistant at the front brakes by pressing the AKTIVLINE® ECO seat system evenly into the steps. The assistant behind holds the back brace, and stabilizes and maintains the wheelchair in the correct position. To prevent tire damage, slowly roll the AKTIVLINE® ECO down the steps. The anti-tipping system must then be properly readjusted.

#### Caution!

Do not hold the AKTIVLINE® ECO seat system by the push bar when negotiating stairs.

Never lift the care system by the foot rest or the foot supports or the arm rests or other attached parts.

Do not lift the care system with the user sitting in it.

### 10.4. Transport on aircraft

The AKTIVLINE® ECO seat system can only be transported without the patient in the cargo hold. There are no concrete instructions for transportation on an aircraft. It is essential to contact the airline beforehand to determine individual transportation options.

## 11. Possible faults

---

In case of faults, please immediately contact your local specialist dealer or contact us at interco Reha.

Your local specialized dealers:

#### interco Group GmbH

Im Auel 50  
53783 Eitorf

Tel.: +49 2243 8807-0  
Fax: +49 2243 8807-29  
E-mail: [info@interco.gmbh](mailto:info@interco.gmbh)  
[www.interco.gmbh](http://www.interco.gmbh)

If you have questions, you are also welcome to visit our Internet homepage at [www.interco.gmbh](http://www.interco.gmbh).

## 12. Material/cleaning

---

The AKTIVLINE® ECO medical device must always be stored in a dry place to prevent possible corrosion. Please keep to ambient temperatures of -10 °C to +40 °C. Sand or other dirt particles should not be allowed to penetrate the product. This could damage the wheel bearings or mechanics. Use in salty water is not permissible. The cushions can be removed for cleaning and care if this option was indicated when ordering.

### 12.1. Materials used

#### 12.1.1. Powder-coated aluminum

All AKTIVLINE® ECO undercarriages are covered with a durable plastic coating. Some attachments and the undercarriage are made of aluminum and then finished with a plastic coating. Powder varnishes are used for coating. Powder varnishes are solvent-free coating materials, composed of different resin systems, pigments as well as special additives.

#### 12.1.2. Fabrics

##### *Synthetic leather*

interco uses easy-care designer synthetic leather for Dynaline restraints. With its outstanding quality characteristics, this synthetic leather is stylish and ideal for all upholstery applications and coverings. Especially excellent features include long-term durability, a comfortable seating experience, resilience to wear, good hold, and tear resistance. In addition, the synthetic leather is very elastic, resistant to urine, blood and disinfectants, free of azo compounds, CFC, formaldehyde, as well as free of PCTs, PCPs and PCBs.

##### *Spacer fabric*

The spacer fabrics used by interco have two outer textile layers, which are connected by a number of threads and spaced at a specific distance; due to their special structure, they are also called 3D textiles. The thermoregulating pile layer wicks moisture away and shields against cold. 3D fabrics yield to pressure like a spring. The material features prolonged durability and pronounced elasticity. Due to the special 3D structure, there is a continuous insulating layer of air between the outer layers. The spacer fabrics feature high contracting force, breathability, and compression elasticity. They compensate well for temperature fluctuations and the material is lightweight, tear-resistant and easy to care for.

##### *Neoprene*

The neoprene used by interco is a foam based on chloroprene rubber (CR). Unique material properties are achieved in the compound due to a specific structural mixture of various types of rubber as the base substance and numerous other components. These include, in particular, excellent insulation properties, high dynamic load capacity, resistance to weather and ozone, good protection against aging, resistance to wear and tear, pronounced elasticity, and impact resilience, good adhesion properties, resistance to sea water and many chemicals, resistance to fungi and bacteria, and low gas permeability. Neoprene is also environmentally friendly, very well tolerated by the skin, and impermeable to water.

##### *Trevira CS (fabrics)*

Trevira CS fabrics are easy to clean, dry quickly and keep their shape. Unlike other fabrics, they require considerably less water for washing, less detergent, and less energy. Drying times are also short thanks to low water absorption. All in all, these are valuable economic benefits. Trevira CS is incomparably superior to flame-retardant cotton fabrics and is easy to disinfect in hospital settings, for example. Flame-retardant Trevira CS fabrics meet all major international fire safety standards.

All materials used are flame-retardant.

### 12.2. Cleaning of materials

#### 12.2.1. Powder-coated aluminum

All AKTIVLINE® ECO undercarriages are covered with a durable plastic coating. All parts can be cleaned with soapy water and a soft sponge. Avoid using aggressive or abrasive cleaning agents, as these can damage the surface of the plastic coating.

#### 12.2.2. Fabrics

You can play a role in maintaining the appearance of the valuable covers of your seating system by cleaning, vacuuming or carefully brushing them regularly before stains set in, by wiping up spilled liquids right away with an absorbent cloth - the longer you wait to do so, the more difficult it will be to clean. Usually it is helpful to apply a small amount of mild liquid detergent for delicate fabrics to a cloth and dab it on the spill (never rub vigorously). Important: Never apply detergent directly to the cover; always apply it to a cleaning cloth. Wipe gently with clean water, if necessary.

Please note: bleach or whiteners should never be applied to the material. The material should be exposed to as little direct sunlight as possible.

All of the fabrics for our covers are washable at 30°C. Use a detergent for delicate fabrics or a gentle dry-cleaning process to clean the covers.

Never dry the covers in the dryer! Iron at a low temperature [single dot]. If you should wish to wash the seat shell cover frequently, then it is advisable to have an auxiliary cover to put on the seat shell while washing the other cover. Ask your local specialized dealers on this topic. Please clean accessories with a conventional mild household soap solution.

### 12.3. Disinfection of materials

#### 12.3.1. Powder-coated aluminum

Please use a water-based agent for disinfection, e.g. Terralin, Quartamon Med or Sagrotan Original concentrate. Please observe the manufacturer's instructions for use when doing so. It is best to clean plastic parts with a commercially available plastic cleaner (please observe the product information for the respective cleaning agent).

**Note!**

If it is necessary to perform wipe disinfection of the AKTIVLINE® ECO seat system, then use a suitable disinfection agent according to the DGH disinfection list for surfaces.

All cleaning and disinfection agents must be rinsed thoroughly and the cleaned item dried fully before use. If this is neglected, this can lead to chemical reactions making the material brittle or lead to skin irritations of the patient.

**12.3.2. Fabrics**

Branded disinfection agents may be used according to the instructions of the fabric manufacturer. Please use a water-based agent for disinfection, e.g. Terralin, Quartamon Med or Sagrotan Original concentrate. Please observe the manufacturer's instructions for use when doing so.

**Important!**

Synthetic leather covers/skai covers may not be cleaned with disinfecting agents, as the material gets porous over time and cracks.

**13. Maintenance /inspection by the user**

Prior to being handed over to the patient, the AKTIVLINE® ECO sea system is checked for functionality and freedom of faults by trained and qualified employees. To ensure safe use of the AKTIVLINE® ECO seat system over its entire service life and to ensure a perfect technical condition, it must be maintained and inspected regularly according to the following maintenance and service plans.

Incorrect settings by the user or laypersons bear the risk of injuries or crushing, which in a worst case scenario could lead to the patient's death!

During normal use, the AKTIVLINE® ECO seat system should be checked regularly by the user for damage and the perfect working condition of all components.

**Note!**

The adjustment of settings is only permitted for qualified and trained personnel!

**Important!**

Only original parts may be used for service and maintenance work.

**13.1. Maintenance plan seating unit**

Checking / components	Daily	Weekly	Monthly
Visual inspection of the seating system, the attached components and accessories for visible damage	×		
Visual inspection of all fastening elements and screws	×		
Fastening/retightening of all fastening elements and screws		×	
Oil and grease moving parts			×

**13.2. Maintenance plan undercarriage**

Checking / components	Daily	Weekly	Monthly
Visual inspection of the undercarriage, the attached components and accessories for visible damage	×		
Check connection to seat system	×		
Check working condition of the brakes	×		
Fastening/retightening of all fastening elements and screws		×	
Check the chassis	×		
Oil and grease moving parts			×

**14. Service plan**

Service must be conducted every six months by a qualified and authorized specialist dealer in accordance with the service plans. The service is subject to charge by our authorized specialist dealers, so that the assumption of costs by the health insurer should be clarified in advance. The AKTIVLINE ECO seat systems are CE marked products, which basically allow use during the

growth phase of the patient, however, they do require a regular check of settings with regard to changes of the growing patient.

This check and potentially necessary adjustment of settings are restricted to qualified and trained specialist personnel only and are part of the regular service.

## 14.1. Service plan complete AKTIVLINE® ECO system

### Note!

Service should be conducted at six-monthly intervals! The manufacturer cannot grant warranty if the AKTIVLINE® ECO seat system has not been checked regularly (at least once per year) as part of service.

Pos.	Test / Execution	OK	Adjustments necessary
1	Cleaning and disinfection	<input type="checkbox"/>	<input type="checkbox"/>
2	Checking of all settings of the AKTIVLINE® ECO seat system in accordance with the CURRENT size of the patient	<input type="checkbox"/>	<input type="checkbox"/>
3	Check of setting of the head rest: <ul style="list-style-type: none"> <li>▶ Does the head rest lie against the patient's occiput (standard version)?</li> <li>▶ Does the bottom edge of the shell-shaped head rest lie against the base of the head (shell shape)?</li> <li>▶ Do the side elements exert pressure on the patient's head?</li> <li>▶ Does the patient have sufficient freedom of movement?</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
4	Check the safety distances according to DIN EN 12182:2012 using the following table	<input type="checkbox"/>	<input type="checkbox"/>

### Safety distances according to DIN EN 12182:2012 for children

To avoid	Between moving parts		Between fixed parts	
Trapping fingers	< 4 mm	or > 25 mm	< 5 mm	or > 12 mm
Trapping feet	< 25 mm	or > 120 mm	< 25 mm	or > 45 mm
Trapping the head	< 60 mm	or > 300 mm	< 60 mm	or > 250 mm
Trapping the genitals	< 8 mm	or > 75 mm	< 8 mm	or > 75 mm

### Safety distances according to DIN EN 12182:2012 for adults

To avoid	Between moving parts		Between fixed parts	
Trapping fingers	< 8 mm	or > 25 mm	< 8 mm	or > 25 mm
Trapping feet	< 35 mm	or > 120 mm	< 35 mm	or > 100 mm
Trapping the head	< 120 mm	or > 300 mm	< 120 mm	or > 250 mm
Trapping the genitals	< 8 mm	or > 75 mm	< 8 mm	or > 75 mm

### 14.1.1. Service plan seating unit

Pos.	Area	Test			
		Setting/function	OK	Damage/distortion	Exchange/replacement
1	<b>Performance test</b>	Back recline	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Fixing unit	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Arm rest	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Head rest	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Abduction pommel	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Therapy tray table	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
2	<b>Screw connections</b>	Screws tightened	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Rust spots	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3	<b>Appearance (visual inspection)</b>	Paintwork	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Flock coating	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Cleanliness	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Fabric seams (position / appearance)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Rating plate	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### 14.1.2. Service plan undercarriage

Pos.	Area	Test	Setting/function			
			OK	Damage/ distortion	Exchange/ replacement	
1	Main frame	Main frame	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Anti-tipping system	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
2	Trapezoidal adapter		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
3	Front wheel	Wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Tire pressure	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Running performance of the wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Front wheel fork	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
4	Rear wheel	Mounting to frame	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Tire pressure	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Running performance of the wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
5	Brakes	Quick-release axles	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Drum brake	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Performance test, level 1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Performance test, level 2	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
6	Adjustability	Performance test, level 3	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Bowden cable	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Gas-pressure springs	Damages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
			Screw connections	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Seat adjustment	Performance test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Folding mechanism	Folding mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Back angle adjustment		Locking mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	Push bar	Performance test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
7		Screw connections	Height adjustment mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Locking mechanism		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	Size and length		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	Completeness		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
8	Options	Damages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Loctite secured	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Covering caps available	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Wheel protector	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
9	Appearance (visual inspection)	Spoke protector	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Arm rests	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Centered foot rest	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Paintwork	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	Cleanliness	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Rating plate	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

#### Important!

Every change / replacement of relevant components is to be documented carefully.

## 15. Reuse/durability

---

Your product can be reused under certain conditions. It is CE marked and must be reconditioned for reuse by the manufacturer to meet all essential safety and performance requirements. This includes that it is properly adapted to the needs of the new user. In particular, individual fittings can only be reinstated after professional adaptation to the new patient. Accompanying documents, as well as these operating instructions, are part of the product and must be handed over to the new user when it is to be reused again.

The durability and service life of the product are not limited and may vary. They depend on the conditions of use and the environment. In addition, daily use and the degree of disability can affect the durability and service life of the product. Regular inspection and the performance of all necessary service and maintenance work can positively influence the durability and service life of the product.

## 16. Disposal/environmental compatibility

---

The issue of environmental protection is a high priority at interco Group.

We only process environment-friendly materials and specifically avoid manufacturing processes which are harmful to the environment. All packaging materials can be recycled. The AKTIVLINE® ECO seat system can

be dismantled into separate material components to comply with optimal recycling and disposal. Please make inquiries with your local waste management company about the regulations concerning disposal and recycling. All components of the product are to be disposed of according to the respective country-specific environmental regulations.

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Eitorf, 2021-01-12

  
Michael Markwald  
Geschäftsführer  
Chief Executive Officer

## 18. Warranty

In addition to our General Terms and Conditions, we grant a two-year warranty for the AKTIVLINE® ECO seat shell system product supplied by us. Unusable or defective parts will be repaired or replaced free of charge during the warranty period. The warranty shall exclude design-related modifications.

We accept no liability in the case of damage incurred as a result of combining our product with third-party products of any kind, which may pose significant risk in certain circumstances. Unless said products have been expressly approved by our company.

The warranty also excludes any defects occurring as the result of normal wear, excessive stress, violent and willful damage, and inappropriate use. Warranty obligations shall cease to apply if original interco replacement parts are not used for repairs and if repairs and modifications are not performed by interco. It is not permissible for any parties other than the manufacturer to make additions or modifications to

your AKTIVLINE® ECO seat system and any cases that should become known can be subject to legal action. Likewise, the warranty obligations of interco Group shall cease to apply in such cases and any liability claims directed at the manufacturer shall no longer be valid.


The continuous use of your AKTIVLINE® ECO seat system may result in loose screws, nuts, and spokes. Retighten these at periodic intervals or arrange to have this done by your local dealer or by us.

### Note!

The manufacturer cannot grant warranty if the AKTIVLINE® ECO seat system has not been checked regularly (at least once per year) as part of service.

The manufacturer cannot grant warranty if third party components (not manufactured by interco) are used without authorization by the manufacturer.

### Extended warranty

Upon submission in this form we will extend the warranty on your AKTIVLINE® ECO to cover a total of 2 years. Please fill out this warranty form carefully and completely and return it to us. Thank you very much.

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Status as per February 2021 | Subject to technical alterations without prior notification

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Management system certified as per

**DIN EN ISO 13485:2016**

Reg.-Nr. 73 105 6351

TÜV Hessen

**DIN EN ISO 9001:2015**

Reg.-Nr. 73 100 6351

TÜV Hessen

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© interco Group GmbH, Eitorf, 2021