# AKTIVLINE® Individual

## Aluminium-Version

Aluminum version | Aluminiumsversjon | Aluminium-versie

Individuelle Ausführung/Sonderanfertigung gemäß Verordnung (EU) 2017/745 über Medizinprodukte  
Individual design/custom-made device according to Regulation (EU) 2017/745 on medical devices  
Individuell utførelse/spesialmodell i henhold til forskriften (EU) 2017/745 om medisinsk utstyr  
Individuele aanpassing/speciaal model volgens regelgeving (EU) 2017/745 over medische apparaten

Gebrauchsanweisung  
Instructions for use  
Bruksanvisning  
Gebruikshandleiding

**interco**®  
Aktivline!

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*The medical device listed below is a custom-made device that has been manufactured in accordance with the aforementioned written prescription and corresponds to the designs and features specified by the prescriber. This medical device is intended exclusively for the person named below in order to meet their individual condition and individual needs.*

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*The dynamic seating system AKTIVLINE® is used for impaired sitting, for severe functional and/or structural damage to the skeletal muscles and for spastic posture. It enables the controlled change between different postures via physiological pivot points, promotes breathing, circulation, digestion and the mobility of the body joints. It has a preventive effect on the development of pressure ulcer and other injuries. The system accompanies the patient in a controlled manner during the stretching movement and brings the patient back to the starting position without losing position. According to §33 German SGB V paragraph 1, the AKTIVLINE seating system serves to compensate the disability as far as possible.*

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*The above-mentioned person MUSTERMANN, Max (m / f / d) / patient no. must sit in an ergonomically optimal seating position in order to participate in social life and to be mobile. For this purpose, the above-mentioned custom-made device was prescribed and manufactured in the specific version in order to compensate for the deficits and to promote existing residual activities.*

Eitorf, JJJJ-MM-TT

Michael Markwald  
CEO Chief Executive Officer

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**We declare under the sole responsibility of the manufacturer that the above mentioned product meets all applicable safety and performance requirements according Annex I & XIII of the EU regulation 2017/745 of the European Parliament and of the Council of 5 April 2017 on Medical devices and will be classified as medical device class 1 (annex VIII, rule 1) for transient and short time use.**

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# Contents

---

1.	<b>Introduction</b> .....	41	12.2.	Cleaning of materials.....	65
1.1.	Foreword.....	41	12.2.1.	Powder-coated aluminum.....	65
1.2.	Product description.....	41	12.2.2.	Fabrics.....	65
1.3.	Product overview.....	42	12.3.	Disinfection of materials.....	65
1.4.	Delivery.....	42	12.3.1.	Powder-coated aluminum.....	65
1.5.	Environmental conditions.....	42	12.3.2.	Fabrics.....	66
2.	<b>Medical requirements and indications</b> .....	43	13.	<b>Maintenance/inspection by the user</b> .....	66
2.1.	Indications.....	43	13.1.	Maintenance plan seating unit.....	66
2.2.	Scope of indications.....	43	13.2.	Maintenance plan undercarriage.....	66
2.3.	Intended use.....	43	14.	<b>Service plan</b> .....	67
3.	<b>Rating plate</b> .....	44	14.1.	Service plan complete AKTIVLINE® system.....	67
4.	<b>Safety instructions</b> .....	44	14.1.1.	Service plan seating unit.....	68
5.	<b>Caution symbols</b> .....	46	14.1.2.	Service plan undercarriage.....	68
6.	<b>Equipment/Options</b> .....	47	15.	<b>Reuse/durability</b> .....	69
6.1.	Standard equipment with various interco undercarriages.....	47	16.	<b>Disposal/environmental compatibility</b> .....	69
6.2.	All equipment and options.....	48	17.	<b>EU declaration for custom made products</b> .....	70
6.2.1.	Lower extremities.....	48	18.	<b>Warranty</b> .....	71
6.2.2.	Pelvis.....	48			
6.2.3.	Upper body.....	48			
6.2.4.	Upper extremities.....	48			
6.2.5.	Head.....	49			
6.2.6.	Undercarriages.....	49			
6.2.7.	Custom items.....	49			
6.3.	Optional accessories.....	49			
6.3.1.	Lower extremities.....	49			
6.3.2.	Pelvis.....	49			
6.3.3.	Upper body.....	49			
6.3.4.	Upper extremities.....	49			
6.3.5.	Head.....	49			
6.3.6.	Undercarriages.....	49			
6.3.7.	Custom items.....	49			
7.	<b>Assembly/adjustment/operation</b> .....	49			
7.1.	Standard equipment.....	49			
7.2.	Options.....	52			
7.2.1.	Lower extremities.....	52			
7.2.2.	Pelvis.....	53			
7.2.3.	Upper body.....	54			
7.2.4.	Upper extremities.....	57			
7.2.5.	Head.....	58			
7.2.6.	Undercarriage.....	59			
7.2.7.	Custom items.....	61			
7.3.	Optional accessories.....	62			
7.3.1.	Lower extremities.....	62			
7.3.2.	Pelvis.....	62			
7.3.3.	Upper body.....	62			
7.3.4.	Upper extremities.....	62			
7.3.5.	Head.....	62			
7.3.6.	Undercarriages.....	62			
7.3.7.	Custom items.....	63			
8.	<b>Combination options</b> .....	63			
9.	<b>Technical data</b> .....	63			
10.	<b>Transport/transport lock/ climbing stairs</b> .....	63			
10.1.	General information on transport.....	63			
10.2.	Transport with an AMF force application point system.....	64			
10.3.	Climbing stairs.....	64			
10.4.	Transport on aircraft.....	64			
11.	<b>Possible faults</b> .....	64			
12.	<b>Material/cleaning</b> .....	65			
12.1.	Materials used.....	65			
12.1.1.	Powder-coated aluminum.....	65			
12.1.2.	Fabrics.....	65			

# 1. Introduction

## 1.1. Foreword

Information  
Date of last update:  
November 2022

Revision status:  
AKTIVLINE-INDIVIDUAL\_Gebrauchsanweisung\_  
interco\_multiling\_06

Please read these instructions for use carefully.

- ▶ Observe all safety instructions.
- ▶ If you have any questions, please contact us, interco Group, or your specialist dealer.

Congratulations on your purchase of a quality product by interco Group!

We have described how to handle the individual components in these instructions for use in order to make it easier for you to use the seat system. Please be aware that these guidelines also provide descriptions of components that may not be included with your product. It is important to adapt the product to meet the needs of each individual patient. The list of accessories and their corresponding descriptions are also a little more detailed in order to meet these many demands.

Please read the sections concerning your particular model carefully and keep these instructions for use in a safe place so that you can refer to them later. Please follow all the safety instructions and symbols to prevent any type of injury. If you should have any questions about your product, please ask us or your local qualified medical supplies dealer.

Please observe the following in particular:

- ▶ Parents, accompanying persons, nurses and family members who operate this product should be instructed in the use of the product by authorized special personnel using these instructions for use.
- ▶ This product was customized to the needs of the user. Subsequent changes to the product may only be made by trained specialist personnel. We recommend regular checks of the product and the product adaptation to ensure optimal long-term operation.
- ▶ This product may only be combined with the combination options listed later on. interco will not accept any liability for combinations with medical devices and/or accessories from other manufacturers.

- ▶ Technical modifications may only be made with the consent of the manufacturer.
- ▶ The manufacturer will not accept liability for modifications, changes to parts of the frame of the product and unapproved repairs to this product.
- ▶ The manufacturer reserves the right to technical changes to the model described in these instructions for use.

## 1.2. Product description

We are very pleased that you have chosen the AKTIVLINE® seat system, a high-quality custom product made by interco Group.


The AKTIVLINE® is based on a completely new concept of care for the severely disabled. With the dynamic arrangement of its seat shell on the undercarriage, the AKTIVLINE® has the capacity to divert the uncontrolled muscle spasms suffered by spastic patients, for example, because the base of its seat system is able to open upwards. This results in a significant reduction in the stresses experienced by the patient. Other positive impacts on the patient were able to be observed in connection with this new concept of care. In addition, you will find that the AKTIVLINE® model is considerably sturdier than conventional fixed/rigid care seat shells. The integrated dynamic range within the system also reduces wear on the materials. If you have an understanding of the statutory regulations in the field of healthcare, you will notice that the AKTIVLINE® does not bear the CE mark. This is due to the fact that each and every AKTIVLINE® is specially built for individual use by one single patient. In other words, each AKTIVLINE® is a "one-off piece" that optimally adjusts to the needs of its owner. It is our sincere hope that your AKTIVLINE® brings you great pleasure and makes your daily life easier!

### Note!

Please note that this medical device is approved exclusively for the patient named on the rating plate.

### 1.3. Product overview

- 1 Push bar
- 2 Complete back
- 3 Top frame
- 4 Undercarriage
- 5 Back unit
- 6 Seating unit
- 7 Foot rest system

### 1.4. Delivery

The product is delivered pre-assembled in a crate or on a pallet, in some cases with detached rear wheels. If necessary, attachments such as head rests, arm rests and/or therapy tray tables are enclosed separately. For correct assembly and adjustment of the attachments please refer to Chapter 7 "Assembly/adjustment/operation". Please read all notes and instructions carefully before assembling.

As a rule, scope of delivery consists of:

- ▶ An undercarriage, roadworthy or for indoor use
- ▶ A seat system consisting of a seat and back unit, connected via physiological pivot points
- ▶ Additional accessories (customized to the needs of the user) depending on the order
- ▶ instructions for use with warranty card

### Important!

Assembly and adjustments may only be carried out by qualified and trained specialist personnel.

### 1.5. Environmental conditions

The AKTIVLINE® medical device must always be stored in a dry place to prevent possible corrosion. Please keep to ambient temperatures of 10 °C to +40 °C. Sand or other dirt particles should not be allowed to penetrate the product. This could damage the wheel bearings or mechanics. Use in salty water is not permissible.

## 2. Medical requirements and indications

### 2.1. Indications

Definition of indications according to the medical aids catalog of the statutory health insurers dated 26.11.2003: seat aids serve to compensate pronounced incorrect sitting postures and/or sitting posture instabilities. They are to enable permanent sitting in a physiological posture free from complaints. The seat system can be the basis for all further treatment/rehabilitation measures.

Seat shells allow patients suffering from considerably reduced or lacking stability of the torso or with pronounced deformities of the torso to sit correctly and without stress. In children with lacking stability of the torso, timely and appropriate provision of a seat shell can counteract the formation of poor posture and deformities.

Seat shells stabilize body posture in the therapeutically desired position through the full body design of these aids. This way they facilitate the patient's perception and contact to his/her surroundings, improve vital functions such as breathing and cardiovascular circulation, reduce spasms and allow active use of arms and hands for autonomous mobility and activities.

### 2.2. Scope of indications

In any case, apart from the indications and contraindications listed here, the individual symptom picture must be taken into account, accompanied by a regularly trained medical device consultant.

- ▶ ICD-10-G80 (Infantile cerebral palsy)
- ▶ ICD-10-G81.1 (Spastic hemiparesis and hemiplegia)
- ▶ ICD-10-G82.1-G82.4 (Spastic paraparesis and paraplegia/ spastic tetraparesis and tetraplegia)
- ▶ ICD-10-P11 (Other birth injuries of the central nervous system)
- ▶ ICD-10-Q00-Q07 (Congenital malformations of the nervous system)
- ▶ ICD-10-Q65-Q79 (Congenital malformations and deformities of the musculoskeletal system)
- ▶ ICD-10-G10 (Huntington's disease)
- ▶ ICD-10-G95.83 (Spinal spasticity of striated muscles)
- ▶ ICD-10-G24 (Dystonia)
- ▶ ICD-10-R25.2, R25.8 (Spasms and convulsions of muscles/Other and unspecified abnormal involuntary movements)
- ▶ ICD-10-R29.2, R29.8 (Abnormal reflexes/Other and unspecified symptoms involving the nervous and musculoskeletal systems)
- ▶ ICD-10-G25.5 (Other chorea)

The provision of a conventional rigid seat shell is often not indicated due to poor acceptance or difficult to realize. Particularly the high forces exerted by conventional seat shells are often not tolerated by patients. Furthermore, a conventional seat shell increases the risk for:

hip luxation, painful hips, shortening of the muscles, muscle contraction, reduced joint mobility, secondary damage due to immobility.

The dynamic AKTIVLINE® seat shell system allows the automatic response to different body postures and thus promotes breathing, circulation, digestion, mobility of the body's joints and acts prophylactically against pressure areas.

### Contra-indications and exclusions

- ▶ ICD-10-Q78 (vitreous bone)
- ▶ ICD-10-H81 (vestibular/perceptual disorder)
- ▶ ICD-10-R42 (vertigo and staggering)
- ▶ ICD-10-F98.4 (rocking movements)
- ▶ ICD-10-F41 (anxiety disorders)
- ▶ ICD-10-R63.4, R63.5 (abnormal weight fluctuation)
- ▶ ICD-10-M25.5 (joint pain)
- ▶ ICD-10-M25.65 (joint stiffness in pelvis)

If the permissible patient weight is exceeded, please note that the AKTIVLINE may require a reinforced frame and/or reinforcement of other components.

### 2.3. Intended use

The delivered medical device is a custom-made device according to Regulation (EU) 2017/745 on medical devices, Article 2 No. 3 and Annex XIII. The delivered medical device is an individual supply, which was specially made for the patient named on the written prescription. The designs and features listed and required on the written prescription are complied with in order to meet exclusively the individual condition and needs of the named patient. The supplied accessories as well as the attached options may not be modified without the approval and prescription of an authorized person, who is entitled to issue prescriptions based on his professional qualifications according to the national legislation. Any modifications interfere with the specified purpose and are not permitted.


The medical device is marked with the patient's name and/or acronym and may only be used by this patient.

The AKTIVLINE dynamic supply system is used in cases of impaired sitting, severe functional and/or structural damage to the skeletal musculature and spastic malpositions. It enables controlled change between different postures via physiological pivot points, promotes breathing, circulation, digestion and mobility of body joints. It has a preventive effect on the development of pressure sores and other injuries. The system accompanies the patient in a controlled manner into the stretching movement and back to the starting position, without loss of position.

According to §33 SGB V paragraph 1, the AKTIVLINE fitting serves to compensate for the disability as far as possible.

### 3. Rating plate

The rating plate is located on the rear connecting tube of the base frame. The rating plate provides the following information:

- A Type designation
- B Manufacturer's data
- C Year of manufacture
- D Serial number
- E Observe warnings/safety instructions
- F Read instructions for use before use
- G REF no.
- H Patient's name
- I Patient's acronym
- J Maximum load

### 4. Safety instructions

You are responsible for the safety of the patient. Your patient's safety could be at risk if you do not follow these guidelines. Nevertheless, it is not possible to cover all potential circumstances and unforeseeable situations. Common sense, caution, and discretion are factors that this product cannot provide; they must be exercised by the person operating the product. Such persons must understand all of the instructions provided. You should explain how to operate the product (and the associated accessories) to everyone who will be using it.

Please contact us or your local specialist dealer if the instructions are unclear and a more detailed explanation is required or if you have any other questions.

#### Meaning of the safety instructions

##### Warning!

Warning against possible severe damage due to accidents or injuries.

##### Caution!

Warning against possible damage due to accidents or injuries.

##### Note!

Warning for possible technical damage.

##### Caution!

*Risk of injury due to operating errors or lacking knowledge of the instructions for use*

- ▶ Please read the instructions for use carefully first.
- ▶ Familiarize yourself with the handling and functions of the product and practice appropriate handling procedures with the product prior to initial use. The safety of the user could be at risk if these guidelines are not followed.

- ▶ The maximum load for the AKTIVLINE® is 70 kg.
- ▶ Do not use the foot rest when getting into or out of the wheelchair.

##### Warning!

###### *Unattended parking*

- ▶ We would like to expressly point out that seat shell occupants must be secured properly at all times by means of suitable restraint systems (straps).
- ▶ Never leave the user unattended, not even when the user has been securely strapped in and the brakes have been set.
- ▶ In order to avoid the risk of injury to other persons, never leave the wheelchair with its occupant in the care of untrained individuals.

##### Warning!

###### *Incorrect handling of packaging material*

- ▶ Risk of suffocation! Keep the packaging materials out of reach of children.

##### Caution!

###### *Incorrect center of gravity / shift in center of gravity*

- ▶ Test the effects of shifts in the center of gravity on the behavior of the wheelchair on inclines, declines and lateral slopes, for example, or while negotiating obstacles only with the reliable support of an assistant. It is definitely advisable for inexperienced wheelchair users to employ an anti-tipping device.
- ▶ Fundamentally, the seat system should only be placed on horizontal, flat surfaces. If it should become necessary to set it down on a slope, then please make sure that the seat has been placed in an upright position. There is a risk of tipping backwards if the occupant is in the supine position on an incline.

- ▶ Ensure that the seat should only be reclined if the anti-tipping system is deployed (see chapter "Reclining the seat via two gas-pressure springs").
- ▶ Make sure that the patient does not lean too far out of the seat shell when reaching for objects (behind, to the side or in front of the wheelchair) because there is a risk of tipping or rolling over due to the shift in the center of gravity.
- ▶ The anti-tipping system is a device intended to prevent the patient from inadvertently tipping backwards. Under no circumstances is the anti-tipping system to be used in place of rollers, for example to transport a person in the wheelchair with the rear wheels removed.
- ▶ Do not hang heavy shopping bags or similar on the push bar or the levers! Risk of tipping.

### Warning!

#### *Negotiating and using stairs*

- ▶ Stairs are only permitted to be negotiated with the assistance of accompanying persons. Other facilities for this purpose, such as ramps or elevators, are to be used wherever available. If such facilities are not available, then the assistants must carry the wheelchair over smaller obstacle.
- ▶ Carrying techniques are explained in detail in chapter 10, it is essential to follow these.
- ▶ Make sure that the accompanying persons grip the seat shell wheelchair on firmly mounted or welded components only.
- ▶ **Do not** lift with the push bar, foot rest, arm rests or similar parts, this can result in a definite risk to the user.
- ▶ Never use escalators with a user in the wheelchair.

### Caution!

#### *Regular check of brakes and wheels*

- ▶ Make sure that the clamp levers and quick-releases on the height-adjustable push bars are firmly tightened.
- ▶ Use the brakes to secure the seat shell undercarriage on uneven surfaces or when transferring the patient (to a motor vehicle, for example).
- ▶ The braking action of the knee lever brake, as well as the general handling of the wheelchair, depend on the air pressure in the tires. The seat shell wheelchair is considerably easier and safer to maneuver with properly inflated rear wheels and equal air pressure in both tires on one axle. Therefore, check the wheels for proper air pressure before using the wheelchair. The correct air pressure value is printed on the tire and/or listed in the air pressure table in chapter 9. Knee lever brakes are only effective with suitable air pressure and proper adjustment (max. 10 mm clearance, subject to technical modifications without notice).
- ▶ Please make sure that the Bowden cables for operating the drum brakes and the trigger for the gas-pressure springs are not protruding into the spokes.
- ▶ Make sure that the wheels have adequate tread depth.

### Warning!

#### *Self-initiated modification of the settings*

- ▶ Adhere to the settings made by the specialist dealer. Only the settings described in chapter 7 may be adjusted by yourself.
- ▶ If you experience difficulties with the settings, please contact your specialist dealer or us.
- ▶ Discuss all changes to the settings properly with the specialist personnel or therapist.
- ▶ It is not permissible for any parties other than the manufacturer to make additions or modifications to your AKTIVLINE® seat system and any cases that should become known can be subject to legal action. Likewise, the warranty obligations of interco Group shall cease to apply in such cases and any liability claims directed at the manufacturer shall no longer be valid.

### Caution!

#### *Public road traffic*

- ▶ Please note that you must observe traffic regulations when traveling on public streets.
- ▶ In dark conditions wear bright clothing wherever possible in order to be more visible. We also recommend installing an active light source.

### Caution!

#### *Risk of injury*

- ▶ In order to prevent hand injuries do not reach between the rear wheel and the knee lever brake when operating the wheelchair.
- ▶ Be aware of your fingers during operation and all adjustment and assembly work (risk of injury).

### Caution!

#### *When using the product*

- ▶ Make sure that screws are re-tightened after all adjustments.
- ▶ Always make sure that the seat shell is securely connected to the undercarriage.
- ▶ Always consult your local specialist dealer or the manufacturer in case of questions associated with the AKTIVLINE®.
- ▶ Carry out regular visual inspections of all connections on the AKTIVLINE®. Consult your local specialist dealer in case of doubt.
- ▶ Keep in mind that the AKTIVLINE® is a customized care system that has been designed and approved for use by a single person only. The name of this person is included on the rating plate of the AKTIVLINE®. The AKTIVLINE® is intended for use by this person only and is not permitted to be used by anyone else.

### Caution!

#### *General cautionary measures*

- ▶ The AKTIVLINE® is only approved for its intended use. Avoid actions, such as driving the wheelchair against an obstacle (stair, curb) without using the brake or down from elevated surfaces.

- ▶ Use caution when handling fire (e.g. a burning cigarette)! Seat and back coverings could catch fire.
- ▶ Make sure to keep children from playing with the wheelchair or any parts thereof.
- ▶ There is a risk of injury due to trapping and crushing fingers.
- ▶ Children must not play next to, underneath or in front of the wheelchair when it is in motion and/or the seat or back unit is in use.

**Important!**


- ▶ Prevent interventions by children and/or other persons!

## 5. Caution symbols



**Note:** to adjust the angle, open the operating lever with at least 3 turns!

**Location:** push bar left

**Note:** when adjusting the angle of the seat shell ensure adequate protection against tipping!

**Location:** push bar left



**Note:** climbing stairs with patient is prohibited!

**Location:** push bar left


**Note:** do not use the top frame as transport lock!

**Location:** top frame, undercarriage left and right


**Note:** this wheelchair must not be used a seat in automobiles!

**Location:** bottom frame, undercarriage left and right


**Note:** this wheelchair is fitted with a force application point system and can be used a seat in motor vehicles if the belt is applied correctly.

**Location:** bottom frame, undercarriage left and right


**Note:** loosen clamp lever with at least two turns. Then re-tighten the clamp lever firmly!

**Location:** head rest bracket

**Note:** loosen clamp lever with at least two turns. Then re-tighten the clamp lever firmly!

**Location:** arm rest brackets

## 6. Equipment / Options

### 6.1. Standard equipment with various interco undercarriages

A custom configured AKTIVLINE® consists of a seat and back unit which is combined using physiological pivot points, as well as a SIMPLY® or SIMPLY LIGHT® aluminum undercarriage. The seat and back unit are fitted with good anatomical basic padding as well as adduction guides and lateral side supports. The various accessories are selected and added according to individual requirements. The AKTIVLINE® seat and back unit in the above basic design can also be fitted to an electric wheelchair provided by the customer (only approved models) as well as a ROOMY NE® indoor undercarriage including electrical motor.

AKTIVLINE® dynamic aluminum seat system on aluminum undercarriage, consisting of:

- ▶ 12"/16"/22"/24" wheels with drum brakes and quick-release axles (SIMPLY/SIMPLY LIGHT)
- ▶ pivoting supporting rollers
- ▶ parking brake handles
- ▶ height- and angle-adjustable push bar
- ▶ trapezoid-shaped adapter receptacle and adapter
- ▶ infinitely adjustable via sliding bars
- ▶ tilt-in-space mechanism via two gas-pressure springs
- ▶ system restraint to 90° sitting position
- ▶ lateral side supports
- ▶ adduction guides
- ▶ anatomical molding for seat and back
- ▶ Locking of the AKTIVLINE® mechanism in the pelvic region incl. pelvis adjustment unit

An extensive range of accessories is available and is selected with care according to the needs of the patient and accompanying person by the specialist personnel.

## 6.2. All equipment and options

### 6.2.1. Lower extremities

<b>Foot rest system</b>	One-piece foot rest Split foot rest
<b>Calf support</b>	One-piece Split
<b>Foot rest retraction unit</b>	
<b>Foot rest mechanism locking option</b>	
<b>Option of adjusting foot angle</b>	
<b>Foot pedals</b>	
<b>Heel stops</b>	
<b>Foot fixation</b>	Foot guide bar Foot gaiters Foot gaiters with forefoot restraint Foot straps with Velcro fleece connection Foot straps with deflection roller and padded instep Foot straps with clip buckle and padded instep Foot straps of leather Foot straps with padded instep and incremental adjustment
<b>Knee support</b>	With 2 individual press buckles With 2 press buckles, in two-piece construction Cover for knee support

### 6.2.2. Pelvis

<b>Hip belt</b>	2-point design with neoprene cushion and press buckle 4-point design with neoprene cushion and press buckle With deflection roller
<b>Sling seat</b>	Standard T-shaped Leg-spreading
<b>Abduction pommel</b>	Velcro attachment on the cover Integrated in knee support Adjustable, removable

### 6.2.3. Upper body

<b>Chest belt</b>	Can be attached with Velcro With flap With deflection roller
<b>Chest harness made of neoprene</b>	Standard version With zipper Vario design Removable design
<b>Strap guide plates</b>	
<b>Reclining supports with Velcro fastening</b>	
<b>Lateral side supports</b>	Foldable with swing-away unit Height and width-adjustable Embraces the trunk, including hinges and safety belt
<b>Lumbar spine support can be attached with Velcro</b>	With Velcro attachment on top of the cover With Velcro attachment under the cover
<b>Back extension, attachable by slipping on</b>	

### 6.2.4. Upper extremities

<b>Arm rests</b>	Height and angle-adjustable, ALULINE Height-adjustable, made of steel
<b>Upper arm supports</b>	
<b>L-shaped arm rests</b>	With integrated upper arm restraint Without integrated upper arm restraint
<b>Lower arm restraint</b>	With expander Without expander
<b>Hand cuff restraints</b>	
<b>Therapy tray table</b>	Standard version made of wood Semi-circular design made of wood Acrylic glass design Semi-circular design made of Neopolen Made of Neopolen, dynamically responsive to movement

Tabletop pad, attachable with Velcro  
 Cover for therapy tray table  
 Swing-away unit  
 Slip-on mechanism

### 6.2.5. Head

#### Mounting for head rest bracket

**Head rest bracket** Rigid design with one fastening element  
 Rigid design with two fastening elements  
 Laterally adjustable design with one fastening element  
 Laterally adjustable design with two fastening elements  
 ALULINE, with depth adjustment

**Head rest cushion** Shell-shaped  
 Standard shape  
 Dynamic head rest  
 Head rest, can be attached with Velcro  
 Cover for head rest

### 6.2.6. Undercarriages

- a. SIMPLY
- b. SIMPLY LIGHT
- c. ROOMY NEW EDITION
- d. Electric wheelchair

**AKTIVLINE® seat system, separable from the undercarriage**

**AKTIVLINE® undercarriage, foldable**

**AMF force application point system**

**Anti-tip**

**Solid tires** Rear wheels  
 Front wheels

**Wheel protectors**

**Spoke protectors**

**Quick-release axles for front wheels**

**Locking of the front wheels**

**Knee lever brake**

**Supporting rollers** Standard  
 Can be swiveled aside by foot pressure

**Stabilizer bar for AKTIVLINE® top frame, able to be fixed in position**

### 6.2.7. Custom items

**Seat depth regulation in 5 cm increments via infinitely adjustable seat plate**

**Dynamic back extension**

**Locking of the AKTIVLINE® mechanism in the pelvic region via hydraulic thrust absorber**

**Spare plates, AKTIVLINE®**

### 6.3. Optional accessories

#### 6.3.1. Lower extremities

- ▶ Rubber pads on foot rest

#### 6.3.2. Pelvis

- ▶ Pelvic bracket

#### 6.3.3. Upper body

- ▶ Custom-made chest harnesses

#### 6.3.4. Upper extremities

- ▶ Custom-made arm restraints

#### 6.3.5. Head

- ▶ Head rest bracket, reinforced version
- ▶ Head rest bracket with ball joint for head rest made of polyurethane foam with ball joint

#### 6.3.6. Undercarriages

- ▶ Split push bar
- ▶ Custom-made push bar
- ▶ Grip ring covers
- ▶ Distance block/ widened wheel base
- ▶ Brake lever extension

#### 6.3.7. Custom items

- ▶ Holder for respirator
- ▶ Maximum hip opening angle
- ▶ AKTIVLINE® rucksack
- ▶ Hook for tube-fed nutrition

## 7. Assembly/adjustment/operation

### 7.1. Standard equipment

#### Parking brakes/drum brakes

SIMPLY and SIMPLY LIGHT are equipped with two independently operating parking drum brakes. These can be activated on both rear wheels. The brake levers are equipped with a three-level ratchet mechanism, which provides for varying levels of braking action.

#### Activating the parking brake

To activate the parking brakes, proceed as follows:

- ▶ Pull both brake handles on the push bar. Make sure that the brake handle does not return to the starting position (locking mechanism will click audibly into place); otherwise the parking brake has not been activated.
- ▶ Check the braking action by attempting to push the SIMPLY or SIMPLY LIGHT.

**Important!**

If the braking action is no longer sufficient, then have the brakes readjusted by your medical supplies dealer.

- ▶ You can also use the parking brake for better braking control on descents. In such cases, engage the parking brake, as needed, into the first or second ratchet position only.


*Releasing the parking brake*

To release the parking brake, proceed as follows: deactivate the locking mechanism by activating the ratchet lever and allowing the brake lever to glide back into the starting position. You may need to discharge the ratchet lever by pulling slightly on the brake lever before the brake can be released.

**Height- and angle-adjustable push bar**

*Settings and adjustment options for the push bar*

The push bar has various adjustment options. This makes it possible to achieve optimal customization to the person being pushed.

The following adjustment options are possible:

- ▶ Adjusting the angle of the top piece of the push bar
- ▶ Adjusting the height of the push bar unit
- ▶ Adjusting the inclination angle of the push bar unit relative to the vertical

*Adjusting the angle of the top piece of the push bar*

To adjust the angle of the push bar, proceed as follows:

- ▶ Activate the parking brakes.
- ▶ Loosen the screw connections on the right and left ratchet joints on the push bar by turning the clamp levers counterclockwise until the top piece of the push bar can be moved.
- ▶ Move the push bar into the desired position.
- ▶ Hold the push bar in the desired position with one hand and secure it in place by tightening the corresponding clamp lever with the other hand.
- ▶ Secure the screw connection on the other side as well by tightening the clamp lever.
- ▶ When tightening the ratchet joints make sure that the teeth engage properly.
- ▶ If necessary, release the parking brakes.

*Adjusting the height of the push bar unit*

Adjustment of the height of the push bar is always necessary when adapting the undercarriage to suit the person pushing the chair or if the inclination angle of the seat shell has been changed.

To adjust the height of the push bar, proceed as follows:

- ▶ First, loosen the quick-release fasteners on the fork tubes by opening the two quick-release levers outwards.
- ▶ Then adjust the push bar to the desired height. To do this, pull or press the push bar to the desired height.
- ▶ Secure the push bar in position by tightening the quick-release connectors. To do so, tilt the lever on the quick-release fasteners back in.
- ▶ If the clamping pressure is not sufficient (push bar not secured in place), then turn the screw across from the quick-release fastener clamp lever (with lever in open position) counter-clockwise, until sufficient clamping pressure is achieved by closing the clamp lever. Repeat this process, as necessary.

*Adjusting the inclination angle of the back unit (including the push bar) relative to the vertical*

In addition to the adjustment options already described above, it is also possible to swivel the entire push bar back by up to 30°. This can become necessary, for example, if a seat shell is being used with a seat angle of more than 90°. This enables the gap between the back of the seat shell and the push bar to be as uniformly small as possible. This adjustment option is not intended for individual adjustment to the person pushing the chair. To adjust the inclination angle of the push bar, proceed as follows:

- ▶ Loosen the two clamp levers located in front of the pivot point for the push bar folding mechanism by turning them counterclockwise.
- ▶ Apply (axial) pressure to the clamp lever and press the coach bolts out of the guide on the ratcheting element until it is possible to freely adjust the inclination adjustment angle.
- ▶ Move the push bar into the desired position.

### Important!

Apply pressure to the rounded head in order to press the square head of the screw into the ratcheting element at the required inclination angle. 

- ▶ Check that the carriage bolt is correctly seated in the ratcheting element.
- ▶ Secure the unit in this position by turning both clamp levers clockwise to tighten them. It is also possible to maintain the inclination angle of the back unit, including the push bar unit, when the push bar is folded down (e.g. during transport).

### Trapezoid-shaped adapter receptacle and adapter

The seat shell is mounted to the undercarriage by means of an adapter system. This system consists of a trapezoidal adapter, which is mounted to the underside of the shell, and a trapezoid-shaped adapter receptacle which is connected to the undercarriage. The adapter is securely locked using two spring-loaded catch mechanisms which snap into the adapter bracket from underneath. Always check to ensure that both spring-loaded catch mechanisms are securely locked on the adapter system.

If you should discover that the spring-loaded catch mechanisms have become loose, then these should be re-tightened using an appropriate tool until the shell is securely positioned on the undercarriage.

### Tilt-in-space mechanism via two gas-pressure springs

SIMPLY and SIMPLY LIGHT are equipped with a gas-pressure spring adjustment for the inclination/tilting of the seat shell.

To adjust the tilt of the seat, proceed as follows:

- ▶ Activate the parking brakes.
- ▶ In case of only one release lever: pull the gas-pressure spring release lever mounted to the right of the push bar frame.
- ▶ In case of two release levers: pull both the gas-pressure spring release levers mounted to the right and left of the push bar frame simultaneously.
- ▶ Adjust the inclination angle of the seat shell by pulling the push bar up or pressing it down while pulling the release lever for both gas-pressure springs.
- ▶ Let go of the release lever to lock the seat shell into the adjusted inclination position.
- ▶ Deactivate the parking brake, as appropriate.

### Important!

Never pull the release lever without a firm grip on the push bar if a patient is seated in the wheelchair. Prior to inclination of the seat shell, ensure that both anti-tipping wheels always face in direction of the driving surface. In a further undercarriage variant, the tilt of the seat shell is adjusted with a telescopic unit. First the mounting of the unit is loosened and the telescopic tube moved either upwards or downwards, depending on how the tilt is to be adjusted. Then the unit is refastened. 

### Dynamic sitting

#### Locking of the AKTIVLINE® mechanism in the pelvic region incl. pelvis adjustment unit

The entire motion cycle of the shell is controlled via the physiological pivotal points of the hip and knee. The seat and back unit are mounted on sliding bars and allow opening of the hip angle to correspond to the patient's motion sequences. The knee-angle and foot rest units permit simultaneous stretching tendencies for the patient. A readjusting spring system both in the seat and the leg system returns the patient automatically into the physiological starting position.

#### Locking of the pelvic mechanism

Activate the pelvis mechanism (mobility in the hip region) by pulling the snap locks on the back of the AKTIVLINE® below the seat shell from the groove and rotating them through 90°. However, sometimes, for example during meals or transport, it is desirable to keep the patient in a passive condition. To deactivate the pelvis mechanism, rotate the snap locks through 90° and return them to their starting position until they engage and the pelvis mechanism is firmly locked.

### Important!

In case of a fault in the mobility of the seat shell, always check the pelvis mechanism first and then the foot rest mechanism!

- ▶ Pull slightly to check whether the quick pin has engaged and whether stretching of the legs is again restrained by the gas-pressure spring.

### Important!

Make sure that you do not stand or knee directly in front of the foot rest during assembly. Sudden stretching movements by the patient could result in your being kicked or injured unexpectedly.

## 7.2. Options

### 7.2.1. Lower extremities

#### Foot rest system

The foot rest system of the AKTIVLINE® must be customized to the lower leg length of the patient. During adjustment it must be ensured that the patient's feet are comfortably positioned on the foot rest. The foot rest system of the AKTIVLINE® is available in both individual as well as single-piece design, a split version must always be selected for a difference in leg length of 40 mm.

In order to adjust the foot rest system to the lower leg length of the patient, release the screw on the clamp block of the foot rest and shift the foot rest up or down, as necessary, on the centered tubing of the foot rest system; this also applies to the foot rest system with individual foot supports. Then tighten the screw after the adjustment is complete.

#### Calf support

Calf supports are already fitted ex-works.

#### Foot rest retraction unit

The foot rest retraction unit is located at the foot rest suspension and mounted under the seating unit. The stretching movement of the legs is restricted by a gas-pressure spring, i.e. the gas-pressure spring restrains stretching of the legs through resistance.

If you wish to release the foot rest retraction unit to allow free movement for the patient's legs, then proceed as follows:

- ▶ The mounting of the foot rest retraction unit is fitted to the foot rest tube.
- ▶ Release the quick pin by pressing the red button and keeping it depressed. During pressing, pull the pin from the opening.
- ▶ This now hangs loosely attached to a suspension.
- ▶ If you wish to refit the foot rest retraction unit, hold the gas-pressure spring into the eyelet, which is attached to the middle tube of the foot rest, with one hand.
- ▶ Then press the red button of the quick pin and guide this into the opening/mounting intended for this purpose.

#### Foot rest mechanism locking option

Activate the foot rest mechanism (mobility in the leg area) by turning the clamp lever on the outside of the seat shell counterclockwise until the foot rest can be moved. Make sure not to turn the clamp lever more than necessary, otherwise there is a risk of the lever popping out of the thread. To deactivate the foot rest mechanism, turn the clamp lever clockwise until the foot rest is firmly locked again.

### Important!

In case of a fault in the mobility of the seat shell, always check the pelvis mechanism first and then the foot rest mechanism!

#### Option of adjusting foot angle

The foot angle adjustment is located on the connection between the foot rest and the center tube of the foot rest. You can adjust this as follows:

- ▶ Turn the clamp lever counterclockwise to loosen it.
- ▶ Now move the foot rest into the desired position.
- ▶ Hold this position with one hand.
- ▶ Now turn the clamp lever clockwise and tighten firmly.

**Tip!**

There is a trick to using the clamp lever if you don't have enough room to turn the lever the full 360°. Turn the clamp lever as far as possible in the desired direction until it won't turn anymore and then pull it straight out. This unlocks the catch mechanism on the clamp lever and you can now turn back the lever without influencing the previously selected adjustment grade.

**Foot pedals**

The foot pedals are bolted in place onto the foot rest, these are already fitted ex-works and need not be adjusted further. However, as slipping of the legs is unavoidable for foot pedals, fixation of the feet is generally recommended.

**Heel stops**

Heel stops are already permanently fitted ex-works and no longer require any adjustment. Heel stops prevent slipping of the feet backwards.

**Foot fixation**

Foot restraints support the dynamic range of the AKTIVLINE® foot rest system. They should be gently placed around the ankles to achieve a snug fit without exerting any pressure. Tensioning of the straps is only effective if they bend the foot to 90°. This requires leg and foot guidance which enables a right-angled angle. The strap runs across the ankle joint and both ends of the strap are screwed to the foot rest. The images illustrate how the foot restraints are to be applied.

**Foot straps**

The foot straps provide for the proper positioning of the feet in order to support proper positioning of the pelvis and upper body. The foot straps relieve pressure on the ankles and improve the placement of the feet.

**Knee support**

Next to pelvic fixation, the knee support is an important part in terms of adjustment for the correct functioning of the AKTIVLINE® if a knee support is necessary. With the aid of the knee support, the patient is given the opportunity of utilizing the dynamics of the foot rest system. The dynamics of the foot rest system are supported by use of the foot restraints already described above (see section on previous page).

Observe the following instructions for the adjustment of the knee support:

- ▶ The knee support should be fitted such that it lies just below the knee cap on the patella tendon. As

the knee cap moves up and down during bending and stretching of the leg, you must ensure that the knee cap is not restricted in its movement by the knee support.

- ▶ The knee support needs to fit closely to the patella tendon but without applying any pressure.
- ▶ Make sure that the locking mechanism of the knee support engages properly and ensures a secure fit of the knee support. Depending on the model, the knee support is either fitted flexibly with a belt webbing (with one side permanently fixed to the shell and on other side with a deflection roller and fleece and Velcro) or rigidly with adjustment plates and hinges (with one side permanently fixed to the shell and the other side with a press connection), or both sides fitted with press connections. The height and depth of the knee support can be adjusted to the patient using the adjustment rails.
- ▶ To do this, loosen the screws on the right and left and adjust the knee support to the desired position, then retighten all the screws again.

**7.2.2. Pelvis****Hip belt**

The AKTIVLINE® comes equipped with a pelvic restraint, which is of decisive importance for the seat system. Only proper adjustment of the pelvic restraint enables the patient to move the AKTIVLINE® correctly, thus achieving the desired effects.

Various versions of the pelvic restraint are customized to meet the needs of the patient. There is a choice of 2-point or 4-point hip belts. As a rule, both versions are bolted in place onto the right and left of the shell or attached to the bottom shell edge with a snap buckle and clamping plate for length adjustment. Using the snap buckle and clamping plate, the length can be ad-

justed individually to the size of the patient, to do this, proceed as follows:

- ▶ Open the snap buckle on one side, adjust the belt to the desired length and reclose the snap buckle. Now proceed in the same sequence on the opposite side.
- ▶ Or: open the clamping plate with an Allen key, adjust the belt to the desired length, and then securely screw the clamping plate back again. Then repeat this procedure on the opposite side.

The hip belt itself is fitted with one, or, if necessary, two tensioning straps. These enable pressure to be applied to the pelvis, as desired or tolerated. To do this, first close the hip belt and then pull the flap, or both flaps simultaneously, to adjust the desired pressure. Please observe the following instructions when applying the pelvic restraint:

When applying the pelvic restraint, the patient must be seated as far back as possible in the seat shell, i.e. the patient's pelvis must be in contact with the back of the seat shell. In order to ensure this, it may be helpful to recline the seat shell using the undercarriage tilt-in-space mechanism until the patient slides comfortably into the desired position.

#### Important!

The patient must always be secured and under observation at all times.

#### Pelvic harness

In combination with the AKTIVLINE®, the DYNALINE® pelvic harnesses provide pelvic restraint and feature materials, which are very well tolerated by the skin and highly elastic, ensuring the greatest possible range of motion.

The standard and spread form of the pelvic harness consists of a broader part on which the patient is to sit with his/her posterior. This broad section is fixated at the back in the seat shell. In the front, the plate splits into two narrower straps. The two straps exiting from the side are pulled upwards between the patient's legs and then outwards/downwards over the groin. The straps must be spread over the thighs across their entire width.

#### Important!


If positioned incorrectly, the belt cuts on the inside and causes pain.

A T-shaped pelvic harness is positioned and attached similar to a hip belt on the right and left of the patient's pelvis, the middle part is attached between the legs. Ensure good positioning of the pelvic harness as this could otherwise lead to pressure areas.

#### Abduction pommel

An integrated abduction pommel is located under the cover and cannot be adjusted. This is always required in the case of a foot rest with individual foot supports. An adjustable or removable abduction pommel can be selected for a single-piece foot rest, proceed as follows: To loosen the abduction pommel, turn the clamp lever - mounted at the front of the abduction pommel - counterclockwise, this loosens the abduction pommel and it can simply be pulled out.

Then reinsert the abduction pommel in the intended receptacle and turn the clamp lever clockwise until tight. Then check whether the abduction pommel has its desired position again.

#### Side cushions for the seat, attachable with Velcro

Side cushions with Velcro attachment serve to adapt to growth. These can be removed by a strong pull of the cushion. Zippers allow removal of the covers for washing. Please follow the exact washing instructions provided for washing the relevant materials.

### 7.2.3. Upper body

#### Chest belt

The chest belt is a belt which is attached at each end of the shell below the patient's armpits and is closed at the front. Various variants of chest belts can be selected. The chest belt embraces the upper body (ribcage) and prevents the user from falling off the seat.

#### Chest belt with flap

The chest belt with flap has been mounted in the proper position during adjustment of the seat shell. It is possible to perform minor adjustments to the length using the hook and loop fastener.

To open the belt, proceed as follows:

- ▶ Pull the front flap upwards.
- ▶ Then pull the overlapping ends of the chest belt apart.
- ▶ You can then lay the two ends of the chest belt aside and remove the patient from the seat shell, as necessary.
- ▶ Repeat these steps in reverse order when applying the belt to secure the patient in the seat shell.
- ▶ When applying the belt, check to ensure that it is not constricting the patient; this is especially im-

portant to ensure that the patient's breathing is not adversely affected.

#### *Chest belt with deflection roller*

The chest belt with deflection roller has been mounted in the proper position during adjustment of the seat shell. It is possible to perform minor adjustments to the length using the hook and loop fastener.

To open the belt, proceed as follows:

- ▶ Pull the overlapping ends of the chest belt apart.
- ▶ Now pull the loose end of the chest belt through the metal loop.
- ▶ You can then lay the two ends of the chest belt aside and remove the patient from the seat shell, as necessary.
- ▶ Repeat these steps in reverse order when applying the belt to secure the patient in the seat shell.

#### *Chest harness*

The chest harness is used to restrain the patient's entire upper body. The chest harness is tailored individually to the needs of the patient. To avoid pressure areas, the fasteners are always sewn to the harness ends. This makes slipping of the fasteners to the side impossible.

The lower end of the harness runs across the two spines (the two upper, anterior iliac spines) to the side of the pelvis. The middle part of the chest harness is guided across the torso at the front and ends at the height of the sternum. The two upper harness ends pass along the side of the neck over the shoulders.

The chest harness is an optimal correction aid to achieve symmetrical alignment of the pelvis and the torso. At the same time it continues straightening the pectoral girdle and the neck. By placing the chest harness symmetrically on the spines (straps pull to straighten the pelvis), the pull upwards helps to straighten the shoulders and to symmetrically stretch the spine, which keeps the head erect.

#### **Important!**

Always adjust the chest harness, working from the hard back of the shell outwards and using the tighteners located there. Restraining the patient in the seat shell.

**Precondition:** place the patient in the seat shell and lay the chest harness over the back edge of the shell.

To restrain the patient, proceed as follows:

- ▶ Position the patient correctly in the seating shell with the pelvis on the inner back of the shell.
- ▶ Lay the chest harness facing forwards. The head is then positioned between the two shoulder straps.
- ▶ The lower end of the chest harness is equipped with two tighteners for fastening the straps.
- ▶ Open the flap on one of the black strap tighteners (if this version applies).
- ▶ Now push the belt webbing through the seat cushion and the back cushion into the slot on the corresponding side between the moving flap and the fixed part.
- ▶ Pull the belt webbing as far as it will go.
- ▶ Hold the belt webbing in this position.

- ▶ Then close the flap on the tightener until it lies flat against the fixed part.
- ▶ Or use the plug lock, if available, to close the strap.
- ▶ To open the lock, pull the flap on the tightener upwards and pull out the belt webbing. Or open the plug lock with two fingers.
- ▶ Repeat this procedure on the other side.
- ▶ Repeat these steps in reverse order to release restraint by the chest harness.

Some chest harnesses are fitted with tensioning straps, proceed here as follows.

- ▶ Close the clip buckles or plug locks on both sides.
- ▶ Now pull on the tensioning straps on the right and left to adjust the desired pressure.
- ▶ Always start first in the pelvic region and then in the region of the shoulders.

Through correct guidance of the slots in the seat, the strap is pulled outwards/downwards and thus stabilizes the pelvis in central position. The pull direction of the straps should be an angled extension of the harness upwards/to the back (erection of shoulders).

To adjust the length of the straps of the chest harness, proceed as follows:

- ▶ Adjust the lower ends of the chest harness first.
- ▶ Open the lower tighteners on the hard shell back or under the seat. To do so, simply pull hard on the free end of the belt webbing and the tightener will open automatically.
- ▶ Pull the belt webbing to the desired length on one side.
- ▶ Hold the belt webbing in this position.
- ▶ Then close the flap on the tightener until it lies flat against the fixed part of the tightener.
- ▶ Repeat this procedure to adjust the other side.
- ▶ Repeat these steps to adjust the upper ends of the chest harness.
- ▶ Adjust these via the upper tighteners on the hard back of the seat shell. The procedure is identical as for the lower ends of the chest harness.

#### **Important!**

Ensure correct strap arrangement. Straps arranged too closely around the neck lead to a risk of strangulation. If the straps are arranged too loosely around the shoulders, the straps could slip off the patient's shoulders; the patient is then no longer secured and could fall out of the chair, resulting in possible injury.

#### **Important!**

The patient must always be secured and under observation at all times.

#### *Strap guide plates*

Strap guide plates are used to optimally adjust the chest harness for each patient. These allow individual adjustment of the strap guide. The strap guide plates are already mounted ex-works to the upper, inner edges of the back.

#### **Important!**

If strap guide plates are present, always use these for correct strap guidance. To

do this, pull the straps on the right and left through the provided slot in the center of the plates and then lock the strap in the clip buckle. Make sure that the straps are not twisted or knotted.

The correct height of the straps is to be adjusted as follows:

- ▶ Loosen the two screws on the strap guide plates with an Allen key.
- ▶ Push the plate upwards or downwards into the best position.
- ▶ Tighten all screws again with an Allen key.
- ▶ Check strap arrangement on the patient.
- ▶ Readjust, as necessary.

**Important!**

Ensure correct strap arrangement. Straps arranged too closely around the neck lead to a risk of strangulation. If the straps are arranged too loosely around the shoulders, the straps could slip off the patient's shoulders; the patient is then no longer secured and could fall out of the chair, resulting in possible injury.

*Reclining supports with Velcro fastening*

Reclining supports are shaped to the patient, starting from the top. These pass from the edge of the shell over the collar bone. The pressure of reclining supports passes from the upper ribcage parallel to the sternum further outwards, upwards and to the back. This pressure helps the patient when erecting the pectoral belt. Reclining supports are screwed firmly to the back of the shell and are rotated sideways using the lock, this allows for easy transfer of the user.

*Lateral side supports*

Lateral side supports support the torso symmetrically, or displaced the spine accordingly, for example, in case of an S-scoliosis. Various versions of lateral side supports allow individual adaptation to the patient's needs. For example, permanently fitted lateral side supports are customized to the patient's need on delivery. Adjustable lateral side supports can be fitted directly on site. The mounting for the adjustable lateral side supports is on the hard edge of the rear back part.

It is best to adjust the lateral side supports with the patient seated in the chair. In order to achieve an ideal level of support, adjust one side first and then the other.

To do so, proceed as follows:


- ▶ Loosen the two screws on the rear plates of the side supports.

- ▶ Now adjust the lateral side supports to the patient's measurements. The lateral side supports can be adjusted both for height as well as width.
- ▶ Tighten all screws again.
- ▶ Proceed in the same manner on the opposite side.
- ▶ Check the fit of the lateral side supports and make any necessary changes.

*Lumbar spine support with Velcro attachment*

Support of the lumbar spine leads to better erection of the torso. Various variants are possible: for one, an integrated lumbar spine support which is already integrated into the back part on delivery, or alternatively, a lumbar spine support with Velcro attachment. Lumbar spine support with Velcro attachment can be both on top or underneath the cover.

Lumbar support allows lordotic bending of the spine in the lumbar region. This type of support is particularly suitable for patients who tend towards a kyphotic back posture and require support in torso erection. The thickness of lumbar spine support is determined during preparation of the concept.

To adjust the lumbar spine support correctly, proceed as follows:

- ▶ Lumbar spine support under the cover: open and remove the cover of the back part.
- ▶ Detach the lumbar spine support from the Velcro fleece connection.
- ▶ Place the lumbar spinal support in the correct desired position and attach again with Velcro.
- ▶ Check the fit of the lumbar spine support and make any necessary changes.
- ▶ Reclose the cover of the back cushion.

#### *Back extension, attachable by slipping on*

If the back length no longer provides adequate due to the patient's growth, then slip-on back extension is used. This is slipped onto the rigid back and is located between the back part and the head rest, thus extending the back height by the length of the back extension.

Proceed as follows:


- ▶ Take the head rest from its mounting (see section "Head rest").
- ▶ Then insert the back extension with the hook facing downwards behind the cover of the back part.
- ▶ Now insert the head rest back into the head rest bracket through the back extension into the actual back unit.
- ▶ Then adjust the head rest correctly and secure all clamp levers.

## 7.2.4. Upper extremities

### *Arm rests*

The AKTIVLINE® seat system can be fitted with different arm rests of the ALULINE range. Arm rests give the patient sufficient support and stability and support erection of the torso. Mount the arm rests by releasing the locking levers of the arm rest holders on the right and left sides of the seating shell, insert the arm rests and tighten the locking levers again. To adjust the height, change the insertion depth on the oval tube universal holder.

To adjust the angle of the arm rest, proceed as follows:

- ▶ Open the clamp on one side of the arm rest bracket by turning the clamp lever counterclockwise.
- ▶ The angle of the arm rest to the seat shell can then be adjusted.
- ▶ Hold the arm rest in the desired position with one hand and secure it in place by tightening the clamp lever clockwise with the other hand.
- ▶ This will maintain the arm rest in position.
- ▶ Repeat this procedure on the other side.

### *Upper arm supports*

The upper arm supports serve to limit the range of motion. The patient is no longer able to move his or her arms backwards, between the back rest and arm rest. The standard equipment includes the upper arm supports mounted on the back rest.

In the option with adjustable upper arm supports, these can be adjusted in height and depth via an adjustment rail on the back of the seating shell. To do so, proceed as follows:

- ▶ Loosen the screws on the adjustment rail of the upper arm supports.
- ▶ Slide the upper arm supports into the desired position.
- ▶ Tighten all screws again.

### *L-shaped arm rests*

L-shaped arm rests meet both the function of an arm rest as well as the function of an upper arm support. This type of arm rest follows the stretching of the patient and automatically returns to the starting position of the seating shell. Please observe the adjustment options and the correct positioning, as already described under arm rests and the upper arm supports.

To adjust to a good position, please proceed as follows:

- ▶ Loosen the clamp lever on the arm rests and the screws on the rear side of the upper arm support.
- ▶ Then bring the arm rests into the desired L-shaped position by pushing up or pressing down.
- ▶ Now hold the desired position with one hand and simultaneously retighten the clamp lever of the arm rest.
- ▶ Tighten all the screws on the upper arm support again.
- ▶ Proceed in the same manner on the other side.

### *Lower arm restraint*

Lower arm restraints are combined with arm rests. The arm rests can be adjusted individually in height and inclination, so that the lower arm of the patient can be placed on them and secured with straps. Restraining of the arms lends greater stability to entire upper body. There is also the possibility of a lower arm restraint with expander. In this variant, an expanding rubber is mounted underneath the arm rest and restricts the patient's radius of arm motion. A lower arm restraint with expander is mainly used for the patient's own protection and to avoid injuries to the arms. Mainly for patients who strike out or grab uncontrolled with their arms and can thus knock over objects and injure themselves in the process. Arm restraints can also be used to provide better stability and guidance for operating a control unit on electric wheelchairs.

To apply the the lower arm restraint, proceed as follows:

- ▶ Lower arm restraints with expander rubber: open the Velcro fleece connection on the lower arm restraint and place this around the lower arm of the patient and then close again. **IMPORTANT!** Do not exercise too much force on the patient's lower arm.
- ▶ Lower arm restraints without expander rubber: open the Velcro fleece connection and place the patient's arm comfortably on the arm rest, then close the lower arm restraint around the patient's lower arms and the arm rest.

**Important!**

Do not exercise too much force on the patient's lower arm.

**Important!**

Ensure that the arm restraint is not too tight. Risk of pressure areas.

*Hand cuff restraints*

Hand cuff restraints provide for direct restraint of the hands. They can be wrapped around the hands in order to prevent patients from injuring themselves with the fingers.

**Important!**

Ensure that the arm restraint is not too tight. Risk of pressure areas.

*Therapy tray table*

Various shapes and materials were taken into consideration in the selection of the therapy tray table. All variants used can be fastened differently, nonetheless the AKTIVLINE® has already been fitted with a mounting ex-works and this only needs to be adjusted to meet the patient's needs.

*Variant with slip-on mechanism*

The holders on the therapy tray table, two tubes on the right and left of the table tray, as well as the receptacles for these tubes which are mounted under the arm rests, are already pre-fitted on delivery. To apply the therapy tray table, proceed as follows:

- ▶ Insert the therapy tray table and the two tubes fitted to the therapy tray table resp. into the intended receptacles under the arm rests.
- ▶ Make sure that the therapy tray table is not pushed too close to the patient. Leave the patient sufficient space in the stomach region to avoid sitting trapped during stretching.
- ▶ You can adjust the height of the table via the arm rests.
- ▶ Mount the table in the correct position by turning the wing screws on the mountings under the arm rests clockwise.
- ▶ Repeat these steps in reverse order to detach the therapy tray table. Loosen the wing screws by turning counterclockwise and pull the therapy tray table from the holder away from the patient.

*Variant with swing-away unit*

To attach the therapy tray table, proceed as follows:

- ▶ There is a receiving tube under the table. Insert this in the intended holder under the arm rest.
- ▶ Make sure that the therapy tray table is not pushed too close to the patient. Leave the patient sufficient space in the stomach region to avoid sitting trapped during stretching.
- ▶ You can adjust the height of the table via the arm rests.
- ▶ Mount the table in the correct position by turning the wing screw on the mounting under the arm rest clockwise.
- ▶ Repeat these steps in reverse order to detach the therapy tray table. Loosen the wing screw by turning counterclockwise and pull the therapy tray table from the holder away from the patient.

- ▶ Once you have adjusted the therapy tray table into the correct position, you no longer need to remove it from the holder for transfers and the like.
- ▶ In the variant with swing-away unit, simply fold the table tray to the side. Always pay attention to the patient. Should the radius for swiveling not be sufficient, loosen the wing screw and pull the table tray slight forwards to increase the radius.

*Variant with dynamic co-motion,**with Velcro fleece connection and strap:*

In the variant with dynamic co-motion, a Velcro fleece connection is fitted to the underside of the therapy tray table and the arm rests.

- ▶ Simply place the therapy tray table in its correct position on the arm rests and press gently.
- ▶ Now close the strap which is fitted to the sides of the therapy table tray by pressing the strap ends of the therapy tray table onto the strap ends of the receptacle mounted to the back.
- ▶ Adjust the correct height of the therapy tray table via the arm rests.
- ▶ Repeat these steps in reverse order to detach the therapy tray table, release the straps and simply pull the table tray from the arm rests in your direction.

**7.2.5. Head***Head rest systems*

Once the patient has been restrained in the correct position in the AKTIVLINE® with the hip belt, adjust the head rest. The holder for the head rest is built into the back of the shell on the AKTIVLINE® and was already fitted before delivery.

The AKTIVLINE® seat system can be fitted with different head rest systems of the ALULINE range (straight or shell-shaped). Both can be adjusted in height and depth. Mount the head rest by loosening the securing lever at the head rest holder located in the upper part of the back rest, insert the head rest bracket and retighten the lever. To adjust the height, change the insertion depth on the oval tube universal holder.

To adjust the depth, proceed as follows:

- ▶ Loosen the clamp lever at the upper oval tubing by turning the clamp lever counterclockwise.
- ▶ Shift the tube of the head rest bracket in the universal bracket until the desired position has been reached.
- ▶ Secure it into position by turning the clamp lever clockwise until it is no longer possible to slide the head rest bracket tube.

Head support with continuous angle adjustment:

Release the clamp lever on the top end of the support bar on the head rest to adjust the height. Adjust the inclination of the head rest cushion using the continuous angle on the head rest frame.

**Tip!**

There is a trick to using the clamp lever if you don't have enough room to turn the lever the full 360°. Turn the clamp lever as far as possible in the desired direction until it won't turn anymore and then pull it straight out. This

unlocks the catch mechanism on the clamp lever and you can now turn back the lever without influencing the previously selected adjustment grade.

It is essential to observe the following instructions for the adjustment of the head rest:

- ▶ Standard head rests should be adjusted to fit directly behind the patient's head.
- ▶ Shell-shaped head rests usually sit lower on the head than standard head rests. In other words, the underside of the shell-shaped head rest should ideally sit at the base of the head, meaning the neck region.

#### Standard-shaped head rest:

This head rest is fitted with a core aluminum plate. The head rest can be bent individually to the needs of the patient. In case of pronounced asymmetry, the usual U-shape may sometimes require stronger bending on one side to counteract the sideways inclination of the head. The head rest is already supplied curved ex-works, if you wish to alter the setting, please contact your local specialist dealer.

The shell-shaped head rest virtually positions the patient's head automatically in the center, due to the shaped side and centering in the recess, the head automatically "slips" into this position.

Ensure that the side elements are not bent too much to avoid exerting pressure on the head. Leave the patient enough space of movement to allow him/her to perceive his/her surroundings adequately and ensure that the patient's ears are not trapped or squashed.

#### Warning!

The head rest is always to be adjusted such, that the gap between the head rest and the back has a maximum breadth of two fingers. There is a risk of strangulation if the patient's head slips between the head rest and the back!

#### Caution!

The gap between the head rest (bottom edge) and the upper edge of the back rest must be less than 60 mm (for children) or 120 mm (for adults) or greater than 250 mm according to DIN EN 12182.

#### Warning!

The combination of the AKTIVLINE® with third party head rests (head rest systems not manufactured by interco) is not permitted. There is a risk of injury!


### 7.2.6. Undercarriage

*AKTIVLINE® seat system, separable from the undercarriage*

#### *Attaching the seat shell (if option exists)*

The seat shell is attached to the trapezoid-shaped adapter receptacle on the undercarriage using the trapezoidal adapter. To do this, perform the following work steps and check secure mounting on the trapezoid-shaped adapter receptacle after placing the seat shell. First release the seat shell fixation by pulling out the snap lock and turning it by 90 degrees. Position the back end of the trapezoidal adapter to the front side of the trapezoid-shaped adapter receptacle so that the lock is pressed downwards.

Now push the set shell to the back until this engages audibly. Relock the seat shell fixation by rotating the snap lock through 90 degrees until it engages. Finally, check the seat shell and the seat shell locking system to ensure they are secure.

Now fix the security pin which connects the back unit of the seat shell to the undercarriage.

#### Important!

Always check whether the seat shell is firmly fitted to the trapezoid-shaped adapter receptacle. Check this by pulling the seat shell and simultaneously pressing on the undercarriage. The seat shell may not slip out.


### Removing the seat shell

First, pull out the two snap locks and turn through 90 degrees. Remove the security pin which connects the back unit of the seat shell to the undercarriage. You will now be able to move the seat shell system forward and detach it from the adapter. Before reinstalling the shell, remember to turn the locking pins 90° again in order to reset them in the 'closed' position. Be careful that the security pin connecting the back unit with the undercarriage is again moved into place and that it is correctly set and adjusted.

### AKTIVLINE® undercarriage, foldable

**Folding mechanism for transporting the SIMPLY LIGHT**  
SIMPLY LIGHT is equipped with a folding mechanism in order to make transport easier. The back unit including push bar, is folded forward; it is also possible to remove the back wheels.

- ▶ To do this, gently step on the black cross bar at the rear back part and simultaneously press the push bar forwards.
- ▶ Now you can remove the wheels.

### AMF force application point system

See chapter on "Transport"

### Anti-tip

The anti-tip is already factory-mounted. It is not possible to perform any readjustments.

### Wheel protectors

The wheel protectors cover the upper part of the wheels. They consist of a plastic mudguard which is screwed to an aluminum mounting plate. For fitting to the chassis, the free end of the mounting plate is inserted into the upper slot of the wheel mounting. To mount the wheel protector, proceed as follows:

- ▶ Take the wheel protector and place the groove-shaped part of the plate into the upper groove of the axle block.
- ▶ Insert the wheel protector to the stop. You can feel a slight click when the plate engages in the guide.
- ▶ Repeat this procedure on the other side.
- ▶ For disassembly, pull the wheel protector off in upward direction.

### Spoke protectors

The spoke protectors are already factory-mounted. It is not possible to perform any readjustments. If it is necessary to perform repairs to the spoke protectors, please contact your medical supplies dealer.

### Quick-release axles for front wheels

To remove the front wheels, proceed as follows:

- ▶ Remove the seat shell (if this option is given).
- ▶ Fold the SIMPLY LIGHT undercarriage (if this option is given), but without activating the parking brake.
- ▶ Lift the SIMPLY LIGHT slightly at the front so that the front wheel no longer touches the floor.
- ▶ Press on the rubber cover of the guide tube with your thumb. This activates the locking button located under the cover.
- ▶ Keep the button pressed down and pull the front wheel off with your free hand.
- ▶ Repeat this procedure on the other side.
- ▶ Repeat these steps in reverse order to assemble the front wheels (e.g. after transport). The wheels must engage audibly.

### Locking of the front wheels

Release the lock of the front wheels by pulling the snap lock out of the groove and turning through 90°. To reactivate the lock of the front wheels, turn the snap locks back into their initial position until they engage. The front wheels are then locked and the wheelchair can only be moved straight ahead.

#### *Knee lever brake (air brake)*

Knee lever brakes are used in addition to parking brakes (drum brakes). These are factory-mounted and adjusted with your order. If the braking action is diminished by changes in distance (between wheel and the knee-lever brake), you can readjust this by shifting the sliding nuts on the C-rail.

#### *Anti-tipping wheels*

The rear end of the frame is equipped with two anti-tipping wheels.

#### **Important!**

In normal driving position, the anti-tipping wheels must always point in the direction of the driving surface!

It may be necessary to swivel the anti-tipping wheels upwards when negotiating high curbs and for transport (reduces the size of the wheelchair for transport). To do so, proceed as follows:

- ▶ Activate the parking brakes.
- ▶ Position yourself behind the SIMPLY or SIMPLY LIGHT undercarriage.
- ▶ Pull on one of the anti-tipping brackets first and slide it partly out of the adapter bracket groove.

- ▶ Then move the anti-tipping bracket upwards until it engages in the upper mounting groove.
- ▶ Now allow the anti-tipping bracket to slide back into the groove of the adapter bracket unit. Repeat the same steps with the second anti-tipping wheel.
- ▶ Release the parking brake.
- ▶ Follow the same procedure for returning into the original position.

#### **Note!**

Always check to ensure that the two anti-tipping wheels are always facing in the same direction. This is the only way to ensure effective protection against tipping!

#### *Stabilizer bar for AKTIVLINE® top frame, able to be fixed in position*

The stabilizers serve to provide greater stability through additional fixation of the desired sitting position. Fixation of the stabilizers is neutralized by releasing the two lower clamp levers. Now the two mini activation handles on the push bar can be activated and tilting of the seat executed. If you now pull the lower clamp lever again, this additionally ensure fixed setting of the tilt.

### 7.2.7. Custom items

So


#### *Seat depth regulation in 5 cm increments via infinitely adjustable seat plate*

Seat depth regulation should be adjusted exclusively by your local specialist dealer or by us!

#### *Locking of the AKTIVLINE® mechanism in the pelvic region via hydraulic thrust absorber*

The hydraulic thrust absorber has the task of regulating the dynamic seat movement of the seat shell. The hydraulic thrust absorber is fitted with two adjusting screws for adjusting the stroke of the absorber.

Proceed as follows for adjustment:

- ▶ By loosening both adjustment screws (by turning to the left), the seat shell can be stretched by the patient's uncontrolled muscle tension and returns into its original position after tension eases.
- ▶ By tightening the upper adjustment screw (by turning to the right), this allows stretching of the shell, but an automatic return to the original position is no longer possible. By slowly releasing the lower adjustment screw, the seat shell slides back into its original position.
- ▶ To become familiar with the function of the hydraulic thrust absorber, you should try out the various setting options.

### 7.3. Optional accessories

#### 7.3.1. Lower extremities

##### *Rubber pads*

If a foot rest retraction unit has not been chosen, then rubber pads are used. These are already fitted ex-works and cannot be adjusted. Two rubber pads are located below the seating unit, these prevent the foot rest from hitting against the material in the rest phase after stretching.

#### 7.3.2. Pelvis

##### *Pelvic bracket*

The mounting of the pelvic brackets on the shell edge and the shape of the pelvic bracket exert pressure on both sides in downward/outward direction for the purpose of symmetrical erection of the pelvis.

##### *Please observe the following instructions when using a pelvic bracket*

The pelvic bracket pressure pads must support the spine (both the superior and anterior iliac spine). Make sure that the pelvic bracket pressure pads do not exert pressure on the abdomen or other soft tissue regions! Thread the fastening strap through the buckle on the pelvic bracket and tighten it until the patient is securely seated in the prescribed position (pelvis against the back of the shell). If this makes the patient uncomfortable, then loosen the pelvic bracket somewhat. Keep in mind, however, that the shell can only be opened easily if the pelvis is in constant contact with the pelvic restraint. It may be necessary to gradually increase the patient's tolerance of the restraint by selecting a lower starting pressure and then increasing it over time until the strap is in the correct position.

#### 7.3.3. Upper body

*Custom-made chest harness*  
See Adjusting the chest harness

#### 7.3.4. Upper extremities

*Custom-made hand restraint*  
See Adjusting the hand restraint.  
*Custom-made therapy tray table*  
See Adjusting the therapy tray table.

#### 7.3.5. Head

Head rest bracket, reinforced version  
In case of very strong patients, a reinforced version of the head rest bracket should be selected, to make it more resilient. For information on adjustment, please refer to the section "Adjustment of the head rest".

#### 7.3.6. Undercarriages

*Split or custom-made push bar*  
See Adjustment of the push bar.

##### *Widened wheel base*

Spacer blocks are used to prevent the front wheels from hitting the foot rest under certain conditions during driving. These are already fitted permanently ex-works and cannot or need not be adjusted.

##### *Solid tires*

These tires are solid and already pre-fitted on delivery. Adjustments cannot be made.

### 7.3.7. Custom items

#### Hook for tube-fed nutrition

The hook for tube-fed nutrition is located at the rear cross-section of the push bar. This can be shifted easily to the right or left.

#### Important!

Use the hook exclusively for attaching tube-fed nutrition! Not intended for rucksacks, carrier bags, etc. – risk of tipping.

## 8. Combination options

- 8.1. Combination options with pushing aids
- 8.2. Combination options with stair climbing aids
- 8.3. Combination options with electric wheelchairs

The combination of AKTIVLINE® seating care systems with other products not manufactured by interco is not permitted. interco Group has made arrangements

with some other manufacturers in terms of combining various additional options, for example, pushing aids or stair-climbing aids.

Please check the individual combination options, as a pre-check is always necessary to to custom manufacturing.

## 9. Technical data

#### Important!

The maximum load (workload) is 70 kg for the AKTIVLINE® seat system.

When exceeding the permissible patient weight, we wish to point out that the AKTIVLINE® seat system may possibly require a reinforced frame and/or reinforcement of other components. Discuss this option with an authorized specialist dealer.

As this medical device - the AKTIVLINE® seat system - is custom-made for a named patient, all technical data may vary from the data given below and must be checked for each case. The AKTIVLINE® seat system is adapted exclusively to the patient's needs and custom-made according to the measurements provided for the patient.

#### Measurements and weights

Total weight (complete seat-back unit with undercarriage and options)	from approx. 40 kg to approx. 80 kg
Overall width (incl. undercarriage)	from approx. 620 mm to approx. 800 mm
Overall height (incl. head rest)	from approx. 1000 mm to approx. 1500 mm
Total length (including supporting rollers and foot rest)	from approx. 1000 mm to approx. 1400 mm
Total weight (only seat unit)	from approx. 25 kg to approx. 45 kg

#### Note!

The data given is approximate only, as the AKTIVLINE® seat system is custom-made.

Tire type	Air pressure max. [bar]
Front wheels, Ø 200 x 50 mm	2.5
Rear wheels 12 inch	2.5
Rear wheels 16 inch	2.5
Rear wheels 22 inch	4.0
Rear wheels 24 inch	4.0

#### Note!

Please note that in choosing the appropriate tire pressure you will need to reach a compromise between riding comfort for the patient (low air pressure) and pushing comfort for accompanying persons (high air pressure). It may be necessary to experiment here. The factory-assembled tires are delivered with an air pressure of 2.5 bar.

## 10. Transport/transport lock/climbing stairs

### 10.1. General information on transport

Only the tested and approved wheelchair and occupant restraint system from AMF-Bruns (Art. No. 75020) may be used for transporting the AKTIVLINE® seat system

in motor vehicles for the transportation of the disabled. In addition, a 4-point restraint system from AMF-Bruns (Art. No. 75021) is required for transport in motor vehicles. The AKTIVLINE® seat system thus complies

with the requirements of DIN 75078. For information on how to correctly secure the force application point system in the respective vehicle, please see the separate manufacturer's instructions and the following chapter 10.2 "Transport with an AMF force application point system".

**Warning!**

Never secure the transport lock to the frame or the back of the AKTIVLINE® seat system!

**10.2. Transport with an AMF force application point system**

The AKTIVLINE® seat system was subjected to a crash test as standard version with maximum fittings in accordance with ISO 7176-19. This model is a custom-made version and cannot be tested in a crash test. Nonetheless, the custom-made model is also fitted with a force application point system manufactured by AMF Bruns. Further information is available at [www.amf-bruns.de](http://www.amf-bruns.de).

The wheelchair is marked on four points which indicate the mounting points for the tension belts for transportation with motor vehicles. Only use the provided mounting points for strapping with tension belts. When fitted correctly, the head rest of the AKTIVLINE® seat system is very robust, but cannot replace head rests fitted to motor vehicles.

The AKTIVLINE® seat system has been testing without an additional drive. If an additional electrical drive is to be added, then this accessory is to be disassembled if the AKTIVLINE® seat system is to be transported in a motor vehicle.

**Important!**

Please note that the dynamic system in the pelvic area (see chapter 7.1. Standard equipment – Fixing the pelvic mechanics) must be closed during vehicle transport.

Please also observe the AMF Bruns instructions for use.

If you have any questions about vehicle transport or general questions about the product, please contact your rehabilitation specialist or your interco Group.

**10.3. Climbing stairs**

Stairs or high curbs are only to be negotiated with the assistance of two accompanying persons. Assistants are only permitted to grip the front frame tube and the back brace of the frame, without lifting or carrying the AKTIVLINE® seat system. The anti-tipping device must be retracted and secured. When ascending, the assistant behind the AKTIVLINE® seat system pulls the wheelchair up the step. The assistant at the front stabilizes the position by pressing the AKTIVLINE® seat system evenly into the steps on both sides. When descending, the assistant at the front brakes by pressing the AKTIVLINE® seat system evenly into the steps. The assistant behind holds the back brace, and stabilizes and maintains the wheelchair in the correct position. To prevent tire damage, slowly roll the AKTIVLINE® down the steps. The anti-tipping system must then be properly readjusted.

**Caution!**

Do not hold the AKTIVLINE® seat system by the push bar when negotiating stairs.

Never lift the care system by the foot rest or the foot supports or the arm rests or other attached parts.

Do not lift the care system with the user sitting in it.

**10.4. Transport on aircraft**

The AKTIVLINE® seat system can only be transported without the patient in the cargo hold. There are no concrete instructions for transportation on an aircraft. It is essential to contact the airline beforehand to determine individual transportation options.

## 11. Possible faults

In case of faults, please immediately contact your local specialist dealer or contact us at interco Group.

Your local specialized dealers:

**interco Group GmbH**

Im Auel 50  
53783 Eitorf

Tel.: +49 2243 8807-0  
Fax: +49 2243 8807-29  
E-mail: [info@interco.gmbh](mailto:info@interco.gmbh)  
[www.interco.gmbh](http://www.interco.gmbh)

If you have questions, you are also welcome to visit our Internet homepage at [www.interco.gmbh](http://www.interco.gmbh).

## 12. Material/cleaning

---

The AKTIVLINE® medical device must always be stored in a dry place to prevent possible corrosion. Please keep to ambient temperatures of -10 °C to +40 °C. Sand or other dirt particles should not be allowed to penetrate the product. This could damage the wheel bearings or mechanics. Use in salty water is not permissible. The cushions can be removed for cleaning and care if this option was indicated when ordering.

### 12.1. Materials used

#### 12.1.1. Powder-coated aluminum

All AKTIVLINE® undercarriages are covered with a durable plastic coating. Some attachments and the undercarriage are made of aluminum and then finished with a plastic coating. Powder varnishes are used for coating. Powder varnishes are solvent-free coating materials, composed of different resin systems, pigments as well as special additives.

#### 12.1.2. Fabrics

##### *Synthetic leather*

interco uses easy-care designer synthetic leather for Dynaline restraints. With its outstanding quality characteristics, this synthetic leather is stylish and ideal for all upholstery applications and coverings. Especially excellent features include long-term durability, a comfortable seating experience, resilience to wear, good hold, and tear resistance. In addition, the synthetic leather is very elastic, resistant to urine, blood and disinfectants, free of azo compounds, CFC, formaldehyde, as well as free of PCTs, PCPs and PCBs.

##### *Spacer fabric*

The spacer fabrics used by interco have two outer textile layers, which are connected by a number of threads and spaced at a specific distance; due to their special structure, they are also called 3D textiles. The thermoregulating pile layer wicks moisture away and shields against cold. 3D fabrics yield to pressure like a spring. The material features prolonged durability and pronounced elasticity. Due to the special 3D structure, there is a continuous insulating layer of air between the outer layers. The spacer fabrics feature high contracting force, breathability, and compression elasticity. They compensate well for temperature fluctuations and the material is lightweight, tear-resistant and easy to care for.

##### *Neoprene*

The neoprene used by interco is a foam based on chloroprene rubber (CR). Unique material properties are achieved in the compound due to a specific structural mixture of various types of rubber as the base substance and numerous other components. These include, in particular, excellent insulation properties, high dynamic load capacity, resistance to weather and ozone, good protection against aging, resistance to wear and tear, pronounced elasticity, and impact resilience, good adhesion properties, resistance to sea water and many chemicals, resistance to fungi and bacteria, and low gas permeability. Neoprene is also environmentally friendly, very well tolerated by the skin, and impermeable to water.

##### *Trevira CS (fabrics)*

Trevira CS fabrics are easy to clean, dry quickly and keep their shape. Unlike other fabrics, they require considerably less water for washing, less detergent, and less energy. Drying times are also short thanks to low water absorption. All in all, these are valuable economic benefits. Trevira CS is incomparably superior to flame-retardant cotton fabrics and is easy to disinfect in hospital settings, for example. Flame-retardant Trevira CS fabrics meet all major international fire safety standards.

All materials used are flame-retardant.

### 12.2. Cleaning of materials

#### 12.2.1. Powder-coated aluminum

All AKTIVLINE® undercarriages are covered with a durable plastic coating. All parts can be cleaned with soapy water and a soft sponge. Avoid using aggressive or abrasive cleaning agents, as these can damage the surface of the plastic coating.

#### 12.2.2. Fabrics

You can play a role in maintaining the appearance of the valuable covers of your seating system by cleaning, vacuuming or carefully brushing them regularly before stains set in, by wiping up spilled liquids right away with an absorbent cloth – the longer you wait to do so, the more difficult it will be to clean. Usually it is helpful to apply a small amount of mild liquid detergent for delicate fabrics to a cloth and dab it on the spill (never rub vigorously). Important: Never apply detergent directly to the cover; always apply it to a cleaning cloth. Wipe gently with clean water, if necessary.

Please note: bleach or whiteners should never be applied to the material. The material should be exposed to as little direct sunlight as possible.

All of the fabrics for our covers are washable at 30°C. Use a detergent for delicate fabrics or a gentle dry-cleaning process to clean the covers.

Never dry the covers in the dryer! Iron at a low temperature [single dot]. If you should wish to wash the seat shell cover frequently, then it is advisable to have an auxiliary cover to put on the seat shell while washing the other cover. Ask your local specialized dealers on this topic. Please clean accessories with a conventional mild household soap solution.

### 12.3. Disinfection of materials

#### 12.3.1. Powder-coated aluminum

Please use a water-based agent for disinfection, e.g. Terralin, Quartamon Med or Sagrotan Original concentrate. Please observe the manufacturer's instructions for use when doing so. It is best to clean plastic parts with a commercially available plastic cleaner (please observe the product information for the respective cleaning agent).

**Note!**

If it is necessary to perform wipe disinfection of the AKTIVLINE® seat system, then use a suitable disinfection agent according to the DGH disinfection list for surfaces.

All cleaning and disinfection agents must be rinsed thoroughly and the cleaned item dried fully before use. If this is neglected, this can lead to chemical reactions making the material brittle or lead to skin irritations of the patient.

**12.3.2. Fabrics**

Branded disinfection agents may be used according to the instructions of the fabric manufacturer.

Please use a water-based agent for disinfection, e.g. Terralin, Quartamon Med or Sagrotan Original concentrate. Please observe the manufacturer's instructions for use when doing so.

**Important!**

Synthetic leather covers/skai covers may not be cleaned with disinfecting agents, as the material gets porous over time and cracks.

**13. Maintenance /inspection by the user**

Prior to being handed over to the patient, the AKTIVLINE® seat system is checked for functionality and freedom of faults by trained and qualified employees. To ensure safe use of the AKTIVLINE® seat system over its entire service life and to ensure a perfect technical condition, it must be maintained and inspected regularly according to the following maintenance and service plans.

Incorrect settings by the user or laypersons bear the risk of injuries or crushing, which in a worst case scenario could lead to the patient's death!

During normal use, the AKTIVLINE® seat system should be checked regularly by the user for damage and the perfect working condition of all components.

**Note!**

The adjustment of settings is only permitted for qualified and trained personnel!

**Important!**

Only original parts may be used for service and maintenance work.

**13.1. Maintenance plan seating unit**

Checking / components	Daily	Weekly	Monthly
Visual inspection of the seating system, the attached components and accessories for visible damage	x		
Visual inspection of all fastening elements and screws	x		
Fastening/retightening of all fastening elements and screws		x	
Oil and grease moving parts			x

**13.2. Maintenance plan undercarriage**

Checking / components	Daily	Weekly	Monthly
Visual inspection of the undercarriage, the attached components and accessories for visible damage	x		
Check connection to seat system	x		
Check working condition of the brakes	x		
Fastening/retightening of all fastening elements and screws		x	
Check the chassis	x		
Oil and grease moving parts			x

## 14. Service plan

Service must be performed at least every twelve months (12 months) by the manufacturer or by a qualified authorized dealer in accordance with the according to the service schedules. The service is a chargeable service service, which is why a possible assumption of costs by the by the respective cost unit should be should be discussed in advance. The AKTIVLINE® Individual is a patient-specific patient-specific custom-made products (see 2.3) according to the individual

requirements of the prescription, which in principle allows use during the patient's during the growth phase of the patient. However regular checking of the settings to to the changes in the patient as he or she grows is mandatory, at least every twelve months (12 months). The review and any necessary adjustments are the responsibility of the manufacturer or qualified, authorized and trained authorized and trained personnel and and can be performed as part of regular service.

### 14.1. Service plan complete AKTIVLINE® system

#### Note!

The service should be carried out every twelve months! The manufacturer does not provide any guarantee or warranty, if the AKTIVLINE® Individual is not serviced regularly, at least once a year, by the manufacturer. service at the manufacturer's checked by the manufacturer.

Pos.	Test / Execution	OK	Adjustments necessary
1	Cleaning and disinfection	<input type="checkbox"/>	<input type="checkbox"/>
2	Checking of all settings of the AKTIVLINE® seat system in accordance with the CURRENT size of the patient	<input type="checkbox"/>	<input type="checkbox"/>
3	Check of setting of the head rest: <ul style="list-style-type: none"> <li>▶ Does the head rest lie against the patient's occiput (standard version)?</li> <li>▶ Does the bottom edge of the shell-shaped head rest lie against the base of the head (shell shape)?</li> <li>▶ Do the side elements exert pressure on the patient's head?</li> <li>▶ Does the patient have sufficient freedom of movement?</li> </ul>	<input type="checkbox"/>	<input type="checkbox"/>
4	Check the safety distances according to DIN EN 12182:2012 using the following table	<input type="checkbox"/>	<input type="checkbox"/>

#### Safety distances according to DIN EN 12182:2012 for children

To avoid	Between moving parts		Between fixed parts	
Trapping fingers	< 4 mm	or > 25 mm	< 5 mm	or > 12 mm
Trapping feet	< 25 mm	or > 120 mm	< 25 mm	or > 45 mm
Trapping the head	< 60 mm	or > 300 mm	< 60 mm	or > 250 mm
Trapping the genitals	< 8 mm	or > 75 mm	< 8 mm	or > 75 mm

#### Safety distances according to DIN EN 12182:2012 for adults

To avoid	Between moving parts		Between fixed parts	
Trapping fingers	< 8 mm	or > 25 mm	< 8 mm	or > 25 mm
Trapping feet	< 35 mm	or > 120 mm	< 35 mm	or > 100 mm
Trapping the head	< 120 mm	or > 300 mm	< 120 mm	or > 250 mm
Trapping the genitals	< 8 mm	or > 75 mm	< 8 mm	or > 75 mm

### 14.1.1. Service plan seating unit

Pos.	Area	Test			
		Setting/function	OK	Damage/ distortion	Exchange/ replacement
1	Performance test	Back recline	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Fixing unit	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Arm rest	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Head rest	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Abduction pommel	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Therapy tray table	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
2	Screw connections	Screws tightened	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Rust spots	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3	Appearance (visual inspection)	Paintwork	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Flock coating	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Cleanliness	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Fabric seams (position / appearance)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Rating plate	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### 14.1.2. Service plan undercarriage

Pos.	Area	Test				
		Setting/function	OK	Damage/ distortion	Exchange/ replacement	
1	Main frame	Main frame	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Anti-tipping system	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
2	Trapezoidal adapter		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
3	Front wheel	Wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Tire pressure	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Running performance of the wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Front wheel fork	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
4	Rear wheel	Mounting to frame	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Tire pressure	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Running performance of the wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Quick-release axles	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Drum brake	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
5	Brakes	Performance test, level 1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Performance test, level 2	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Performance test, level 3	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Bowden cable	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
6	Adjustability					
		Gas-pressure springs	Damages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Screw connections	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
		Seat adjustment	Performance test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Folding mechanism	Folding mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
			Locking mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Back angle adjustment	Performance test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Push bar	Height adjustment mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
			Locking mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	7	Screw connections	Size and length	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Pos.	Area	Test			
		Setting/function	OK	Damage/ distortion	Exchange/ replacement
8	Options	Completeness	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Damages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Loctite secured	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Covering caps available	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Wheel protector	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Spoke protector	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Arm rests	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		E-Fix	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		Centered foot rest	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		9	Appearance (visual inspection)	Paintwork	<input type="checkbox"/>
Cleanliness	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>
Rating plate	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>

**Important!**

Every change / replacement of relevant components is to be documented carefully.

## 15. Reuse/durability

The AKTIVLINE® seat shell system is not suitable for reuse. The supplied aid is individually manufactured to patient dimensions and requirements of a regulation and represents a custom-made device according to Regulation (EU) 2017/745 on medical devices, Article 2 No. 3 as well as Annex XIII. Please urgently observe the intended purpose listed in item 2.3.

The service life and thus the product life of the AKTIVLINE® Individual can be specified as min. 2 to max.

5 years under normal conditions. It should be noted that growth and development-related changes of the patient can lead to deviations, as this is an individual fitting / custom-made product. Regular inspection and performance of service and maintenance work can positively influence the service life of your medical device. The degree of disability and the duration of daily use also influence the service life. The product service life can therefore vary greatly from patient to patient in relation to the factors mentioned above.

## 16. Disposal/environmental compatibility

The issue of environmental protection is a high priority at interco Group.

We only process environment-friendly materials and specifically avoid manufacturing processes which are harmful to the environment. All packaging materials can be recycled. The AKTIVLINE® seat system can be

dismantled into separate material components to comply with optimal recycling and disposal. Please make inquiries with your local waste management company about the regulations concerning disposal and recycling. All components of the product are to be disposed of according to the respective country-specific environmental regulations.

## 17. EU declaration for custom made products

**Name des Herstellers / Manufacturer name**  
**Adresse des Herstellers / Manufacturer address**  
**SRN (Single Registration Number)**

interco Group GmbH  
Im Auel 50, Eitorf, 53783, Germany, [www.interco-reha.de](http://www.interco-reha.de)  
DE-MF-000005334

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## 18. Warranty

In addition to our General Terms and Conditions, we grant a two-year warranty for the AKTIVLINE® seat shell system product supplied by us. Unusable or defective parts will be repaired or replaced free of charge during the warranty period. The warranty shall exclude design-related modifications.

We accept no liability in the case of damage incurred as a result of combining our product with third-party products of any kind, which may pose significant risk in certain circumstances. Unless said products have been expressly approved by our company.

The warranty also excludes any defects occurring as the result of normal wear, excessive stress, violent and willful damage, and inappropriate use. Warranty obligations shall cease to apply if original interco replacement parts are not used for repairs and if repairs and modifications are not performed by interco. It is not permissible for any parties other than the manufacturer to make additions or modifica-

tions to your AKTIVLINE® seat system and any cases that should become known can be subject to legal action. Likewise, the warranty obligations of interco Group shall cease to apply in such cases and any liability claims directed at the manufacturer shall no longer be valid.


The continuous use of your AKTIVLINE® seat system may result in loose screws, nuts, and spokes. Retighten these at periodic intervals or arrange to have this done by your local dealer or by us.

### Note!

The manufacturer cannot grant warranty if the AKTIVLINE® seat system has not been checked regularly (at least once per year) as part of service.

The manufacturer cannot grant warranty if third party components (not manufactured by interco) are used without authorization by the manufacturer.

### Extended warranty

Upon submission in this form we will extend the warranty on your AKTIVLINE® to cover a total of 2 years. Please fill out this warranty form carefully and completely and return it to us. Thank you very much.

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© interco Group GmbH, Eitorf, 2022