CE

MINY®

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# Contents

---

1.	<b>Introduction</b>	28
1.1.	Product description	28
1.2.	Product overview	28
1.3.	Installation plan	28
1.4.	Environmental conditions	29
2.	<b>Medical requirements</b>	29
2.1.	Field of application/Target group	29
2.2.	Indications	29
2.3.	Contraindications	29
2.4.	Intended use	29
3.	<b>Product labeling</b>	29
4.	<b>Safety instructions</b>	29
5.	<b>Caution symbols</b>	31
6.	<b>Scope of delivery</b>	32
6.1.	Basic configuration	32
6.2.	Options/Accessories	32
6.2.1.	Lower extremities	32
6.2.3.	Other Accessories	33
7.	<b>Settings/Operation</b>	33
7.1.	Basic configuration	33
7.1.1.	Trapezoidal adaptor plate	33
7.1.2.	Attaching the seat shell	33
7.1.3.	Adjusting the inclination of the seat shell (tilt-in-space mechanism)	34
7.1.4.	Adjustment options for the push bar	34
7.1.5.	Folding mechanism for transporting the MINY®	35
7.1.6.	Drum brake-enabled back wheels	35
7.1.7.	Checking the tyre pressure	36
7.1.8.	Tilting mechanism and protection	36
7.2.	Options/Accessories	37
7.2.1.	Lower extremities	37
8.	<b>Compatibility/Combination options</b>	37
9.	<b>Technical data</b>	38
9.1.	Overall dimensions of MINY®	38
9.2.	Tyre pressure table	38
10.	<b>Transport/transport lock/climbing stairs</b>	39
10.1.	General information on transport	39
10.2.	Transport with an AMF force application point system	39
10.3.	Transport on aircraft	39
10.4.	Climbing stairs	39
11.	<b>Maintenance plan</b>	40
11.1.	Cleaning	40
11.1.1.	Materials used	40
11.1.2.	Cleaning of materials	40
11.1.3.	Disinfection of the product	41
11.2.	Maintenance	41
11.3.	Service	42
11.3.1.	Seating system	42
12.	<b>Reuse/durability</b>	43
13.	<b>Disposal/environmental compatibility</b>	43
14.	<b>Possible faults</b>	43
15.	<b>EU Declaration of Conformity</b>	44
16.	<b>Warranty policy</b>	45

## Information

Date of last update:  
February 2021

Revision status:  
MINY\_Gebrauchsanweisung\_interco\_de-en\_01

Please read these operating instructions carefully! Observe all safety instructions. If you have any questions, please contact your specialist dealer or interco.

## Note

As a rule, all interco products consist of a seat shell and/or an undercarriage. To simplify matters, only the interco product is mentioned in the following.

In these operating instructions we use the masculine form exclusively for the male/female users of an interco product and its male/female operators (parents, male/female caregivers, male/female therapists).

This way we wish to maintain the readability and comprehensibility of the texts. At the same time, we would like to expressly point out that this naturally includes the female form and other gender identities.

## Foreword

Dear Customer,


Thank you very much for choosing a quality product from interco. This is intended to contribute to a sustainable improvement in the everyday life and quality of life of the user. To ensure that you can always handle

the product safely and use it for a long time, we have described what you should bear in mind during use in these operating instructions.

It is important that you carefully read and understand all sections that apply to your interco product. Keep the operating instructions for future reference. Please follow all the safety instructions and symbols to prevent any type of injury. If you should have any questions about your product and its use, please ask your local qualified medical supplies dealer or contact us directly. We would also be pleased to inform you about further customization and adaptation options for your interco product.

Please observe the following in particular:

- ▶ Parents, accompanying persons, nurses and family members who operate this product should be instructed in the use of the product by authorized specialist personnel using these operating instructions.
- ▶ Your interco product was customized to the needs of the user. Subsequent changes to the product may only be made by trained specialist personnel. We recommend regular checks of the product and the product adaptation to ensure optimal long-term operation.
- ▶ Your interco product is only approved for combination with one of the products listed later on. interco will not accept any liability for combinations with medical devices and/or accessories from other manufacturers.
- ▶ Technical modifications may only be made with the written consent of interco.
- ▶ interco will not accept liability for modifications, changes to parts of the frame of the product and unapproved repairs to this product.
- ▶ interco reserves the right to technical changes to the model described in these operating instructions.

# 1. Introduction

## 1.1. Product description

MINY® is the street undercarriage for the interco seat unit range and offers a variety of adjustment possibilities. The inclination of the seat unit can be continuously adjusted via gas spring. Pneumatic tyres and drum brakes supplied as standard ensure a high degree of safety and comfort. At the front, MINY® features swivel wheels so that it is extremely manoeuvrable and particularly easy to push.

## 1.2. Product overview

MINY® can be folded up easily and quickly for transportation and it fits into any car boot.

The MINY® street undercarriage is suitable for all interco seat units featuring a trapezoidal adaptor. MINY® has been especially designed for outdoor use. Due to the pivoting front wheels the device is also easy and convenient to use indoors. For safety reasons, driving with grip tires should be avoided when knee brakes (optional) are missing.

- ① Front wheels
- ② Rear wheels
- ③ Lower frame
- ④ Upper frame
- ⑤ Trapezoidal seat unit adaptor
- ⑥ Back
- ⑦ Handle

## 1.3. Installation plan

### Authorization for delivery

MINY® may only be delivered to the end customer by authorized specialist dealers, medical device consultants or by the medical device consultants of interco Group.

### Scope of delivery

MINY®'s basic features are:

- ▶ aluminum base frame, slightly abducted version
- ▶ tilt-in-space mechanism using gas-pressure springs by approx. 35°
- ▶ trapezoidal adapter for interco seat shells
- ▶ adjustable seat height and seat depth
- ▶ depth- and angle-adjustable back
- ▶ back can be folded forward for transport
- ▶ different wheel sizes with pneumatic tires or solid tires
- ▶ steering wheels available in 200x50 or 7"
- ▶ Colour according to RAL chart
- ▶ Operating Guidelines

The mounted adapter receptacle, trapezoid-shaped is standard on all undercarriages by Interco GmbH. This makes it possible to mount any Interco seat shell with a trapezoidal adapter on various undercarriages by Interco. The adapter receptacle, trapezoid-shaped mounted on the undercarriage essentially consists of 3 functional groups:

- ▶ Mounting plate with the guide bushings compatible with the seat shell interface
- ▶ Spring-loaded catch mechanism for securing the seat shell to be mounted and keeping it in place
- ▶ Receptacle unit (centered) for mounting and attaching the center leg rest system with incrementally adjustable holder
- ▶ Side mounting units for assembling and attaching the ALULINE leg rest system with incrementally adjustable holder

### Setup/assembly

The product is delivered fully pre-assembled to the end customer, thus no components need to be assembled on-site. If options are installed at a later date, the assembly instructions are enclosed with the option.

The medical product advisor of the specialist dealer must ensure that the instructions are carried out correctly. (see below) The medical device consultant of interco Group does not have to be present during the installation of the product. However, he can be consulted for support in case of questions.

### Note!

Please make sure that the individual adjustment of the adaptor system has been carried out correctly. You can do so by checking that the adaptor rests flat on the front edge of the adaptor plate and the locking pin is fully engaged.

### Instruction of the customer

All functions of MINY® which are explained in the operating instructions must be demonstrated to the customer. The operating instructions remain with the customer. All adjustments to the product that directly affect the user may only be carried out by medical device consultants or in consultation with these.

### 1.4. Environmental conditions

The MINY® medical device must always be stored in a dry place to prevent possible corrosion. Please keep to ambient temperatures of -10 °C to +40 °C. Sand or other dirt particles should not be allowed to penetrate the product. Use in salty water is not permissible.

## 2. Medical requirements

### 2.1. Field of application/Target group

If there is an indication for a seat shell  
▶ To accommodate the existing seat shell

### 2.2. Indications

If there is an indication for a seat shell  
▶ To accommodate the existing seat shell

### 2.3. Contraindications

There are no known contraindications.

### 2.4. Intended use

When supplying children with seat units and modular seating systems, undercarriages or chassis' for indoor and/or outdoor use are compulsory. They increase the mobility and make the transport of the insured person possible. Therefore, the purpose of all chassis' – in particular MINY® – always is the reception of a seat unit.

## 3. Product labeling

- A Type designation
- B Country of manufacture
- C Manufacturer's data
- D Serial number
- E Observe warnings/safety instructions
- F Read operating instructions prior to use
- G Seat width
- H Maximum weight

Please always state the serial number and size of the MINY® when ordering spare parts or accessories!

## 4. Safety instructions

Please note that you are responsible for the safety of the user. This could be at risk if you do not follow the safety guidelines. Nevertheless, it is not possible to cover all potential circumstances and unforeseeable situations in these guidelines. Therefore common sense, caution and prudence are required when handling the product. Every person operating the product should know and understand all the instructions. Furthermore, you should explain the operation of the product and its accessories to every person using the product.

Please contact your local specialist dealer or interco directly if the instructions are unclear and a more detailed explanation is required or if you have any other questions.

### Meaning of the safety instructions

#### Warning!

Warning against possible severe damage due to accidents or injuries.

#### Caution!

Warning against possible damage due to accidents or injuries.

#### Note!

Warning for possible technical damage.

### Product-specific safety instructions

#### Warning!

*Risk of injury from unattended parking*

- If the product is left unattended with user, the user or other persons may be injured.
- ▶ Please make absolutely sure that the person sitting in the product is properly secured at all times by a suitable restraint system (harness).
  - ▶ Never leave the user sitting in the product unattended, not even when the user has been

securely strapped in and the brakes have been set.

- ▶ Never leave the product alone with the user without expert supervision to exclude the risk of injury to the user or other persons.

### Warning!

*Risk of suffocation due to the packaging material*

Incorrect handling of packaging material may cause a risk of suffocation.

- ▶ Keep the packaging materials out of reach of children.

### Warning!

*Risks when climbing stairs*

When using stairs, there may be a risk of injury to the user of the product or the product may be damaged. Therefore only lift the product without the user sitting in it.

- ▶ Stairs are only permitted to be negotiated with the assistance of accompanying persons.
- ▶ If facilities such as access ramps or lifts are available for this purpose, these must be used. In the absence of such facilities, two accompanying persons can carry the product over obstacles without the user.
- ▶ Techniques for climbing stairs are explained in detail in Section 10; it is imperative that you comply with these.
- ▶ Make sure that the accompanying persons only grip the product on firmly mounted or bonded components.
- ▶ Do not lift the product by the foot rest, head rest, arm rests or similar parts.

### Warning!

*Risks of self-initiated modification of the settings*

If settings are changed autonomously, the function of the product and the safety of the user may be affected.

- ▶ Adhere to the settings made by the specialist dealer. You may only adjust settings by yourself as described in Section 7.
- ▶ If you experience problems with settings, please contact your specialist consultant or interco directly.
- ▶ Discuss all changes to the settings properly with the specialist personnel or therapist.
- ▶ It is not permissible for any parties other than interco to make additions or modifications to the product and any cases that should become known can be subject to legal action. In addition, the warranty obligations of interco Group shall cease to apply in such cases and any liability claims directed at the manufacturer shall no longer be valid.

### Caution!

*Risk of injury due to operating errors or lacking knowledge of the operating instructions*

If the instructions in the operating instructions are disregarded, application errors may occur and the safety of the user may be impaired. Please read these operating instructions carefully first.

- ▶ Familiarize yourself with the functions of the product and its handling prior to initial use. Practice handling the product before use.
- ▶ Observe the maximum load capacity of the product (see Section 9: Technical Data).
- ▶ Do not use the foot rest when getting into or out of the product.

### Caution!

*Risk of tipping due to incorrect center of gravity / shift in center of gravity*

If the product's center of gravity is shifted significantly, for example on slopes and inclines or if the user leans out of the product too far, the product may tip or roll over.

- ▶ Test the effects of shifts in the center of gravity on the behavior of the product, for example, on declines, inclines, and lateral slopes or while negotiating obstacles only with the reliable support of an assistant.
- ▶ Always stand the product on horizontal, flat surfaces only. If it is unavoidable to park it on an incline, put the seat in an upright position. There is a risk of the product tipping backwards on inclines in the lying position.
- ▶ When grasping objects (lying in front of, to the side or behind the device) make sure that the user does not lean out of the seat too far.
- ▶ Do not hang heavy shopping bags or similar on the push bar or the levers.

### Caution!

*Risks in public road traffic*

Dangerous situations can occur when driving on public roads.

- ▶ Please note that you must observe traffic regulations when traveling on public streets.
- ▶ In dark conditions wear bright clothing wherever possible in order to be more visible. We also recommend installing an active light source.

### Caution!

*Risk of injury to hands and fingers*

Your hands or fingers can get caught during adjustment and assembly work.

- ▶ Always be aware of your fingers during operation and all adjustment and assembly work.

### Caution!

*Dangers due to loosened components during use*

Screws and nuts on the product may become loose when the product is used. This can result in risks of accidents and injuries for the user and the operators.

- ▶ Check all screws and nuts regularly and retighten them.
- ▶ Alternatively, you can also instruct your specialist dealer or interco to do this.

### Caution!

*Dangers when making adjustments to the product and connections*

If you do not carry out adjustments to the product properly, the stability and safety of the product may be affected.

- ▶ Ensure that the screws are tightened again after all adjustment and assembly work.
- ▶ Always make sure that the seating system is securely connected to the undercarriage.
- ▶ If you have a question regarding the connection of the seating system and the undercarriage, always contact your specialist dealer or interco directly.
- ▶ Carry out a regular visual inspection of all connections of the product. Consult your local specialist dealer in case of doubt.

### Caution!

*General risks due to careless handling of the product*

The product is only approved for its intended use. Improper or careless use can result in safety and injury risks for the user and other persons.

- ▶ Avoid actions, such as driving against an obstacle (stair, curb) without using the brake or down from elevated surfaces.
- ▶ Use caution when handling fire (e.g. a burning cigarette)! Upholstery and coverings could catch fire.
- ▶ Make sure to keep children from playing with the therapy chair or any parts thereof. There is a risk of injury due to trapping and crushing hands and fingers.
- ▶ Children must not play next to, underneath or in front of the seating system when the seating system is in motion and/or the seat or back unit is in use.

## 5. Caution symbols

**Note:** For adjusting the angle, open control lever at least 3 turns!

**Place of mounting:** push handle on the left

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**Note:** Assure appropriate tilting protection when adjusting seat angle!

**Place of mounting:** push handle on the right

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**Note:** All screw connections, that were loosened by changing the settings, must be re-tightened!

**Place of mounting:** Both upper frames outside

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**Note:** This wheelchair is completed with a special wheelchair restraining system. It is allowed to use this wheelchair as seat in vehicles, if the security belts are used according to the specifications.

**Place of mounting:** Lower frame outside

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**Note:** It is not permitted to use this wheelchair as a car seat.

**Place of mounting:** Lower frame outside

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**Note:** Do not use on stairs!

**Place of mounting:** Push handle on the left

**Note:** Do not carry!

**Place of mounting:** Push handle on the right

## 6. Scope of delivery

### 6.1. Basic configuration

MINY® is a seat shell undercarriage with

- ▶ Aluminum base frame, slightly abducted version
- ▶ Tilt-in-space mechanism using gas-pressure springs by approx. 35°
- ▶ Trapezoidal adapter for interco seat shells
- ▶ Adjustable seat height and seat depth
- ▶ Depth- and angle-adjustable back
- ▶ Back can be folded forward for transport
- ▶ Tilting mechanism and protection
- ▶ Different wheel sizes with pneumatic tires or solid tires
- ▶ Steering wheels available in 200x50 or 7"
- ▶ Colour according to RAL chart
- ▶ Operating guidelines

### 6.2. Options/Accessories

#### 6.2.1. Lower extremities

Article	Size	Art. no.
ALULINE leg rest system with incrementally adjustable holder for leg rest system Foot rest, complete version, foot rest powder coated in black, includes corner guards and foot rest overlay for all interco undercarriages	size 0-3	46013
ALULINE leg rest system with incrementally adjustable holder for leg rest system Foot rest with individual foot supports, foot rest powder coated in black, includes corner guards and foot rest overlay for all interco undercarriages	size 0-3	46016
Centered leg rest	size 0-3	77050-0 77050-1 77050-2 77050-3
Calf support, complete version		46030

Article	Size	Art. no.
Black synthetic leather cover for calf support, complete version		46031
Calf support, single version		46120
Black synthetic leather cover for calf support, single version		46121
Plastic foot rest riser, black		46530
Heel stop for leg rest system, complete version		46551
Heel stop for leg rest with individual foot supports		46552
Heel stop cushion without cover		46510
Black synthetic leather cover for heel stop cushion		46511
Foot rest cushion without cover		46520
Black synthetic leather cover for foot rest cushion		46521
Foot straps with hook and loop fastener with padded instep		46410
Foot straps made of leather with padded instep		46411
Foot straps with deflection roller and padded instep with hook and loop fastener		46412
Foot straps with incremental adjustment and padded instep		46419
Foot straps with clip buckle and padded instep		46420
Foot shells made of aluminum	size 1-3	46210
Foot gaiters without forefoot restraint made of neoprene	size 1-3	46414
Foot gaiters with forefoot restraint made of neoprene	size 1-3	46416

### 6.2.3. Other Accessories

Article	Size	Art. no.
Trapezoidal adapter to be mounted underneath the seat shell incl. mounting material		40110

## 7. Settings/Operation

Every interco product is customized to the individual needs of the user and configured accordingly. To do justice to all possible product designs, the list of options and accessories as well as their descriptions is therefore somewhat more extensive.

Therefore please bear in mind that this Section may also describe components that possibly do not exist on your product.

### 7.1. Basic configuration

#### 7.1.1. Trapezoidal adaptor plate

The mounted adapter receptacle, trapezoid-shaped is standard on all undercarriages by Interco GmbH. This makes it possible to mount any Interco seat shell with a trapezoidal adapter on various undercarriages by Interco. The adapter receptacle, trapezoid-shaped mounted on the undercarriage essentially consists of 3 functional groups (Fig. 1):

- ▶ Mounting plate with the guide bushings compatible with the seat shell interface (Fig. 1/1)
- ▶ Spring-loaded catch mechanism for securing the seat shell to be mounted and keeping it in place (Fig. 1/2)
- ▶ Receptacle unit (centered) for mounting and attaching the center leg rest system with incrementally adjustable holder (Fig. 1/3)
- ▶ Side mounting units for assembling and attaching the ALULINE leg rest system with incrementally adjustable holder (Fig. 1/4)


Fig. 1

#### Important!

Please make sure that the interface system has been adjusted correctly in each individual case. This can be determined by ensuring that the seat shell interface is flush against the front edge on the adapter receptacle, trapezoid-shaped, and the spring-loaded catch mechanism snaps completely into place!

#### 7.1.2. Attaching the seat shell

The seat shell is attached to the adapter receptacle, trapezoid-shaped (Fig. 1) on the undercarriage with the trapezoidal adapter and secured in place with the two spring-loaded catch mechanisms.

- ▶ First, apply the parking brakes to prevent MINY® from accidentally rolling away.

- ▶ Then extend the spring-loaded catch mechanisms (Fig. 2) and turn them 90 degrees.
- ▶ Set the back end of the seat shell interface on the front side of the adapter receptacle, trapezoid-shaped; then slide the seat shell until it stops at the back and secure it with the spring-loaded catch mechanisms (rotate through 90 degrees so you hear it snap into place).
- ▶ Finally, check the seat shell and the seat shell locking system to ensure they are secure.

#### Removing the seat shell

First, extend the spring-loaded catch mechanisms and turn them 90 degrees. Then slide the seat shell forward and remove it.

#### 7.1.3. Adjusting the inclination of the seat shell (tilt-in-space mechanism)

##### Caution!

Due to inclination adjustment of the seat shell the stability of the MINY® is only guaranteed if the tilt stop system is extended.

The MINY® is equipped with a gas-pressure spring for the inclination of the seat shell. In order to adjust the inclination of the seat shell, proceed as follows:

- ▶ Activate the parking brakes
- ▶ Pull the gas-pressure spring release lever (Fig. 3) mounted to the right of the push bar frame.
- ▶ Adjust the inclination of the seat shell by pulling the push bar up or pressing it down while pulling on the gas-pressure spring release lever.
- ▶ Let go of the release lever to lock the seat shell into the adjusted inclination angle.
- ▶ Deactivate the parking brake, as appropriate

##### Important!

Never pull the release lever without a firm grip on the push bar if a patient is seated in the wheelchair.

#### 7.1.4. Adjustment options for the push bar

The push bar has various adjustment options. This makes it possible to achieve optimal customization to the person being pushed. The following adjustment options are possible:

##### *Adjusting the angle of the top piece of the push bar (with incremental adjustment only)*

In order to adjust the angle of the push bar (Fig. 4), proceed as follows:

- ▶ Activate the parking brakes
- ▶ Loosen the screw connections on the right and left ratchet joints on the push bar by turning the clamp levers counter-clockwise until the top piece of the push bar becomes loose.
- ▶ Move the push bar into the desired angle.
- ▶ Hold the push bar in the desired position with one hand and secure it in place by tightening the corresponding clamp lever with the other hand.
- ▶ Secure the screw connection on the other side as well by tightening the clamp lever.
- ▶ When tightening the ratchet joints make sure that the teeth are engaging properly.
- ▶ Deactivate the parking brake, as appropriate


Fig. 2

Fig. 3


Fig. 4

Fig. 5

##### *Adjusting the height of the push bar unit*

The height of the push bar always needs to be adjusted when it is necessary to customize the MINY® to the person being pushed or if the inclination of the seat shell has been changed.

In order to adjust the height of the push bar, proceed as follows:

- ▶ First, loosen the quick-release fasteners (Fig. 5) on the fork tubes by opening the two quick-release levers outwards.
- ▶ Then adjust the push bar to the desired height. To do this, pull or press the push bar to the desired height.

- ▶ Secure the push bar in position by tightening the quick-release connectors. To do this, tilt the lever on the quick-release fasteners back in.
- ▶ If the clamping pressure is not sufficient (push bar not secured in place), then turn the screw across from the quick-release fastener clamp lever (with lever in open position) counter-clockwise, until sufficient clamping pressure is achieved by closing the clamp lever. Repeat this step, as necessary.

*Adjusting the inclination of the back unit with push bar*  
(Here, the back unit, including push bar, is folded back.)

In addition to the adjustment options described above, it is also possible to swivel the push bar backwards by up to 30° (in 5° increments). This can become necessary, for example, if a seat shell is being used with a seat angle of more than 90°. This enables the gap between the back of the seat shell and the push bar to be as uniformly small as possible. This adjustment option is not intended for individual customization to the person being pushed.

In order to adjust the inclination of the push bar, proceed as follows:

- ▶ Loosen the spring-loaded locking system for the folding mechanism on the push bar by pressing down on the central cable system (Fig. 6) connecting the two perforated plates to the bottom of the two back posts and simultaneously swiveling the push bar backwards.
- ▶ Then it is possible to shift the back unit including the push bar backwards (up to 30 degrees max.).

#### 7.1.5. Folding mechanism for transporting the MINY®

The MINY® is equipped with a folding mechanism in order to make transport easier. The back unit including push bar is folded forward; it is also possible to remove the back wheels. This greatly condenses the size of the wheelchair for transport.

##### *Folding up the MINY®*

In order to fold up the MINY®, proceed as follows:

- ▶ Remove the seat shell (see Section 7.1.2.)
- ▶ Loosen the height adjustment mechanism on the push bar (see section 7.1.4.)
- ▶ Loosen the screw connections on the push bar angle adjustment mechanisms on the right and left side of the push bar (see section 7.1.4.)
- ▶ Loosen the spring-loaded locking system for the folding mechanism on the push bar by pressing down on the central cable system (Fig. 6) connecting the two perforated plates to the bottom of the two back posts and simultaneously swiveling the push bar forward until the back together with the push bar comes into contact with the top frame of the MINY®.
- ▶ Then remove the back wheels. To do this, raise the MINY® slightly so that the back wheel to be removed is no longer in contact with the ground.
- ▶ Grasp the spokes on the back wheel and use your thumb to press the button located at the center of the wheel axle.
- ▶ Keep the button pressed down and pull the wheel out of the axle receptacle. Make sure that the brake has been deactivated on the wheel being removed!
- ▶ Repeat the same steps with the second back wheel.
- ▶ The MINY® is now ready for transport.

Fig. 6

Fig. 7

##### *Unfolding the MINY®*

In order to unfold the MINY®, proceed as follows:

- ▶ First, check to ensure that the MINY® parking brake is deactivated.
- ▶ Pick the MINY® up by the back end.
- ▶ Grasp the spokes on one of the back wheels and use your thumb to press the button located on the wheel axle.
- ▶ Keep the button pressed down and slide the wheel back into the axle receptacle on the frame.
- ▶ Release the button and check to ensure that the wheel is secure. Make sure that the wheel can no longer be pulled off!
- ▶ Attention: Do not press the button while performing this check.
- ▶ Repeat the same steps with the other back wheel.
- ▶ Fold the push bar unit backward until it audibly clicks into place.
- ▶ Then secure the pivoting top piece of the push bar into position, as described above (see section 7.1.4.).
- ▶ Then adjust the push bar to the desired height. To do this, pull the push bar to the desired height.
- ▶ Secure the push bar in position by tightening the quick-release connectors, as described above (see section 7.1.4.).
- ▶ If the clamping pressure is not sufficient (push bar still loose), then proceed as described above (see section 7.1.4.).
- ▶ Mount the seat shell

#### 7.1.6. Drum brake-enabled back wheels

The MINY® is equipped with two parking drum brakes on the two back wheels, which can be activated independently of each other. The brake levers are equipped with a three-level ratchet mechanism, which provides for varying levels of braking action.

##### *Activating the parking brake*

In order to activate the parking brake, proceed as follows:

- ▶ Engage both brake handles on the push bar (Fig. 7).

Make sure that the brake handle does not return to the starting position (locking mechanism will click audibly into place); otherwise the parking brake has not been activated.

- ▶ Check to make sure that the brake is indeed activated by attempting to push the MINY®.

**Caution!**

If the braking action is not sufficient, then have the brakes readjusted by your medical supplies dealer or readjust them yourself, as described below in this section.

You can also use the parking brakes for better control of the wheelchair on descents. In such cases, engage the parking brake, as needed, by pressing it into the first or second ratchet position only.

*Deactivating the parking brake*

In order to deactivate the parking brake, proceed as follows:

- ▶ Deactivate the locking mechanism by activating the ratchet lever and allowing the brake lever to glide back into the starting position.
- ▶ You may need to release the ratchet lever further by pulling slightly on the brake lever before the brake can be deactivated.

*Adjusting the parking brake*

**Tool:** Size 8 spanner

**Prerequisite:** Insufficient braking power of the parking brake; MINY® can be pushed even though the parking brake is applied.

To adjust the parking brake, proceed as follows:

- ▶ Make sure that the parking brake is not applied; if necessary, release the parking brake (see section 7.1.6.)
- ▶ The Bowden cables of the parking brakes are routed through an adjusting bolt before the brake backing plates, which allows an adjustment of the braking power.
- ▶ By turning this bolt anti-clockwise, the brake cable is brought under greater tension, pre-tensioning the brake. The adjusting bolt is secured by a nut
- ▶ Before adjusting the brake, release the securing nut by turning it to the left (depending on the version, a size 8 spanner may be required)
- ▶ Next, turn the adjusting bolt anti-clockwise until it becomes more difficult to turn the wheel. To check, lift the MINY® on one side and turn the respective wheel with your other hand
- ▶ Then, turn the adjusting bolt clockwise until the wheel which has to be adjusted runs freely again
- ▶ Turn the counter nut located on the thread of the adjusting nut clockwise as far as it will go
- ▶ If necessary, repeat the procedure for the other wheel.

**Important!**

In order to check the adjustment of the brake, pull the brake lever several times before finally fastening the counter nut. If necessary, you have to repeat this operation, as some experience is required to carry out the adjustment.

Fig. 8

**7.1.7. Checking the tyre pressure**

MINY® features pneumatic tyres as standard. In order to keep rolling resistance to an absolute minimum, you should check the tyre pressure at regular intervals. The tyres are standardly equipped with automobile valves, so you can check the pressure by using the tyre gauge at your local petrol station. The required filling pressure is printed on the tyres and is also specified in Section 9 ("Technical data") of these instructions.

To check the tyre pressure, please proceed as follows:

- ▶ Apply the parking brake (see section 7.1.6.)
- ▶ Unscrew the valve cap from the valve by turning it anti-clockwise
- ▶ Press the tyre gauge nozzle firmly onto the valve
- ▶ Correct the tyre pressure to the specified value
- ▶ Remove the tyre gauge nozzle from the valve
- ▶ Screw the valve cap onto the thread of the valve (clockwise)
- ▶ Repeat the operation until you have checked all tyres
- ▶ If necessary, release the parking brake (see section 7.1.6.)

**Caution!**

Depending on the version, the front and rear wheels require different tyre pressures.

**7.1.8. Tilting mechanism and protection**

Tilting mechanism 2 in 1 (tilting protection and support). Sliding and lockable with locking pins.

- ▶ Pull out the spring-loaded catch mechanism and push in the tilt stop system in order to use it as a tilt assist device (Fig. 8).
- ▶ Make sure you hear the spring-loaded catch mechanism snap into place and check to make sure it is secure.

**Important!**

In this position there is no protection against tilting!

To prevent tilting, the tilt stop system must have snapped into place in the extended position.

**Important!**

The tilt assist function is restricted because it can move in again!

*Use as a tilt assist device*

Pull back the push bars slightly and at the same time step on the tilt assist device. As a result you can tilt the undercarriage onto the rear wheels without any great effort. In this way you can easily overcome obstacles when traveling forward or backward.

**Important!**

The user must be physically capable of holding the undercarriage together with the seat shell and occupant securely!

**Important!**

Always pay attention that the tilting protection wheels together with the two small wheels point towards the ground when a patient is sitting in the system.

**7.2. Options/Accessories**

**7.2.1. Lower extremities**

*Centered leg rest system*

The centered leg rest is factory-mounted with your order. Subsequent changes are possible for individual customization.

**Adjustment options**

The centered leg rest system has various adjustment options for optimal customization:

**a) Vertical inclination adjustment range**

In order to adjust the centered leg rest, proceed as follows:

- ▶ Activate the parking brakes
- ▶ Loosen the screw connections on the centered leg rest by turning the clamp levers counter-clockwise until the top piece becomes loose.
- ▶ Move the centered foot rest into the desired angle.
- ▶ Hold the centered foot rest in the desired angle with one hand and secure it in place by tightening the corresponding clamp lever (Fig. 14) with the other hand.
- ▶ Deactivate the parking brake, as appropriate

**b) Foldable back**

The foot plate is folded up and back by hand (Fig. 15).

**c) Height and depth adjustment**

Disengage the clamping lever (Fig. 16) on the clamping block of the foot rest and set the appropriate length of lower leg. To adjust the depth (thigh length) undo the screw connection (Fig. 17) below the seat adapter. Tighten all the screws again when adjustment has been completed.

**8. Compatibility/Combination options**

Undercarriage	MINY® SB40 with 22" / 24" wheels	MINY® SB40 with 12" / 16" wheels	MINY® SB35 with 22" / 24" wheels	MINY® SB35 with 12" / 16" wheels	MINY® SB30 with 22" / 24" wheels	MINY® SB30 with 12" / 16" wheels
ORTHOLIGHT® size 1	✓	✓	✓	✓	✓	✓
ORTHOLIGHT® size 2	✓	✓	✓	✓	✓	✓
ORTHOLIGHT® size 3	✓	✓	✓	✓	-	-
ORTHOLIGHT® individual max. SB 32 cm	✓	✓	-	✓	-	-
ORTHOLIGHT® size 4	✓	✓	-	-	-	-
ORTHOLIGHT® size 5	✓ without wheel cover	✓ max. 60 kg	-	-	-	-
ORTHOLIGHT® size 6	-	-	-	-	-	-
ORTHOLIGHT® individual from SB 32 cm up to max. 41 cm	-	-	-	-	-	-
MAGICLIGHT® size 1	✓	✓	✓	✓	✓ without wheel cover	✓
MAGICLIGHT® size 2	✓	✓	✓	✓	✓ without wheel cover	✓
MAGICLIGHT® size 3	✓	✓	✓	✓	✓	✓
MAGICLINE® size 1	-	-	✓	✓	✓	✓
MAGICLINE® size 2	✓	✓	✓	✓	-	✓
MAGICLINE® size 3	✓	✓	✓	✓	-	-
ERGOLINE® according to patient's measures	✓ up to seat width 45 cm without wheel cover	✓ up to seat width 45 cm	✓ up to seat width 40 cm	✓ up to seat width 40 cm	✓ up to seat width 35 cm without wheel cover	✓ up to seat width 35 cm

## 9. Technical data

### Caution!

The maximum load of MINY® with attached seat unit amounts to 60 kg.

### 9.1 Overall dimensions of MINY®

	SB 35	SB 35	SB 40
Weight, including wheels in the basic configuration (without leg rest)	approx. 16,5 kg	approx. 17,0 kg	approx. 18,0 kg
Frame width without attachments (measured in abduction)	44 cm	49 cm	54 cm
Seat width	30 cm	35 cm	40 cm
Seat height adapter min. - max.	23-30 cm	23-30 cm	23-30 cm
Seat height from foot rest min.-max.	25-39 cm	15-39 cm	7-39 cm
Overall width (with wheel camber)	50 cm (+ 5 cm)	55 cm (+ 5 cm)	60 cm (+ 5 cm)
Max. length without push bar (24" wheels)	75 cm (78 cm)	75 cm (78 cm)	75 cm (78 cm)
Height / length folded	52/79 cm	52/79 cm	52/79 cm
Height to top adapter plate min. - max.	46 - 48 cm	46 - 48 cm	46 - 48 cm
Height of push bar min. - max.	98 - 110 cm	98 - 110 cm	98 - 110 cm
Seat shell inclination angle	-5° to 30°	-5° to 30°	-5° to 30°
Leg rest holder to seat area (adapter plate) - incremental adjustment	0° to 90°	0° to 90°	0° to 90°
Gripping wheel diameter	48 cm (22"-wheels) · 53 cm (24"-wheels)		
Turning radius	approx. 92 cm (22"-wheels)		
Max. load	60 kg	60 kg	60 kg
Wheel camber	6°	6°	6°

### 9.2 Tyre pressure table

	Tyre pressure (bar) maximum
Front swivel wheels	2.5
Rear wheels 12"	2.5
Rear wheels 16"	2.5
Rear wheels 22"	4.0
Rear wheels 24"	4.0

## 10. Transport/transport lock/climbing stairs

### 10.1. General information on transport

Only the tested and approved wheelchair and occupant restraint system from AMF-Bruns may be used for transporting the MINY® undercarriage in motor vehicles for the transportation of the disabled. In addition, a 4-point restraint system from AMF-Bruns is required for transport in motor vehicles. For information on how to correctly secure the force application point system in the respective vehicle, please see the manufacturer's instructions and the following chapter 10.2. „Transport with an AMF force application point system“.

#### Warning!

*Dangers due to incorrect fastening of the transport safety device*

If you attach the safety device for transport in a vehicle to the wrong places on the product, it cannot be transported safely.

- ▶ Never secure the transport safety device to the frame or other unsuitable basic components!

### 10.2. Transport with an AMF force application point system

The MINY® undercarriage was subjected to a crash test as standard version with maximum fittings in accordance with ISO 7176-19. Some models are custom-made and cannot be tested in a crash test. Nonetheless, the custom-made model is also fitted with a force application point system manufactured by AMF Bruns. Further information is available at [www.amf-bruns.de](http://www.amf-bruns.de).

The undercarriage is marked on four points which indicate the mounting points for the tension belts for transportation with motor vehicles. Only use the provided mounting points for strapping with tension belts. When fitted correctly, the head rest of the seat shell is very robust, but cannot replace head rests fitted to motor vehicles.

The MINY® undercarriage has been tested without an additional drive. If an additional electrical drive is to be added, then this accessory is to be disassembled if the MINY® undercarriage is to be transported in a motor vehicle.

### 10.3. Transport on aircraft

When traveling by air, the product may only be transported without the user in the cargo hold of the aircraft. There are no further instructions or regulations for transport in an aircraft. We therefore strongly recommend that you contact the airline before traveling by air to clarify the transport options in each individual case.

#### Caution!

*Transport on aircraft*

The user must never be in the device when it is transported in the cargo hold of an aircraft!

### 10.4. Climbing stairs

Stairs or high curbs can only be negotiated with the assistance of two accompanying persons. Assistants are only permitted to grip the front frame tube and the back brace of the frame, without lifting or carrying the product when doing so.

When ascending, the assistant behind the device pulls the device up the step. The assistant in the front stabilizes the position by pressing the device equally on both sides into the step.

When descending, the assistant at the front brakes by pressing the device evenly into the steps. The assistant behind holds the back brace, and stabilizes and maintains the device in the correct position. To prevent damage, slowly roll the product down the steps. Then properly readjust the anti-tipping system (if installed).

#### Caution!

*Risks when climbing stairs*

If you hold or lift the product by its attachments when climbing stairs, these may tear off and the user may be injured.

- ▶ Do not hold the product by the push bar when negotiating stairs.
- ▶ Never lift the product by the foot rest, the wheels or the arm rests or other attached parts.
- ▶ Do not lift the product with the user sitting in it.

## 11. Maintenance plan

---

### 11.1. Cleaning

Your interco product is made of high-quality and durable materials which are easy to clean. Please observe the following instructions for long durability of the product:

- ▶ Always keep the product dry, at temperatures between -10 and +40 °C. This prevents corrosion.
- ▶ Sand or other dirt particles should not be allowed to penetrate the product. This could damage the mechanics or - if applicable - the wheels.
- ▶ Never use the product in salt water.
- ▶ The cushions can be removed for cleaning and care.

#### 11.1.1. Materials used

interco uses the following materials in the manufacturing of medical devices and accessories:

- ▶ Metals (aluminum, steel, stainless steel)
- ▶ Plastics (POM, PA6, PTEG)
- ▶ Beech wood (for therapy tray tables)
- ▶ Foams (frame foam, polyurethane foam, hard foam, azote foam, neoprene)
- ▶ Upholstery fabrics (mesh material, Trevira CS fabrics, artificial leather, Dartex)
- ▶ Powder coatings
- ▶ DIN and standard parts

All materials used are classified as harmless to health.

#### Powder-coated aluminum

interco undercarriages and some add-on parts are made of aluminum and feature a resistant plastic coating. Powder varnishes are used for coating. Powder varnishes are solvent-free coating materials, composed of different resin systems, pigments as well as special additives.

#### Fabrics

interco uses different fabrics for covers and fixation aids. All materials used are flame-retardant.

#### Synthetic leather

interco uses easy-care designer synthetic leather for Dynaline restraints. This synthetic leather is long-lasting, hard-wearing, very elastic and tear-resistant. It provides a good grip and offers a comfortable seating experience. In addition, it is resistant to urine, blood and disinfectants, free of azo compounds, CFC, formaldehyde, as well as free of PCTs, PCPs and PCBs.

#### Spacer fabrics/mesh material

The spacer fabrics used have two outer textile layers, which are connected by a number of spacer threads. This 3D structure creates an insulating layer of air between the textile layers, which balances temperature fluctuations effectively. Furthermore, the spacer fabrics are very breathable, pressure-elastic, light and tear-resistant.

#### Supra-tex

The SUPRA-tex upholstery fabric consists of highly tear-resistant polyamide. It is impregnated to be water-repellent, PU-coated on the back and flame-retardant.

#### Neoprene

The neoprene used is a foam based on chloroprene rubber (CR). It insulates very well, is waterproof and very well tolerated by the skin. In addition, the neoprene is tear and abrasion resistant, very elastic, resistant to fungi and bacteria and eco-friendly.

#### Trevira CS (fabrics)

Among others, interco uses Trevira CS fabrics for seat covers. These are easy to clean, dry quickly and keep their shape. In addition, Trevira CS can be disinfected under hospital conditions.

#### Dartex incontinence cover fabric

interco uses the incontinence-resistant upholstery fabric Dartex for incontinence-resistant seat shells. The high-quality fabric meets all quality and safety requirements and is easy to wash off, sufficiently lightfast and flame-retardant. It does not cause allergies, is free of formaldehyde and azo dyes.

### 11.1.2. Cleaning of materials

#### Powder-coated aluminum

Clean undercarriages and attachments made of powder-coated aluminum with mild soapy water and a soft sponge. Do not use aggressive or abrasive agents. These can damage the surface of the coating.

#### Fabrics

Clean the fabric covers regularly before stains set in. You can vacuum or carefully brush the textiles. Spilled liquids should be wiped up immediately with an absorbent cloth. You can also put a mild detergent solution on a cloth, dab the stain with it and then rub off gently with clear water.

#### Note!

Bleach or whiteners should never be applied to the fabrics / covers. Never dry the covers in the dryer! Only iron at a low temperature.

If you wish to wash the cover of the seating system more often, a second cover for changing is recommended. Ask your local specialized dealers on this topic.

#### Synthetic leather

Clean synthetic leather covers with a cloth and commercially available mild soap solution.

#### Spacer fabrics/mesh material

Wash removable covers made of spacer fabrics/mesh material at a maximum of 30 °C in the washing machine with the delicates program.

#### Supra-tex

interco uses Supra-tex for the arm rests. Clean these regularly with a cloth and commercially available mild soapy water.

### Neoprene

Wash removable fixation aids made of neoprene at a maximum of 30 °C in the washing machine. Use a mild liquid detergent or a commercially available neoprene detergent. Use a gentle wash cycle and do not spin-dry the neoprene covers. Then hang up the covers wet to dry.

### Trevira CS (fabrics)

Removable covers made of Trevira CS can be dry cleaned. Washing under easy-care conditions at 60°C is possible. Do not dry and do not iron!

### Dartex

Test washed by the manufacturer at temperatures between 40 and 95 °C.

### Accessory components made of metal, plastic, wood

Please clean accessory components with a conventional mild household soap solution.

#### 11.1.3. Disinfection of the product

Please use a water-based agent for disinfection of the product, e.g. Terralin, Quartamon Med or Sagrotan Original concentrate. Please observe the manufacturer's instructions for use when doing so. When disinfecting the fabrics, also follow the instructions of the fabric manufacturer regarding the use of brand name disinfectants.

Rinse off all detergents and disinfectants thoroughly and leave the cleaned object to dry completely before use. If this is neglected, this can lead to chemical reactions making the material brittle or lead to skin irritations of the patient.

#### Important!

*No disinfectants for synthetic leather/skai*  
Synthetic leather covers/skai covers may not be cleaned with disinfecting agents, as the material gets porous over time and cracks.

#### Note!

##### Wipe disinfection

If wipe disinfection is necessary, use a disinfectant that is suitable for surface disinfection according to the DGH disinfectant list.

#### 11.2. Maintenance

Your interco product is checked by trained and qualified employees for proper functioning and freedom from defects before being handed over to the user.

To ensure safe use of the product over its entire service life (according to the manufacturer from at least 2 years to 8 years maximum) and to ensure a perfect technical condition, it must be maintained and inspected regularly according to the following maintenance (every 6 months) and service (every 12 months) plans.

#### Note!

##### Changing settings

The changing of settings (other than those described in Section 7) is only permitted for qualified and trained personnel!

Incorrect settings by the user, operator or laypersons bear the risk of injuries or crushing, which in a worst case scenario could lead to the patient's death! During normal use, you should regularly check the interco product for damage and the perfect working condition of all components.

#### Important!

##### Original spare parts

Only original parts may be used for service and maintenance work.

Checking / components	Daily	Weekly
Visual inspection of the seating system and the attached components and accessories for visible damage	x	
Check of connection between undercarriage and seating system	x	
Visual inspection of all fastening elements	x	
Visual inspection of all accessories / options	x	
Retightening of all fastening elements		x
Retightening of all screw connections of the accessories / options		x

### 11.3. Service

Test / Execution	OK	Adjustments necessary
Cleanliness	<input type="checkbox"/>	<input type="checkbox"/>
Check of all settings of the MINY® undercarriage	<input type="checkbox"/>	<input type="checkbox"/>
Check of setting of the foot rest: ▶ Is there enough free space at the back of the knee?	<input type="checkbox"/>	<input type="checkbox"/>
Check the safety distances according to DIN EN 12182:2012 using the following table	<input type="checkbox"/>	<input type="checkbox"/>

#### Safety distances according to DIN EN 12182:2012 for children

To avoid	Between moving parts	Between fixed parts
Trapping fingers	<4 mm or >25 mm	<5 mm or >12 mm
Trapping feet	<25 mm or >120 mm	<25 mm or >45 mm
Trapping the head	<60 mm or >300 mm	<60 mm or >250 mm
Trapping the genitals	<8 mm or >75 mm	<8 mm or >75 mm

#### Safety distances according to DIN EN 12182:2012 for adults

To avoid	Between moving parts	Between fixed parts
Trapping fingers	<8 mm or >25 mm	<8 mm or >25 mm
Trapping feet	<35 mm or >120 mm	<35 mm or >100 mm
Trapping the head	<120 mm or >300 mm	<120 mm or >250 mm
Trapping the genitals	<8 mm or >75 mm	<8 mm or >75 mm

#### Note!

*Expiry of the manufacturer's warranty if service is not conducted*

The product must be checked regularly – at least once a year – as part of a (payable) service.

Otherwise interco Group cannot provide any warranty.

#### 11.3.1. Seating system

Check of setting / function	OK	Damage/ distortion	Exchange/ replacement
<b>Main frame</b>			
Main frame	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Anti-tipping system	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Trapezoidal adapter</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Front wheel</b>			
Wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tire pressure	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Running performance of the wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Front wheel fork	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Mounting to frame	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Rear wheel</b>			
Wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tire pressure	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Running performance of the wheels	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Quick-release axles	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Drum brake	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Brakes</b>			
Performance test, level 1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Performance test, level 2	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Performance test, level 3	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Bowden cable	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Adjustability</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Check of setting / function		OK	Damage/ distortion	Exchange/ replacement
<b>Gas-pressure springs</b>	Damages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Screw connections	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Seat adjustment</b>	Performance test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Folding mechanism</b>	Folding mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Locking mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Back angle adjustment</b>	Performance test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Push bar</b>	Height adjustment mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Locking mechanism	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Screw connections</b>	Size and length	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Completeness	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Damages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Loctite secured	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Covering caps available	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Options</b>	Wheel protector	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Spoke protector	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	E-Fix	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Centered foot rest	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Appearance (visual inspection)</b>		<input type="checkbox"/>	<input type="checkbox"/>	
<b>Paintwork</b>		<input type="checkbox"/>	<input type="checkbox"/>	
<b>Cleanliness</b>		<input type="checkbox"/>	<input type="checkbox"/>	
<b>Rating plate</b>		<input type="checkbox"/>	<input type="checkbox"/>	

## 12. Reuse/durability

Your product can be reused under certain conditions. It is CE marked and must be reconditioned for reuse by the manufacturer to meet all essential safety and performance requirements. This includes that it is properly adapted to the needs of the new user. In particular, individual fittings can only be reinstated after professional adaptation to the new patient. Accompanying documents, as well as these operating instructions, are part of the product and must be handed over to the

new user when it is to be reused again. The durability and service life of the product are not limited and may vary. They depend on the conditions of use and the environment. In addition, daily use and the degree of disability can affect the durability and service life of the product. Regular inspection and the performance of all necessary service and maintenance work can positively influence the durability and service life of the product.

## 13. Disposal/environmental compatibility

Environmental protection is important. Therefore interco does not work with hazardous substances and tries to integrate environmental aspects into its daily work and thinking and to pass these approaches on to our suppliers if necessary. All employees are frugal in using the resources provided to us.

All packaging materials are recyclable. interco participates in the Dual System so that environment-friendly disposal of transport packaging is ensured. The product can be broken down into its individual material components at the end of its useful life. Dispose of these raw materials, depending on their type, in accordance with the applicable environmental protection regulations.

## 14. Possible faults

All serious incidents and faults occurring in connection with the product and its use must be reported immediately to the specialist dealer. Furthermore, you can inform the Compliance Officer of interco and the competent authority in your country.

Your Compliance Officer at interco:  
Ute Markwald (E-mail: [msib@interco-reha.de](mailto:msib@interco-reha.de))

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Eitorf, 2021-01-12

  
Michael Markwald  
Geschäftsführer  
Chief Executive Officer

## 16. Warranty policy

In addition to our General Terms and Conditions, we grant a two-year warranty on all products supplied by us. Unusable or defective parts will be repaired or replaced free of charge during the warranty period. The warranty shall exclude design-related modifications.

interco Group does not assume any liability for damages resulting from the combination of our products with third-party products of any kind, unless our company has given its express approval for such a product combination.

Furthermore, warranty also excludes any defects occurring as the result of normal wear, excessive stress, violent and willful damage, as well as inappropriate use. Warranty obligations shall cease to apply if original interco replacement parts are not used for repairs and if repairs and modifications are not performed by interco. The installation of additional parts or changes to the interco products by others than the manufacturer are not permitted and can be subject to legal

action if they become known. Likewise, the warranty obligations of interco Group shall cease to apply in such cases and any liability claims directed at the manufacturer shall no longer be valid.


Certain screws and nuts on your interco product may become loose when the product is used. Re-tighten these at periodic intervals or arrange to have this done by your local dealer or by interco.

### Note!

Expiry of the manufacturer's warranty  
The product must be checked regularly – at least once a year – as part of a (payable) service. Otherwise interco Group cannot provide any warranty.

interco cannot grant warranty if third party components (not manufactured by interco) are used without authorization by interco.

### Warranty extension

On receipt of this form, we will grant an extended warranty for your MINY® for a total of two years. Please fill out this warranty card carefully and completely and return it to us. Thank you very much.

Funding agency/health insurance fund: \_\_\_\_\_

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
Status as per February 2021 | Subject to technical alterations without prior notification

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Management system certified as per

**DIN EN ISO 13485:2016**

Reg.-Nr. 73 105 6351

TÜV Hessen

**DIN EN ISO 9001:2015**

Reg.-Nr. 73 100 6351

TÜV Hessen

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© interco Group GmbH, Eitorf, 2021
**interco Group GmbH**