ORTHOLIGHT®

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E-Mail*: _____

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interco GmbH
Im Auel 50
53783 Eitorf

**Operating instructions
ORTHOLIGHT® custom-made**

English

Contents

1.	Introduction	32
1.1.	Product description	32
1.2.	Product overview	32
1.3.	Installation plan	32
1.4.	Environmental conditions	32
2.	Medical requirements	33
2.1.	Field of application/Target group	33
2.2.	Indications	33
2.3.	Contraindications	33
2.4.	Intended use	33
3.	Product labeling	33
4.	Safety instructions	33
5.	Caution symbols	35
6.	Scope of delivery	35
6.1.	Basic configuration	35
6.2.	Options/Accessories	36
6.2.1.	Lower extremities	36
6.2.2.	Pelvis	36
6.2.3.	Upper body	37
6.2.4.	Upper extremities	38
6.2.5.	Head	38
7.	Settings/Operation	39
7.1.	Basic configuration	39
7.1.1.	Attachment and removal of the seat shell via adapter plate	39
7.1.2.	Positioning of the patient	39
7.1.3.	Fixing the patient	40
7.1.4.	Back angle adjustment	40
7.1.5.	Seat depth extension	40
7.2.	Options/Accessories	41
7.2.1.	Lower extremities	41
7.2.2.	Pelvis	41
7.2.3.	Upper body	42
7.2.4.	Upper extremities	43
7.2.5.	Head	43
8.	Compatibility/Combination options	44
9.	Technical data	45
10.	Transport/transport lock/climbing stairs	45
10.1.	General information on transport	45
10.2.	Transport with an AMF force application point system	45
10.3.	Transport on aircraft	45
10.4.	Climbing stairs	45
11.	Maintenance plan	46
11.1.	Cleaning	46
11.1.1.	Materials used	46
11.1.2.	Cleaning of materials	46
11.1.3.	Disinfection of the product	47
11.2.	Maintenance	47
11.3.	Service	48
11.3.1.	Seating system	48
12.	Reuse/durability	49
13.	Disposal/environmental compatibility	49
14.	Possible faults	49
15.	EU Declaration of Conformity	50
16.	Warranty policy	51

Information

Date of last update:
11/2020

Revision status:
ORTHOLIGHT-NACH-MASS_Gebrauchs-
anweisung_interco_de-en_01

Please read these operating instructions carefully! Observe all safety instructions. If you have any questions, please contact your specialist dealer or interco.

Note

As a rule, all interco products consist of a seat shell and/or an undercarriage. To simplify matters, only the interco product is mentioned in the following.

In these operating instructions we use the masculine form exclusively for the male/female users of an interco product and its male/female operators (parents, male/female caregivers, male/female therapists).

This way we wish to maintain the readability and comprehensibility of the texts. At the same time, we would like to expressly point out that this naturally includes the female form and other gender identities.

Foreword

Dear Customer,


Thank you very much for choosing a quality product from interco. This is intended to contribute to a sustainable improvement in the everyday life and quality

of life of the user. To ensure that you can always handle the product safely and use it for a long time, we have described what you should bear in mind during use in these operating instructions.

It is important that you carefully read and understand all sections that apply to your interco product. Keep the operating instructions for future reference. Please follow all the safety instructions and symbols to prevent any type of injury. If you should have any questions about your product and its use, please ask your local qualified medical supplies dealer or contact us directly. We would also be pleased to inform you about further customization and adaptation options for your interco product.

Please observe the following in particular:

- ▶ Parents, accompanying persons, nurses and family members who operate with this product should be instructed in the use of the product by authorized specialist personnel using these operating instructions.
- ▶ Your interco product was customized to the needs of the user. Subsequent changes to the product may only be made by trained specialist personnel. We recommend regular checks of the product and the product adaptation to ensure optimal long-term operation.
- ▶ Your interco product is only approved for combination with one of the products listed later on. interco will not accept any liability for combinations with medical devices and/or accessories from other manufacturers.
- ▶ Technical modifications may only be made with the written consent of interco.
- ▶ interco will not accept liability for modifications, changes to parts of the frame of the product and unapproved repairs to this product.
- ▶ interco reserves the right to technical changes to the model described in these operating instructions.

1. Introduction

1.1. Product description

ORTHOLIGHT is designed for children and adolescents who are unable to take up or maintain a seating position independently and whose movement is reduced to a minimal residual activity, for example because the interaction of the pelvis and the spinal column presents functional deficits, due to the deformation of the support apparatus or due to insufficiencies of the musculoskeletal system.

For customized construction, the seat and the back have been individually manufactured for these patients from aluminum according to the patient's measurements. ORTHOLIGHT is an aluminum seat shell con-

sisting of a seat and a back section, connected by a hinge. ORTHOLIGHT is available with removable padding in multiple sections in the standard design or in a custom-made design. ORTHOLIGHT can be fitted with an optional stepless manual back adjustment. Depending on the patient's clinical picture and individual requirements, a wide range of special accessories, such as a variety of immobilization aids, padding and supports for the head, torso, hips, arms, legs and feet is also available.

Using the trapezoidal adapter, ORTHOLIGHT can be combined with our room undercarriages, ROOMY and ROOMY NEW EDITION, and our street undercarriages, SIMPLY, SIMPLY LIGHT and MINY.

1.2. Product overview

- 1 Complete Seat shell (seat and back)
- 2 Chest harness
- 3 Lateral side supports
- 4 Arm rests
- 5 Lap belt
- 6 Trapezoidal adapter

1.3. Installation plan

Authorization for delivery

ORTHOLIGHT® may only be delivered to the end customer by authorized specialist dealers, medical device consultants or by the medical device consultants of interco GmbH.

Scope of delivery

Basic configuration includes:

- ▶ Two-piece 3mm aluminum standard base
- ▶ Connected with rear hinge
- ▶ Powder coated in various RAL colors
- ▶ With anatomical padding using the sandwich method for custom-made seat and back
- ▶ Back bolted in place at 90°
- ▶ Optional conversion to angle-adjustable back
- ▶ With high-quality fabric cover in various colors

Setup/assembly

The product is delivered fully pre-assembled to the end customer, thus no components need to be assembled on-site. If options are installed at a later date, the assembly instructions are enclosed with the option.

The medical product advisor of the specialist dealer must ensure that the instructions are carried out correctly (see below). The medical device consultant of interco GmbH does not have to be present during the installation of the product. However, he can be consulted for support in case of questions.

Instruction of the customer

All functions of ORTHOLIGHT® which are explained in the operating instructions must be demonstrated to the customer. The operating instructions remain with the customer. All adjustments to the product that directly affect the user may only be carried out by medical device consultants or in consultation with these.

1.4. Environmental conditions

The ORTHOLIGHT® medical device must always be stored in a dry place to prevent possible corrosion. Please keep to ambient temperatures of -10°C to +40°C. Sand or other dirt particles should not be allowed to penetrate the product. Use in salty water is not permissible.

2. Medical requirements

2.1. Field of application/Target group

Excerpt from the GKV Association

These seating systems can be put together from different units. The basic configuration always consists of a seat and a back in the form of a seat unit. Such modular systems can be individually adapted to the handicap due to the large number of different sizes of the individual parts without the need for custom-made products.

Children's seating systems are intended for the use of chassis for the interior and exterior. Some of these products are already equipped with a chassis as standard.

In many cases, accessories such as neck / head rests, systems for securing and supporting positioning, side supports, arm rests and therapy tables are required.

The provision of modular children's seating systems means that a favorable seating position is provided, thus creating support and positioning effects for the insured. These products aim to straighten the pelvis and stretch the spine as well as provide lateral support for the pelvis and thorax.

Due to the many different configuration options, modular children's seating systems are easier to adapt than seat shells to the size development of children without the need for a new product after a short time.

2.2. Indications

- ▶ ICD-10-G00-G99 (diseases of the nervous system)
- ▶ ICD-10-M00-M99 (diseases of the musculoskeletal system and connective tissue)
- ▶ ICD-10-P11 (Other Central Nervous System Birth Injuries)
- ▶ ICD-10-Q00-Q07 (congenital malformations of the nervous system)
- ▶ ICD-10-Q65-Q79 (Congenital malformations and deformities of the musculoskeletal system)
- ▶ ICD-10-R25-R29 (symptoms affecting the nervous system and the musculoskeletal system)


2.3. Contraindications

There are no known contraindications.

2.4. Intended use

If sitting is impaired because of severe functional and/or structural damage to the trunk of trunk- and possibly cervix muscles with low posture, a seating system may be necessary to maintain a sitting position that is as physiological as possible. This seating system consists of an outer shell, suitable upholstery and a cover. Individual seating systems are used where ready-made seating systems cannot be used. Standardized seating systems are used when a ready-made seating system can be used.

3. Product labeling

- A Type designation
- B Country of manufacture
- C Manufacturer's data
- D Serial number
- E Observe warnings/safety instructions
- F Read operating instructions prior to use
- G Size
- H Maximum weight
- I Upholstery

Please always state the serial number and size of the ORTHOLIGHT® when ordering spare parts or accessories!

4. Safety instructions

Please note that you are responsible for the safety of the user. This could be at risk if you do not follow the safety guidelines. Nevertheless, it is not possible to cover all potential circumstances and unforeseeable situations in these guidelines. Therefore common sense, caution and prudence are required when handling the product. Every person operating the product should know and understand all the instructions. Furthermore, you

should explain the operation of the product and its accessories to every person using the product. Please contact your local specialist dealer or interco directly if the instructions are unclear and a more detailed explanation is required or if you have any other questions.

Meaning of the safety instructions

Warning!

Warning against possible severe damage due to accidents or injuries.

Caution!

Warning against possible damage due to accidents or injuries.

Note!

Warning for possible technical damage.

Product-specific safety instructions

Warning!

Risk of injury from unattended parking

If the product is left unattended with user, the user or other persons may be injured.

- ▶ Please make absolutely sure that the person sitting in the product is properly secured at all times by a suitable restraint system (harness).
- ▶ Never leave the user sitting in the product unattended, not even when the user has been securely strapped in and the brakes have been set.
- ▶ Never leave the product alone with the user without expert supervision to exclude the risk of injury to the user or other persons.

Warning!

Risk of suffocation due to the packaging material

Incorrect handling of packaging material may cause a risk of suffocation.

- ▶ Keep the packaging materials out of reach of children.

Warning!

Risks when climbing stairs

When using stairs, there may be a risk of injury to the user of the product or the product may be damaged. Therefore only lift the product without the user sitting in it.

- ▶ Stairs are only permitted to be negotiated with the assistance of accompanying persons.
- ▶ If facilities such as access ramps or lifts are available for this purpose, these must be used. In the absence of such facilities, two accompanying persons can carry the product over obstacles without the user.
- ▶ Techniques for climbing stairs are explained in detail in Section 10; it is imperative that you comply with these.
- ▶ Make sure that the accompanying persons only grip the product on firmly mounted or bonded components.
- ▶ Do not lift the product by the foot rest, head rest, arm rests or similar parts.

Warning!

Risks of self-initiated modification of the settings

If settings are changed autonomously, the function of the product and the safety of the user may be affected.

- ▶ Adhere to the settings made by the specialist dealer. You may only adjust settings by yourself as described in Section 7.
- ▶ If you experience problems with settings, please contact your specialist consultant or interco directly.
- ▶ Discuss all changes to the settings properly with the specialist personnel or therapist.
- ▶ It is not permissible for any parties other than interco to make additions or modifications to the product and any cases that should become known can be subject to legal action. In addition, the warranty obligations of interco GmbH shall cease to apply in such cases and any liability claims directed at the manufacturer shall no longer be valid.

Caution!

Risk of injury due to operating errors or lacking knowledge of the operating instructions

If the instructions in the operating instructions are disregarded, application errors may occur and the safety of the user may be impaired.

Please read these operating instructions carefully first.

- ▶ Familiarize yourself with the functions of the product and its handling prior to initial use. Practice handling the product before use.
- ▶ Observe the maximum load capacity of the product (see Section 9: Technical Data).
- ▶ Do not use the foot rest when getting into or out of the product.

Caution!

Risk of tipping due to incorrect center of gravity / shift in center of gravity

If the product's center of gravity is shifted significantly, for example on slopes and inclines or if the user leans out of the product too far, the product may tip or roll over.

- ▶ Test the effects of shifts in the center of gravity on the behavior of the product, for example, on declines, inclines, and lateral slopes or while negotiating obstacles only with the reliable support of an assistant.
- ▶ Always stand the product on horizontal, flat surfaces only. If it is unavoidable to park it on an incline, put the seat in an upright position. There is a risk of the product tipping backwards on inclines in the lying position.
- ▶ When grasping objects (lying in front of, to the side or behind the device) make sure that the user does not lean out of the seat too far.
- ▶ Do not hang heavy shopping bags or similar on the push bar or the levers.

Caution!*Risks in public road traffic*

Dangerous situations can occur when driving on public roads.

- ▶ Please note that you must observe traffic regulations when traveling on public streets.
- ▶ In dark conditions wear bright clothing wherever possible in order to be more visible. We also recommend installing an active light source.

Caution!*Risk of injury to hands and fingers*

Your hands or fingers can get caught during adjustment and assembly work.

- ▶ Always be aware of your fingers during operation and all adjustment and assembly work.

Caution!*Dangers due to loosened components during use*

Screws and nuts on the product may become loose when the product is used. This can result in risks of accidents and injuries for the user and the operators.

- ▶ Check all screws and nuts regularly and retighten them.
- ▶ Alternatively, you can also instruct your specialist dealer or interco to do this.

Caution!*Dangers when making adjustments to the product and connections*

If you do not carry out adjustments to the product

properly, the stability and safety of the product may be affected.

- ▶ Ensure that the screws are tightened again after all adjustment and assembly work.
- ▶ Always make sure that the seating system is securely connected to the undercarriage.
- ▶ If you have a question regarding the connection of the seating system and the undercarriage, always contact your specialist dealer or interco directly.
- ▶ Carry out a regular visual inspection of all connections of the product. Consult your local specialist dealer in case of doubt.

Caution!*General risks due to careless handling of the product*

The product is only approved for its intended use. Improper or careless use can result in safety and injury risks for the user and other persons.


- ▶ Avoid actions, such as driving against an obstacle (stair, curb) without using the brake or down from elevated surfaces.
- ▶ Use caution when handling fire (e.g. a burning cigarette)! Upholstery and coverings could catch fire.
- ▶ Make sure to keep children from playing with the therapy chair or any parts thereof. There is a risk of injury due to trapping and crushing hands and fingers.
- ▶ Children must not play next to, underneath or in front of the seating system when the seating system is in motion and/or the seat or back unit is in use.

5. Caution symbols

Note: All screw connections, that were loosened by changing the settings, must be re-tightened !

Place of mounting: outside, at the bottom of the adduction guides

6. Scope of delivery

6.1. Basic configuration

ORTHOLIGHT® is an passive, light-weight seating system with

- ▶ two-piece 3 mm aluminum standard base
- ▶ connected with rear hinge
- ▶ powder coated in various RAL colors
- ▶ with anatomical padding using the sandwich method for custom-made seat and back
- ▶ back bolted in place at 90°
- ▶ optional conversion to angle-adjustable back
- ▶ with high-quality fabric cover in various colors

6.2. Options/Accessories

6.2.1. Lower extremities

Article	Size	Art. no.
ALULINE leg rest system with incrementally adjustable holder for leg rest system Foot rest, complete version, foot rest powder coated in black, includes corner guards and foot rest overlay for all interco undercarriages	size 0-3	46013
ALULINE leg rest system with incrementally adjustable holder for leg rest system Foot rest with individual foot supports, foot rest powder coated in black, includes corner guards and foot rest overlay for all interco undercarriages	size 0-3	46016
Leg rest, complete version, with fixed foot brackets and foot guide slots	size 1-3	46010
Leg rest, complete version, with foot bracket with continuous angle adjustment and foot guide slots	size 1-3	46012
Leg rest with individual foot supports, with foot bracket with continuous angle adjustment and foot guide slots	size 0-3	46014
Leg rest, with individual foot supports and fixed foot brackets and foot guide slots	size 0-3	46015
Holder for leg rest system as fixed 90° or 110° version	size 1-3	46020
Holder for leg rest system with continuous angle adjustment	size 1-3	46021
Calf support, complete version		46030
Black synthetic leather cover for calf support, complete version		46031
Calf support, single version		46120
Black synthetic leather cover for calf support, single version		46121
Plastic foot rest riser, black		46530
Heel stop for leg rest system, complete version		46551
Heel stop for leg rest with individual foot supports		46552
Heal stop cushion without cover		46510
Black synthetic leather cover for heel stop cushion		46511
Foot rest cushion without cover		46520
Black synthetic leather cover for foot rest cushion		46521
Foot straps with hook and loop fastener with padded instep		46410
Foot straps made of leather with padded instep		46411
Foot straps with deflection roller and padded instep with hook and loop fastener		46412
Foot straps with incremental adjustment and padded instep		46419
Foot straps with clip buckle and padded instep		46420
Foot shells made of aluminum	size 1-3	46210
Foot gaiters without forefoot restraint made of neoprene	size 1-3	46414
Foot gaiters with forefoot restraint made of neoprene	size 1-3	46416

6.2.2. Pelvis

Article	Size	Art. no.
Abduction pommel without cover, adjustable and attachable		44010
Abduction pommel, swing-away, without cover		44011
Abduction pommel, integrated, including cover		44020
Cover for abduction pommel		44025
Hip belt with deflection roller black synthetic eather or neoprene		44170u
2-point hip belt, size 1, with neoprene cushion and small plastic press buckle with one pull		44175-1
2-point hip belt with neoprene cushion and large plastic press buckle with one pull	size 2-3	44175-2
		44175-3
4-point hip belt, size 1, with neoprene cushion and small plastic press buckle with one pull		44176-1

Article	Size	Art. no.
4-point hip belt, with neoprene cushion and large plastic press buckle with one pull	size 2-3	44176-2 44176-3
2-point hip belt with neoprene cushion with two pulls	size 1-3	44177
4-point hip belt with neoprene cushion with two pulls	size 1-3	44178
Pelvic harness standard, leg-spreading or T-shaped in spacer fabrics, synthetic leather or neoprene	size 1-3	44190 44192 44198
PELVI-LOCK 3DS pelvis retraction belt with toothed belt and buckle for seat shells with flexible belt strap fastener	size XS (SB=16-24 cm), size S (SB=20-32 cm), size M (SB=30-40 cm), size L (SB=36-46 cm)	44046 44047 44048 44049
PELVI-LOCK 3DA pelvis retraction belt with toothed belt and concealed buckle for high traction forces for seat shells with rigid plastic belt strap	size XS (SB=16-24 cm), size S (SB=20-32 cm), size M (SB=30-40 cm), size L (SB=36-46 cm)	44050 44051 44052 44053
PELVI-LOCK 3DAS pelvis retraction belt with toothed belt and concealed buckle for high traction forces for seat shells with flexible belt strap fastener	size XS (SB=16-24 cm), size S (SB=20-32 cm), size M (SB=30-40 cm), size L (SB=36-46 cm)	44054 44055 44056 44057
PELVI-LOCK 3DB pelvis retraction belt with toothed belt and red ratchet buckle for seat shells with rigid plastic belt strap	size XS (SB=16-24 cm), size S (SB=20-32 cm), size M (SB=30-40 cm), size L (SB=36-46 cm)	44061 44062 44063 44064
PELVI-LOCK pelvis retraction belt with toothed belt and buckle and rigid plastic belt strap SB= 30-40cm		44060

6.2.3. Upper body

Article	Size	Art. no.
Lateral side support, right/left, fixed version		44210
Swing-away unit for lateral side support, Aluline version		44220
Lateral side support, height- and width-adjustable		44250
Lateral side support, height- and width-adjustable, integrated into the back unit (diagonally adjustable)		44251
Cover for lateral side supports		60081
Lumbar spine support, integrated		44300
Lumbar spine support, with Velcro attachment, without cover		44301
Cover for lumbar spine support, with Velcro attachment		44302
Chest belt with deflection roller	synthetic leather or neoprene	44110u
Chest belt, basic version (with Velcro attachment)	synthetic leather or neoprene size 1-3	44114
Chest belt with flap	synthetic leather or neoprene size 1-3	44115
Standard RV chest harness with push-in buckle, 25 mm, lock-stitched straps on the shoulders and zipper in the middle (includes clip buckle with fastening tab for mounting on seat shells)	for boys or girls in spacer fabrics, synthetic leather or neoprene size 1-4	47130
VARIO chest harness with twist-lock at the pelvis and push-in buckles on the upper shoulder straps (includes clip buckle with fastening tab for mounting on seat shells)	for boys or girls in spacer fabrics, synthetic leather or neoprene size 1-4	47131
REMOVABLE chest harness with removable upper straps with push-in clasp and push-in clasps at the pelvis	for boys or girls in spacer fabrics, synthetic leather or neoprene size 1-4	47132
Standard chest harness with push-in buckle, 25 mm, lock-stitched straps on the shoulders (includes clip buckle with fastening tab for mounting on seat shells)	for boys or girls in spacer fabrics, synthetic leather or neoprene size 1-4	47133

Article	Size	Art. no.
Strap guide plates (includes clip buckle compatible with straps) for mounting on the edges of the seat shell, used to guide upper shoulder strap and enable continuous height adjustment		47200

6.2.4. Upper extremities

Article	Size	Art. no.
ALULINE arm rests, height- and angle-adjustable, including PUR arm cushions (only in combination with adduction guides)		42025
Height-adjustable arm rests (15x15 mm square) including PUR arm cushions (only in combination with adduction guides)		42010
Height- and angle-adjustable arm rests (15x15 mm square) including PUR arm cushions (only in combination with adduction guides)		42020
Arm cushion, trough-shaped or flat, without cover		42030
Cover for arm cushion, trough-shaped		42031
Upper arm supports		44261
Cover for upper arm supports		44262
Therapy tray table made of acrylic, attachable, square-shaped, with upturned front and side bars and two-segment cut-out	size 1-4	42811 42812 42813 42814
Therapy tray table made of beech wood, attachable, square-shaped, with upturned front and side bars and two-segment cut-out	size 1-4	42511 42512 42513 42514
Therapy tray table made of beech wood, attachable, semicircular, with two-segment cut-out	size 1-4	42521 42522 42523 42524
Therapy tray table made of acrylic with edges around, in special measurements		42810
Swing-away unit for therapy tray table		42515
Attachment mechanism for therapy tray table for attachment with two brackets		42517
Tabletop cushion with synthetic leather cover		42610
Therapy tray table made of Neopolen, semicircular, with Velcro attachment, with additional fastening strap, with synthetic leather cover	size 1-4	42541 42542 42543 42544

6.2.5. Head

Article	Size	Art. no.
Head rest holder, angled model		43020
Head rest holder, straight model		43022
Head rest holder with continuous angle adjustment		43021
ALULINE head rest bracket, fixed model, includes 1 x ALULINE universal bracket and support block with one fastening element		43046
ALULINE head rest bracket, fixed model, includes 1 x ALULINE universal bracket and support block with two fastening elements		43048
ALULINE head rest bracket, laterally adjustable model, with one fastening element and support block, includes 1 x ALULINE universal bracket		43049
ALULINE head rest bracket with depth adjustment option		43050

Article	Size	Art. no.
ALULINE head rest bracket, fixed model, for head rests made of PU, compatible with ball and socket system with ball and socket system and support block, includes 1 x ALULINE universal bracket	with one fastening element	43056
ALULINE head rest bracket, fixed model, for head rests made of PU, compatible with ball and socket system with ball and socket system and support block, includes 1 x ALULINE universal bracket	with two fastening elements	43058
Standard head rest cushion with black synthetic leather cover	size 0-2	43010
Head rest cushion, shell-shaped, with black synthetic leather cover	size 0-2	43011
Fabric head rest cover for standard or shell-shaped head rest in addition to synthetic leather cover, easily removable using zipper		43015
Head rest cushion with lateral guides, made of PU, compatible with ball and socket system	size 1-2	43017
Head rest cushion, shell-shaped, made of PU, compatible with ball and socket system	size 1-3	43018
Head rest cushion, individually customized version without cover		43013
Head rest fabric for head rest in SA in addition to synthetic leather cover, easily removable using zipper		43015-SA

7. Settings/Operation

Every interco product is customized to the individual needs of the user and configured accordingly. To do justice to all possible product designs, the list of options and accessories as well as their descriptions is therefore somewhat more extensive.

Therefore please bear in mind that this section may also describe components that possibly do not exist on your product.

7.1. Basic configuration

7.1.1 Attachment and removal of the seat shell via adapter plate

Attaching the seat shell

The seat shell is fixed to the adapter reception plate of the underframe or wheel chair by means of the trapezoidal adapter (Fig. 1). For that purpose perform the subsequent steps and check the firm fit on the adapter plate after the seat shell has been placed. First unlock the seat shell fixation (Fig. 2).



Set the rear end of the seat shell adapter to the front of the adapter plate such that the locking is pressed downward (Fig. 3). Then push the seat shell backwards until it engages audibly. Lock the seat shell fixation again (Fig. 4). Finally, check the seat shell and the seat shell catch for firm fit.

Removing of the seat shell

First loosen the protective device. Open the seat shell catch by pressing the lever upwards. Then you can pull out the seat shell to the front.

Caution!

The removal of the seat shell from the wheel chair may only be performed without patient.

7.1.2 Positioning of the patient

In order to position and properly fix a patient in a seatshell, you should adhere to some basic rules. It is not possible to name all potential discrepancies in this short description. Should you feel insecure when using the appliance, please ask your medical device consultant or your therapist. They can best help you with problems.

When positioning a patient into the seatshell, please proceed as follows:

Caution!

Place the patient only into a seatshell that is fixed onto an undercarriage and has all necessary fixations.

- ▶ Pay attention that the seatshell is properly fixed onto the undercarriage.
- ▶ Pay attention that all breaks of the undercarriage are tightened. Read the operation manual of the undercarriage in order to do so.
- ▶ Now you can place the patient into the seatshell. Mind that all fixations (except for seating pants) are placed in such a way, that the patient does not sit on them. Also make sure that no sharp foreign bodies are placed in the seatshell.
- ▶ The pelvis should be placed as far back into the angle of the seatshell as possible. Only then an optimal support is possible.
- ▶ When applying additional fixations, always start from the pelvis. Only when those sit properly, continue with fixing the upper body and extremities.

7.1.3 Fixing the patient

Body fixations can be done with:

- ▶ Upper body pressure pads and belt
- ▶ Pelvis belt
- ▶ Seating pants

Your personal, local medical device consultant will assist you in the planning of an individual solution, taking account of your personal requirements. Please coordinate with him the correct use of the different fixations. The operation of the fixations is also described in chapter 7.2. and below.

7.1.4 Back angle adjustment

The back angle is serially adjusted at 90° (Fig. 5). Optionally the ORTHOLIGHT can be equipped with a manual continuous back angle adjustment. In order to change the back angle, loosen the screws on the sides with the clamp lever (Fig. 6). Adjust the desired back angle and tighten the screws again.

Caution!

Afterwards inspect levers to ensure that they are secure; re-tighten, as necessary!

**7.1.5 Seat depth extension**

The seat depth of the ORTHOLIGHT can be optionally extended by a total of 5 cm. For this, use the 2 elongated holes on the underside of the seat (Fig. 7). As a result, you will gain an increase of up to 5 cm in the seat depth. Tighten the hexagonal socket screw again after reaching the desired position.

Caution!

Afterwards inspect hexagonal socket screw to ensure that they are secure; re-tighten, as necessary!

Fig. 3

Fig. 4

Fig. 5


Fig. 6


Fig. 7

7.2. Options/Accessories

7.2.1. Lower extremities

Footrest

To attach the footrest, proceed as follows:

- ▶ Loosen the two wing screws on the lower right and left sides of the front (long) side of the trapezoidal adapter plate.
- ▶ Slide in the footrest to the desired depth.
- ▶ Tighten the two wing screws.

Footplate

To attach the footplate, proceed as follows:

- ▶ Loosen both wing screws on the right and left sides of the front (long) side of the trapezoidal adapter plate.
- ▶ Loosen both wing screws on the rear edge of the footplate. Ensure that the thread of the wing screws does not protrude onto the guides.
- ▶ Slide the footplate onto the ends of the footrest and fasten at the desired height with the wing screws.
- ▶ Fasten footrests with the wing screws at the trapezoidal adapter plate to the required depth.

7.2.2. Pelvis

Abduction pommel integrated

The integrated abduction pommel was affixed during manufacturing of the seat shell. Adjustments can not be made after that. Should adjustments become necessary, please contact your medical device consultant.

Abduction pommel adjustable and attachable

The adjustable abduction pommel is affixed at the bottom of the outer side of the hard seat shell. An adjustment is possible. The position of the pommel toward the inner seat can be adjusted. It should be aimed that both knees are centered toward the abduction pommel side.

To adjust the abduction pommel proceed as follows:

- ▶ Loosen the clamp by turning the clamping lever anticlockwise. The lever is placed under the seat at the mounted holder of the abduction pommel.
- ▶ Push the abduction pommel into the correct position.
- ▶ Fix the position by turning the lever clockwise until the abduction pommel can no longer be moved.

Pelvis belt with pressure lock

Requirement: The patient is sitting in the seatshell. The pelvic brace is placed on the pelvis with the insertion slide.

- ▶ Push the insertion slide into the belt lock until it hearably locks while pushing down the iliac crest with one hand.
- ▶ Pull the belt back until the pelvic brace sits in the correct position.

- ▶ Check the correct seating position and the even pressure on the iliac crest. The two large upholsteries of the pelvic brace should sit central and with even pressure pressure on the iliac crest.
- ▶ To open the belt lock push the red Button in the middle of the belt lock. Pay attention to hold the patient at the same time.

Caution!

The adjustment of the belt is usually done once. This adjustment remains the same even after opening the belt lock. You should check the correct adjustment however after every time.

Pelvis belt with reverse rolls

Requirement: The patient is sitting in the seatshell. Place the pelvic belt with the reverse rolls over the pelvis.

- ▶ Pull the remaining end of the belt through the reverse rolls.
- ▶ Pull until the pelvic belt sits in its final position. You should check this by testing the correct seat and the even pressure on the iliac crest.
- ▶ Fix this position by using the affixed velcro fastener.
- ▶ To open the pelvic belt loosen the velcro fastener.

Caution!

The adjustment of the belt is necessary every time. Always check the correct fit after applying.

Restraining pants

Requirement: Before positioning the patient in the seatshell, make sure that the restraining pants lay on the the seat with no pleats. Sit the patient into the seatshell. Proceed as follows:

- ▶ The restraining pants have two ends at the knees, one belt for each leg with a insertion lock.
- ▶ Place the respective belt around the upper leg and push the insertion slide into the insertion adapter, which is on the side of the seatshell.
- ▶ Repeat this procedure on the other side.
- ▶ Check the correct fit of the restraining pants. The belts should support the pelvis without being too tight.

Caution!

Avoid pressure through pleats (danger of pressure marks!)

7.2.3. Upper body

Height and width adjustable side supports

The side supports are serially fixed. Optionally they can be adjusted continuously with a slot hole in the back and support piece (Fig. 8). In order to do so, loosen the screws with a 4 allen key. Move the side supports in the slot holes according to the requirements. Afterwards tighten the screws in the desired position.

Caution!

Afterwards inspect hexagonal socket screw to ensure that they are secure; re-tighten, as necessary!

Retractable side supports

The retractable side supports were already integrated into the side supports during production of the seat-shell. No changes can be undertaken to the measurements of the upholstery afterwards. Should an adaptation become necessary, please contact your medical device consultant.

In order to retract the side supports, please proceed as follows:

- ▶ Grap the front of the upholstery with one hand.
- ▶ With the other hand press the red part of the release lever on the outside of the retraction hinge. Retract the side support to the side while holding the release lever.
- ▶ In order to bring the side supports back into their original position, simply move them back until they audibly lock.

Side supports with no retraction function

The side supports with no retraction function were already integrated into the side pieces of the ORTHOLIGHT back during manufacturing. No changes can be undertaken afterwards. Should an adaptation become necessary, please contact your medical device consultant.

Pressure pads for breast and shoulder

The breast and shoulder pressure pads are used to fix the upper body of the patient. The use is described in the next two paragraphs.

Caution!

Adjustments should always be done from the hard back of the shell through the 4 belt clamps.

Fixation of the respective patient

Requirement: The person is already sitting in the seat-shell, the pressure pads are laying backwards on the shell rim. Please proceed as follows:

- ▶ Bring the shoulder pads to the front. The head is now between the belts.
- ▶ On the bottom part of the breast and shoulder pressure pads you find two belt clams to fix the belts.
- ▶ Open the black clamp.
- ▶ Push the belt in between the movable flap and the fixed part of the respective side that comes out between the seat and the back.
- ▶ Pull the belt as far as possible, hold it and fold the flap of the clamp back into the fixed part.


Fig. 8

- ▶ The belt is now fixed.
- ▶ To open the lock, lift the flap of the clamp and pull the belt.
- ▶ Repeat the process for the second belt.
- ▶ To loosen the fixation of the breast and shoulder pressure pads, proceed in opposite direction.

Adjustment of the breast and shoulder pressure pads

Caution!

Adjustments should always be done from the hard back of the shell through the 4 belt clamps.

Requirement: The person is already sitting in the seat-shell, the pressure pads are fixed.

To adjust the breast and shoulder pressure pads, proceed as follows:

- ▶ At first adjust the lower parts of the pressure pads.
- ▶ Open the lower belt clamps at the hard back of the seatshell by simply pulling the free end of the belt, the clamp will open automatically.
- ▶ Bring the belt into the desired position, hold it and fold the flap of the clamp back into the fixed part.
- ▶ Repeat the adjustment on the other side.
- ▶ Repeat the procedure with the upper ends of the pressure pads. They are adjusted at the hard back of the seatshell. The procedure is identical to the adjustment of the lower ends.

Breast belt with flip-over flap

The breast belt with a flip-over flap was mounted in the correct position when adjusting the seatshell. Small adjustments are possible by using the velcro fastener.

To open the belt proceed as follows:

- ▶ Pull the front flap off from the bottom to the top.
- ▶ Pull the overlapping ends of the breast belt apart.
- ▶ You can lay the ends of the breast belt over both sides to lift the patient out of the seat.
- ▶ To fix the patient in the seatshell proceed in opposite order.

Caution!

Pay attention that after applying, nothing is constricted and breathing is not impaired.

Breast belt with reverse rolls

The breast belt with reverse rolls was mounted in the correct position when adjusting the seatshell. Small adjustments are possible by using the velcro fastener.

To open the belt proceed as follows:

- ▶ Pull the overlapping ends of the breast belt apart.
- ▶ Pull the loose end of the belt through the metal flap.
- ▶ You can lay the ends of the breast belt over both sides to lift the patient out of the seat.
- ▶ To fix the patient in the seatshell proceed in opposite order.

7.2.4. Upper extremities

Arm rests

ORTHOLIGHT® can be equipped with different arm rests (Fig. 9) Fix the arm rest by loosening the clamping lever of the arm rest receptacle on the right and left side of the seat shell, insert the arm rest holder and tighten the lever again. Adjustments are always done by changing the depth of the insertion in the square pipe.

Hollow-form armrest

The hollow-form armrests are used as an alternative to the standard armrests. The only difference is the hollow-form that is used to position the patients' arms in. To change the height and angle, please read chapter 7.2.4. (arm rests).

Therapy table

ORTHOLIGHT® can be equipped with different therapy tables if there are armrests on the seat shell. The brackets of the therapy table are attached under the armrests. The therapy table is inserted into the brackets and can be fixed by tightening the wing screw. By loosening the wing screw, the therapy table can be removed again or, depending on the brackets, swing away.

7.2.5. Head

Headrest systems ALULINE

ORTHOLIGHT can be equipped with different head rest systems from the series ALULINE (straight or shell shaped). Both are adjustable in height as well as depth. Fix the head rest by loosening the clamping lever of the head rest receptacle on the upper part of the back rest, insert the head rest holder and tighten the lever again. Adjustments are always done by changing the depth of the insertion in the square pipe.


To adjust the depth proceed as follows:

- ▶ Turn the fixing lever anticlockwise until the square pipe can be moved.
- ▶ Change the square pipe to the desired position.
- ▶ Turn the black fixing lever clockwise to fix the square pipe.

Head rest cropped

Components for head positioning are affixed to the hard shell with universal holders. Adjustments are always done by changing the depth of the insertion in the square pipe. Please proceed as follows:

- ▶ Turn the black fixing lever anticlockwise until the square pipe can be moved.
- ▶ Change the square pipe to the desired position.
- ▶ Turn the black fixing lever clockwise to fix the

square pipe.

- ▶ Repeat this procedure on the other side.

You can adjust the cropped head rest on two directions, height and depth.

To adjust the height proceed as follows:

- ▶ The universal holder for the head rest is on the back of the seat shell.
- To adjust, follow the instructions above.

To adjust the depth proceed as follows:

- ▶ Loosen the clamp on the upper cross tube by turning the fixing lever anticlockwise.
- ▶ Move the head rest pipe in the upper cross tube until you reach the desired position.
- ▶ Fix the position by turning the fixing lever clockwise until the head rest pipe can not be moved anymore.

Head rest straight

You can adjust the straight head rest on two directions, height and depth. The procedure is the same as described as for the head rest cropped (see above).

Head rest holder with conus adjustment

You can adjust any head rest in two ways with the head rest holder with conus adjustment, in height and in depth.

For height adjustment proceed as follows:

- ▶ The universal holder for the head rest is affixed to the back of the seat shell. The procedure is the same as described as for the head rest cropped (see above).

For depth adjustment proceed as follows:

- ▶ Open the clamp of the head rest holder by turning the fixing lever anticlockwise.

The head rest holder can now be moved.

- ▶ Hold the head rest in the desired position and turn the fixing lever clockwise to fix it.
- ▶ The position of the head rest remains the same now.
- ▶ It is now necessary to readjust the head rest upholstery.
- ▶ Turn the head rest in its fixation spot until the head of the patient rests in its respective upholstery.

Head rest upholstery

Three differently shaped head rests that are affixed with the head rest holders described above are available for the seat shell ORTHOLIGHT. The head rest upholsteries are:

- ▶ Standard
- ▶ Shell shaped

The velcro head rest upholstery are not affixed with the head rest holder, but they are simply attached with the velcro fixation.

Standard, shell shaped, half roll head rest upholstery
 These head rest upholsteries are cropped with the head rest holder, straight or conus, by a screwing mechanism. The basic adjustment of the head rest is done at the holder.

8. Compatibility/Combination options

Undercarriage	Size 1	Size 2	Size 3	individual max. SB 32 cm
SIMPLY SB 45 · 22"/24" wheels	-	-	✓	-
SIMPLY SB 45 · 12"/16" wheels	-	-	✓	-
SIMPLY SB 40 · 22"/24" wheels	✓	✓	✓	✓
SIMPLY SB 40 · 12"/16" wheels	✓	✓	✓	✓
SIMPLY SB 35 · 22"/24" wheels	✓	✓	✓	✓
			without wheel cover	without wheel cover
SIMPLY SB 35 · 12"/16" wheels	✓	✓	✓	✓
SIMPLY LIGHT SB 45 · 22"/24" wheels	-	-	✓	-
SIMPLY LIGHT SB 45 · 12"/16" wheels	-	-	✓	-
SIMPLY LIGHT SB 40 · 22"/24" wheels	✓	✓	✓	✓
SIMPLY LIGHT SB 40 · 12"/16" wheels	✓	✓	✓	✓
SIMPLY LIGHT SB 35 · 22"/24" wheels	✓	✓	✓	✓
			without wheel cover	without wheel cover
SIMPLY LIGHT SB 35 · 12"/16" wheels	✓	✓	✓	✓
ROOMY NEW EDITION · 12"/16" wheels outdoor	✓	✓	✓	✓
ROOMY NEW EDITION · 100/125mm wheels indoor	✓	✓	✓	✓
MINY seat shell wheelchair SB40 · 22"/24" wheels	✓	✓	✓	✓
MINY seat shell wheelchair SB40 · 12"/16" wheels	✓	✓	✓	✓
MINY seat shell wheelchair SB35 · 22"/24" wheels	✓	✓	✓	-
			without wheel cover	
MINY seat shell wheelchair SB35 · 12"/16" wheels	✓	✓	✓	✓
MINY seat shell wheelchair SB30 · 22"/24" wheels	✓	✓	-	-
MINY seat shell wheelchair SB30 · 12"/16" wheels	✓	✓	-	-

Undercarriage	size 4	size 5	size 6	individual from SB 32 cm to max. 41 cm
SIMPLY SB 45 · 22"/24" wheels	✓	✓	✓	✓
			without wheel cover max. 80 kg	
SIMPLY SB 45 · 12"/16" wheels	✓	✓	✓	✓
			max. 80 kg	
SIMPLY SB 40 · 22"/24" wheels	✓	✓	-	-
		without wheel cover		
SIMPLY SB 40 · 12"/16" wheels	✓	✓	-	-
SIMPLY SB 35 · 22"/24" wheels	-	-	-	-
SIMPLY SB 35 · 12"/16" wheels	-	-	-	-
SIMPLY LIGHT SB 45 · 22"/24" wheels	✓	✓	✓	✓
			without wheel cover max. 80 kg	
SIMPLY LIGHT SB 45 · 12"/16" wheels	✓	✓	✓	✓
			max. 80 kg	
SIMPLY LIGHT SB 40 · 22"/24" wheels	✓	✓	-	-
		without wheel cover		
SIMPLY LIGHT SB 40 · 12"/16" wheels	✓	✓	-	-
SIMPLY LIGHT SB 35 · 22"/24" wheels	-	-	-	-
SIMPLY LIGHT SB 35 · 12"/16" wheels	-	-	-	-

Undercarriage	size 4	size 5	size 6	individual from SB 32 cm to max. 41 cm
ROOMY NEW EDITION · 12"/16" wheels outdoor	✓	✓ max. 70 kg	✓ max. 70 kg	✓
ROOMY NEW EDITION · 100/125mm wheels indoor	✓	✓ max. 70 kg	✓ max. 70 kg	✓
MINY seat shell wheelchair SB40 · 22"/24" wheels	✓	✓ without wheel cover	–	–
MINY seat shell wheelchair SB40 · 12"/16" wheels	✓	✓ max. 60 kg	–	–
MINY seat shell wheelchair SB35 · 22"/24" wheels	–	–	–	–
MINY seat shell wheelchair SB35 · 12"/16" wheels	–	–	–	–
MINY seat shell wheelchair SB30 · 22"/24" wheels	–	–	–	–
MINY seat shell wheelchair SB30 · 12"/16" wheels	–	–	–	–

9. Technical data

Caution!

The maximum load (workload) for ORTHOLIGHT® is 20 kg for size 1, 25 kg for size 2, 30 kg for size 3, 35 kg for size 4, 40 kg for size 5 and 45 kg for size 6.

	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Weight without patient (complete with upholstery, cover, arm rest, head rest and abduction pommel)	approx. 7 kg	approx. 7,5 kg	approx. 8 kg	approx. 7 kg	approx. 7,5 kg	approx. 8 kg
Total outer width without armrests	34 cm	37 cm	40 cm	34 cm	37 cm	40 cm
Total outer width with armrests	44 cm	47 cm	50 cm	44 cm	47 cm	50 cm
Back-angle adjustment	30°	30°	30°	30°	30°	30°
Seat width max.	26 cm	29 cm	32 cm	26 cm	29 cm	32 cm
Seat depth	20-25 cm	25-30 cm	30-35 cm	20-25 cm	25-30 cm	30-35 cm
Back height	38 cm	42 cm	46 cm	38 cm	42 cm	46 cm

10. Transport/transport lock/climbing stairs

10.1. General information on transport

The ORTHOLIGHT® product in combination with a TÜV approved interco undercarriage is permissible as a seat in a motor vehicle.

Warning!

Dangers due to incorrect fastening of the transport safety device

If you attach the safety device for transport in a vehicle to the wrong places on the product, it cannot be transported safely.

- ▶ Never secure the transport safety device to the frame or other unsuitable basic components!

recommend that you contact the airline before traveling by air to clarify the transport options in each individual case.

Caution!

Transport on aircraft

The user must never be in the device when it is transported in the cargo hold of an aircraft!

10.2. Transport with an AMF force application point system

This option is only permissible for ORTHOLIGHT® with an interco undercarriage.

10.3. Transport on aircraft

When traveling by air, the product may only be transported without the user in the cargo hold of the aircraft.

There are no further instructions or regulations for transport in an aircraft. We therefore strongly recom-

10.4. Climbing stairs

Stairs or high curbs can only be negotiated with the assistance of two accompanying persons. Assistants are only permitted to grip the front frame tube and the back brace of the frame, without lifting or carrying the product when doing so.

When ascending, the assistant behind the device pulls the device up the step. The assistant in the front stabilizes the position by pressing the device equally on both sides into the step.

When descending, the assistant at the front brakes by pressing the device evenly into the steps. The assistant behind holds the back brace, and stabilizes and

maintains the device in the correct position. To prevent damage, slowly roll the product down the steps. Then properly readjust the anti-tipping system (if installed).

Caution!

Risks when climbing stairs

If you hold or lift the product by its attachments when climbing stairs, these may tear off and the user may be injured.

- ▶ Do not hold the product by the push bar when negotiating stairs.
- ▶ Never lift the product by the foot rest, the wheels or the arm rests or other attached parts.
- ▶ Do not lift the product with the user sitting in it.

11. Maintenance plan

11.1. Cleaning

Your interco product is made of high-quality and durable materials which are easy to clean. Please observe the following instructions for long durability of the product:

- ▶ Always keep the product dry, at temperatures between -10 and +40 °C. This prevents corrosion.
- ▶ Sand or other dirt particles should not be allowed to penetrate the product. This could damage the mechanics or - if applicable - the wheels.
- ▶ Never use the product in salt water.
- ▶ The cushions can be removed for cleaning and care.

11.1.1. Materials used

interco uses the following materials in the manufacturing of medical devices and accessories:

- ▶ Metals (aluminum, steel, stainless steel)
- ▶ Plastics (OPM, PA6, PTEG)
- ▶ Beech wood (for therapy tray tables)
- ▶ Foams (frame foam, polyurethane foam, hard foam, azote foam, neoprene)
- ▶ Upholstery fabrics (mesh material, Trevira CS fabrics, artificial leather, Dartex)
- ▶ Powder coatings
- ▶ DIN and standard parts

All materials used are classified as harmless to health.

Powder-coated aluminum

interco undercarriages and some add-on parts are made of aluminum and feature a resistant plastic coating. Powder varnishes are used for coating. Powder varnishes are solvent-free coating materials, composed of different resin systems, pigments as well as special additives.

Fabrics

interco uses different fabrics for covers and fixation aids. All materials used are flame-retardant.

Synthetic leather

interco uses easy-care designer synthetic leather for Dynaline restraints. This synthetic leather is long-lasting, hard-wearing, very elastic and tear-resistant. It provides a good grip and offers a comfortable seating experience. In addition, it is resistant to urine, blood and disinfectants, free of azo compounds, CFC, formaldehyde, as well as free of PCTs, PCPs and PCBs.

Spacer fabrics/mesh material

The spacer fabrics used have two outer textile layers, which are connected by a number of spacer threads. This 3D structure creates an insulating layer of air between the textile layers, which balances temperature fluctuations effectively. Furthermore, the spacer fabrics are very breathable, pressure-elastic, light and tear-resistant.

Supra-tex

The SUPRA-tex upholstery fabric consists of highly tear-resistant polyamide. It is impregnated to be water-repellent, PU-coated on the back and flame-retardant.

Neoprene

The neoprene used is a foam based on chloroprene rubber (CR). It insulates very well, is waterproof and very well tolerated by the skin. In addition, the neoprene is tear and abrasion resistant, very elastic, resistant to fungi and bacteria and eco-friendly.

Trevira CS (fabrics)

Among others, interco uses Trevira CS fabrics for seat covers. These are easy to clean, dry quickly and keep their shape. In addition, Trevira CS can be disinfected under hospital conditions.

Dartex incontinence cover fabric

interco uses the incontinence-resistant upholstery fabric Dartex for incontinence-resistant seat shells. The high-quality fabric meets all quality and safety requirements and is easy to wash off, sufficiently lightfast and flame-retardant. It does not cause allergies, is free of formaldehyde and azo dyes.

11.1.2. Cleaning of materials

Powder-coated aluminum

Clean undercarriages and attachments made of powder-coated aluminum with mild soapy water and a soft sponge. Do not use aggressive or abrasive agents. These can damage the surface of the coating.

Fabrics

Clean the fabric covers regularly before stains set in. You can vacuum or carefully brush the textiles. Spilled liquids should be wiped up immediately with an absorbent cloth. You can also put a mild detergent solution on a cloth, dab the stain with it and then rub off gently with clear water.

Note!

Bleach or whiteners should never be applied to the fabrics / covers. Never dry the covers in the dryer! Only iron at a low temperature.

If you wish to wash the cover of the seating system more often, a second cover for changing is recommended. Ask your local specialized dealers on this topic.

Synthetic leather

Clean synthetic leather covers with a cloth and commercially available mild soap solution.

Spacer fabrics/mesh material

Wash removable covers made of spacer fabrics/mesh material at a maximum of 30 °C in the washing machine with the delicates program.

Supra-tex

interco uses Supra-tex for the arm rests. Clean these regularly with a cloth and commercially available mild soapy water.

Neoprene

Wash removable fixation aids made of neoprene at a maximum of 30 °C in the washing machine. Use a mild liquid detergent or a commercially available neoprene detergent. Use a gentle wash cycle and do not spin-dry the neoprene covers. Then hang up the covers wet to dry.

Trevira CS (fabrics)

Removable covers made of Trevira CS can be dry cleaned. Washing under easy-care conditions at 60°C is possible. Do not dry and do not iron!

Dartex

Test washed by the manufacturer at temperatures between 40 and 95 °C.

Accessory components made of metal, plastic, wood

Please clean accessory components with a conventional mild household soap solution.

11.1.3. Disinfection of the product

Please use a water-based agent for disinfection of the product, e.g. Terralin, Quartamon Med or Sagrotan Original concentrate. Please observe the manufacturer's instructions for use when doing so. When disinfecting the fabrics, also follow the instructions

of the fabric manufacturer regarding the use of brand name disinfectants.

Rinse off all detergents and disinfectants thoroughly and leave the cleaned object to dry completely before use. If this is neglected, this can lead to chemical reactions making the material brittle or lead to skin irritations of the patient.

Important!

No disinfectants for synthetic leather/skai
Synthetic leather covers/skai covers may not be cleaned with disinfecting agents, as the material gets porous over time and cracks.

**Note!**

Wipe disinfection
If wipe disinfection is necessary, use a disinfectant that is suitable for surface disinfection according to the DGH disinfectant list.

**11.2. Maintenance**

Your interco product is checked by trained and qualified employees for proper functioning and freedom from defects before being handed over to the user. To ensure safe use of the product over its entire service life (according to the manufacturer from at least 2 years to 6 years maximum) and to ensure a perfect technical condition, it must be maintained and inspected regularly according to the following maintenance (every 6 months) and service (every 12 months) plans.

Note!

Changing settings
The changing of settings (other than those described in Section 7) is only permitted for qualified and trained personnel!

Incorrect settings by the user, operator or laypersons bear the risk of injuries or crushing, which in a worst case scenario could lead to the patient's death!

During normal use, you should regularly check the interco product for damage and the perfect working condition of all components.

Important!

Original spare parts
Only original parts may be used for service and maintenance work.

Checking / components	Daily	Weekly
Visual inspection of the seating system and the attached components and accessories for visible damage	x	
Check of connection between undercarriage and seating system	x	
Visual inspection of all fastening elements	x	
Visual inspection of all accessories / options	x	
Retightening of all fastening elements		x
Retightening of all screw connections of the accessories / options		x

11.3. Service

Test / Execution	OK	Adjustments necessary
Cleanliness	<input type="checkbox"/>	<input type="checkbox"/>
Check of all settings of the ORTHOLIGHT® seating system	<input type="checkbox"/>	<input type="checkbox"/>
Check of setting of the head rest:	<input type="checkbox"/>	<input type="checkbox"/>
▶ Does the head rest lie against the patient's occiput (standard version)?		
▶ Does the bottom edge of the head rest lie against the base of the head (shell shape)?		
▶ Do the side elements exert pressure on the patient's head?		
▶ Does the patient have sufficient freedom of movement?		
Check of setting of the foot rest:	<input type="checkbox"/>	<input type="checkbox"/>
▶ Is there enough free space at the back of the knee?		
Check of setting of the arm rests:	<input type="checkbox"/>	<input type="checkbox"/>
▶ Are the shoulders in a neutral / ergonomic position when the elbows are in the resting position?		
Check of setting of the lateral side supports:	<input type="checkbox"/>	<input type="checkbox"/>
▶ Is there enough free space under the armpits?		
▶ Is the upper body guided sufficiently?		
Check the safety distances according to DIN EN 12182:2012 using the following table	<input type="checkbox"/>	<input type="checkbox"/>

Safety distances according to DIN EN 12182:2012 for children

To avoid	Between moving parts	Between fixed parts
Trapping fingers	<4 mm or >25 mm	<5 mm or >12 mm
Trapping feet	<25 mm or >120 mm	<25 mm or >45 mm
Trapping the head	<60 mm or >300 mm	<60 mm or >250 mm
Trapping the genitals	<8 mm or >75 mm	<8 mm or >75 mm

Safety distances according to DIN EN 12182:2012 for adults

To avoid	Between moving parts	Between fixed parts
Trapping fingers	<8 mm or >25 mm	<8 mm or >25 mm
Trapping feet	<35 mm or >120 mm	<35 mm or >100 mm
Trapping the head	<120 mm or >300 mm	<120 mm or >250 mm
Trapping the genitals	<8 mm or >75 mm	<8 mm or >75 mm

Note!

Expiry of the manufacturer's warranty if service is not conducted

The product must be checked regularly – at least once a year – as part of a (payable) service.

Otherwise interco GmbH cannot provide any warranty.

11.3.1. Seating system

Check of setting / function	OK	Damage/ distortion	Exchange/ replacement
Performance test General	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Screw connections Screws tightened	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Rust spots	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Damages	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Covering caps	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Completeness	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Appearance Paintwork	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Cleanliness	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Fabric seams	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Rating plate	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Main frame General	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Seat adjustment Performance test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Check of setting / function		OK	Damage/ distortion	Exchange/ replacement
Back angle adjustment	Performance test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Options	Arm rests: + mounting	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Head rest + mounting	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Abduction pommel	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Therapy tray table + mounting	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Harness + mounting	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Lateral side supports	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Upper arm supports	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Calf support + mounting	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Heel stop + mounting	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

12. Reuse / durability

Your product can be reused under certain conditions. It is CE marked and must be reconditioned for reuse by the manufacturer to meet all essential safety and performance requirements. This includes that it is properly adapted to the needs of the new user. In particular, individual fittings can only be reinstated after professional adaptation to the new patient. Accompanying documents, as well as these operating instructions, are part of the product and must be handed over to the

new user when it is to be reused again. The durability and service life of the product are not limited and may vary. They depend on the conditions of use and the environment. In addition, daily use and the degree of disability can affect the durability and service life of the product. Regular inspection and the performance of all necessary service and maintenance work can positively influence the durability and service life of the product.

13. Disposal / environmental compatibility

Environmental protection is important. Therefore interco does not work with hazardous substances and tries to integrate environmental aspects into its daily work and thinking and to pass these approaches on to our suppliers if necessary. All employees are frugal in using the resources provided to us.

All packaging materials are recyclable. interco participates in the Dual System so that environment-friendly disposal of transport packaging is ensured. The product can be broken down into its individual material components at the end of its useful life. Dispose of these raw materials, depending on their type, in accordance with the applicable environmental protection regulations.

14. Possible faults

All serious incidents and faults occurring in connection with the product and its use must be reported immediately to the specialist dealer. Furthermore, you can inform the Compliance Officer of interco and the competent authority in your country.

Your Compliance Officer at interco:
Ute Markwald (E-mail: msib@interco-reha.de)

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REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 on medical devices, annex IV, annex IX, chapter 1

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16. Warranty policy

In addition to our General Terms and Conditions, we grant a two-year warranty on all products supplied by us. Unusable or defective parts will be repaired or replaced free of charge during the warranty period. The warranty shall exclude design-related modifications.

interco GmbH does not assume any liability for damages resulting from the combination of our products with third-party products of any kind, unless our company has given its express approval for such a product combination.

Furthermore, warranty also excludes any defects occurring as the result of normal wear, excessive stress, violent and willful damage, as well as inappropriate use. Warranty obligations shall cease to apply if original interco replacement parts are not used for repairs and if repairs and modifications are not performed by interco. The installation of additional parts or changes to the interco products by others than the manufacturer are not permitted and can be subject to legal

action if they become known. Likewise, the warranty obligations of interco GmbH shall cease to apply in such cases and any liability claims directed at the manufacturer shall no longer be valid.


Certain screws and nuts on your interco product may become loose when the product is used. Re-tighten these at periodic intervals or arrange to have this done by your local dealer or by interco.

Note!

Expiry of the manufacturer's warranty
The product must be checked regularly – at least once a year – as part of a (payable) service. Otherwise interco GmbH cannot provide any warranty.

interco cannot grant warranty if third party components (not manufactured by interco) are used without authorization by interco.

Warranty extension

On receipt of this form, we will grant an extended warranty for your ORTHOLIGHT® for a total of two years. Please fill out this warranty card carefully and completely and return it to us. Thank you very much.

Funding agency/health insurance fund: _____

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Status as per December 2020 | Subject to technical alterations without prior notification

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© interco GmbH, Eitorf, 2020