

AKTIVLINE® BASIS

HMV no.: in the application process

Order Cost estimate

Name of the orderer

Company

Order/Patient name

Your order no.

Date

Signature/Company stamp

AKTIVLINE dynamic care system, standardized version complete with the MINY aluminum undercarriage

- ▶ seat and back unit made of aluminum with physiological pivot points in the pelvic and knee regions for the performance of dynamic movements
- ▶ anatomical standard upholstery for seat and back units
- ▶ height- and width-adjustable lateral side supports
- ▶ seat width, seat depth, back height all adjustable by approx. 5 cm
- ▶ AFS system (Abduction pommel - leg and Foot rest - System), complete version, swing-away
- ▶ foot rest, flip-up
- ▶ with 12"/16"/22"/24" wheels, including drum brakes with quick release axles in back
- ▶ 200x50 front wheels
- ▶ parking brake handles with three catches, height- and angle-adjustable push bar
- ▶ 2 in 1 anti-tipping system (support roller and anti-tip all in one), retractable beneath the undercarriage
- ▶ seating system removable from the undercarriage, undercarriage foldable for transport
- ▶ adjustable overall seat height and adjustable seat depth on undercarriage



Technical Data

	Size 1	Size 2	Size 3	Size 1	Size 2	Size 3
Wheels	12"/16"	12"/16"	12"/16"	22"/24"	22"/24"	22"/24"
Frame width, undercarriage	35 cm	40 cm	40 cm	35 cm	40 cm	40 cm
Total weight	approx. 27 kg	approx. 28 kg	approx. 30 kg	approx. 27 kg	approx. 28 kg	approx. 30 kg
Seat width	18 - 24 cm	24 - 30 cm	30 - 36 cm	18 - 24 cm	24 - 30 cm	30 - 36 cm
Seat depth	20 - 25 cm	27 - 32 cm	34 - 39 cm	20 - 25 cm	27 - 32 cm	34 - 39 cm
Back height	33 - 38 cm	38 - 43 cm	43 - 48 cm	33 - 38 cm	38 - 43 cm	43 - 48 cm
Exterior width, seat unit without arm rest	28 - 34 cm	34 - 40 cm	40-46 cm	28 - 34 cm	34 - 40 cm	40 - 46 cm
Height, push bar	90 - 105 cm	90 - 105 cm	90 - 105 cm	90 - 105 cm	90 - 105 cm	90 - 105 cm
Overall width, (wheel camber)	55 cm (+ 5 cm)	60 cm (+ 5 cm)	60 cm (+ 5 cm)	55 cm (+ 5 cm)	60 cm (+ 5 cm)	60 cm (+ 5 cm)
Max. length	80 cm	80 cm	80 cm	80 cm with 22" wheels / 85 cm with 24" wheels		
Max. load	60 kg	60 kg	60 kg	60 kg	60 kg	60 kg
Frame colour	Brilliant Silver 08/G050					
Fabric Trevira CS (black) / Spacerfabrics	Orange	Green	Blue	Orange	Green	Blue

Art. no.	39800-1	39800-2	39800-3	39800-1	39800-2	39800-3
Suggested retail price*	9,868.90 €	9,868.90 €	9,980.40 €	9,868.90 €	9,868.90 €	9,980.40 €
Unit	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh

MINY Equipment

MINY aluminum undercarriage:

- ▶ aluminum base frame, slightly abducted version
- ▶ tilt-in-space mechanism using gas-pressure springs by approx. 35°
- ▶ different wheel sizes with pneumatic tires or solid tires
- ▶ steering wheels available in 200x50 or 7"

Rear wheels

- 12" tires 16" tires 22" tires 24" tires
 pneumatic tires solid tires

Front wheels

- 200 × 50 mm wheels 7" wheels
 pneumatic tires solid tires

Art. no.	Article name	Article description	Options	UNIT	Price *	Quantity
----------	--------------	---------------------	---------	------	---------	----------

Head



43046	ALULINE head rest bracket, fixed model,		<input type="checkbox"/> with one fastening element	ST	268.70 €	<input type="text"/>
43048	with support block, includes 1 x ALULINE universal bracket		<input type="checkbox"/> with two fastening elements		316.70 €	<input type="text"/>
43049	ALULINE head rest bracket, laterally adjustable model,			ST	370.20 €	<input type="text"/>
	with one fastening element and support block, includes 1 x ALULINE universal bracket					
43010	Standard head rest cushion		<input type="checkbox"/> size 0	ST	141.60 €	<input type="text"/>
	with black synthetic leather cover		<input type="checkbox"/> size 1			
			<input type="checkbox"/> size 2			
43011	Head rest cushion, shell-shaped,		<input type="checkbox"/> size 0	ST	141.60 €	<input type="text"/>
	with black synthetic leather cover		<input type="checkbox"/> size 1			
			<input type="checkbox"/> size 2			
43015	Fabric head rest cover for standard or shell-shaped head rest			ST	73.60 €	<input type="text"/>
	in addition to synthetic leather cover, easily removable using zipper					



Upper body

44301	Lumbar spine support, with Velcro attachment, without cover			ST	95.90 €	<input type="text"/>
	under the cushion					
44110u	Chest belt with deflection roller		<input type="checkbox"/> synthetic leather	ST	89.20 €	<input type="text"/>
			<input type="checkbox"/> neoprene			
44114	Chest belt, basic version (with Velcro attachment)	<input type="checkbox"/> made of synthetic leather	<input type="checkbox"/> size 1	ST	137.20 €	<input type="text"/>
		<input type="checkbox"/> made of neoprene	<input type="checkbox"/> size 2			
			<input type="checkbox"/> size 3			
44115	Chest belt with flap	<input type="checkbox"/> made of synthetic leather	<input type="checkbox"/> size 1	ST	155.00 €	<input type="text"/>
		<input type="checkbox"/> made of neoprene	<input type="checkbox"/> size 2			
			<input type="checkbox"/> size 3			

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh

Art. no.	Article name Article description	Options	UNIT	Price *	Quantity
47130	Standard RV chest harness with push-in buckle, 25 mm, lock-stitched straps on the shoulders and zipper in the middle (includes clip buckle with fastening tab for mounting on seat shells)	<input type="checkbox"/> for boys <input type="checkbox"/> for girls <hr/> <input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3 <input type="checkbox"/> size 4	<input type="checkbox"/> spacer fabrics <input type="checkbox"/> synthetic leather <input type="checkbox"/> neoprene	ST 214.10 € 189.60 € 281.00 €	<input type="text"/>
47131	VARIO chest harness with twist-lock at the pelvis and push-in buckles on the upper shoulder straps (includes clip buckle with fastening tab for mounting on seat shells)	<input type="checkbox"/> for boys <input type="checkbox"/> for girls <hr/> <input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3 <input type="checkbox"/> size 4	<input type="checkbox"/> spacer fabrics <input type="checkbox"/> synthetic leather <input type="checkbox"/> neoprene	ST 220.80 € 196.30 € 281.00 €	<input type="text"/>
47132	REMOVABLE chest harness with removable upper straps with push-in clasp and push-in clasps at the pelvis	<input type="checkbox"/> for boys <input type="checkbox"/> for girls <hr/> <input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3 <input type="checkbox"/> size 4	<input type="checkbox"/> spacer fabrics <input type="checkbox"/> synthetic leather <input type="checkbox"/> neoprene	ST 214.10 € 196.30 € 281.00 €	<input type="text"/>
47133	Standard chest harness with push-in buckle, 25 mm, lock-stitched straps on the shoulders (includes clip buckle with fastening tab for mounting on seat shells)	<input type="checkbox"/> for boys <input type="checkbox"/> for girls <hr/> <input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3 <input type="checkbox"/> size 4	<input type="checkbox"/> spacer fabrics <input type="checkbox"/> synthetic leather <input type="checkbox"/> neoprene	ST 209.60 € 161.70 € 267.60 €	<input type="text"/>



Upper extremities

42025	ALULINE arm rests, height- and angle-adjustable, including PUR arm cushions (only in combination with adduction guides)		PA	310.00 €	<input type="text"/>
42811		<input type="checkbox"/> size 1		322.30 €	
42812	Therapy tray table made of acrylic, attachable,	<input type="checkbox"/> size 2	ST	334.50 €	<input type="text"/>
42813	square-shaped, with upturned front and side bars and two-segment cut-out	<input type="checkbox"/> size 3		345.70 €	
42814		<input type="checkbox"/> size 4		352.40 €	
42511		<input type="checkbox"/> size 1		262.10 €	
42512	Therapy tray table made of beech wood, attachable,	<input type="checkbox"/> size 2	ST	262.10 €	<input type="text"/>
42513	square-shaped, with upturned front and side bars and two-segment cut-out	<input type="checkbox"/> size 3		274.30 €	
42514		<input type="checkbox"/> size 4		274.30 €	
42515	Swing-away unit for therapy tray table		ST	142.70 €	<input type="text"/>
42517	Attachment mechanism for therapy tray table for attachment with two brackets		ST	174.00 €	<input type="text"/>
42610	Tabletop cushion		ST	107.00 €	<input type="text"/>
42611	Synthetic leather cover for tabletop cushion		ST	75.80 €	<input type="text"/>

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh

Art. no.	Article name Article description Options	UNIT	Price *	Quantity
----------	--	------	---------	----------



Pelvis

44170u	Hip belt with deflection roller	<input type="checkbox"/> black synthetic leather <input type="checkbox"/> neoprene	ST	90.30 €	<input type="text"/>
44175-1	2-point hip belt, size 1, with neoprene cushion and small plastic press buckle with one pull		ST	179.50 €	<input type="text"/>
44175-2 44175-3	2-point hip belt with neoprene cushion and large plastic press buckle with one pull	<input type="checkbox"/> size 2 <input type="checkbox"/> size 3	ST	179.50 €	<input type="text"/>
44176-1	4-point hip belt, size 1, with neoprene cushion and small plastic press buckle with one pull		ST	190.70 €	<input type="text"/>
44176-2 44176-3	4-point hip belt with neoprene cushion and large plastic press buckle with one pull	<input type="checkbox"/> size 2 <input type="checkbox"/> size 3	ST	190.70 €	<input type="text"/>
44177	2-point hip belt with neoprene cushion with two pulls	<input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3	ST	184.00 €	<input type="text"/>
44178	4-point hip belt with neoprene cushion with two pulls	<input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3	ST	195.10 €	<input type="text"/>
44046 44047 44048 44049	PELVI-LOCK 3DS pelvis retraction belt with toothed belt and buckle for seat shells with flexible belt strap fastener	<input type="checkbox"/> size XS (SB=16-24 cm) <input type="checkbox"/> size S (SB=20-32 cm) <input type="checkbox"/> size M (SB=30-40 cm) <input type="checkbox"/> size L (SB=36-46 cm)	ST	304.40 €	<input type="text"/>
44050 44051 44052 44053	PELVI-LOCK 3DA pelvis retraction belt with toothed belt and concealed buckle for high traction forces for seat shells with rigid plastic belt strap	<input type="checkbox"/> size XS (SB=16-24 cm) <input type="checkbox"/> size S (SB=20-32 cm) <input type="checkbox"/> size M (SB=30-40 cm) <input type="checkbox"/> size L (SB=36-46 cm)	ST	299.34 €	<input type="text"/>
44054 44055 44056 44057	PELVI-LOCK 3DAS pelvis retraction belt with toothed belt and concealed buckle for high traction forces for seat shells with flexible belt strap fastener	<input type="checkbox"/> size XS (SB=16-24 cm) <input type="checkbox"/> size S (SB=20-32 cm) <input type="checkbox"/> size M (SB=30-40 cm) <input type="checkbox"/> size L (SB=36-46 cm)	ST	299.34 €	<input type="text"/>

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh

Art. no.	Article name Article description Options	UNIT	Price *	Quantity
44061	PELVI-LOCK 3DB			
44062	pelvis retraction belt with toothed belt and			
44063	red ratchet buckle			
44064	for seat shells with rigid plastic belt strap			
	<input type="checkbox"/> size XS (SB=16-24 cm) <input type="checkbox"/> size S (SB=20-32 cm) <input type="checkbox"/> size M (SB=30-40 cm) <input type="checkbox"/> size L (SB=36-46 cm)	ST	304.40 €	<input type="text"/>



Lower extremities

46520	Foot rest cushion without cover	Dimensions in cm <input type="text"/>	ST	75.80 €	<input type="text"/>
46521	Black synthetic leather cover for foot rest cushion		ST	53.50 €	<input type="text"/>
46410	Foot straps with hook and loop fastener with padded instep		PA	35.70 €	<input type="text"/>
46412	Foot straps with deflection roller and padded instep with hook and loop fastener		PA	89.20 €	<input type="text"/>
46419	Foot straps with incremental adjustment and padded instep		PA	113.70 €	<input type="text"/>
46420	Foot straps with clip buckle and padded instep		PA	113.70 €	<input type="text"/>
46210	Foot shells made of aluminum	<input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3	PA	102.60 €	<input type="text"/>
46414	Foot gaiters without forefoot restraint made of neoprene	<input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3	PA	165.00 €	<input type="text"/>
46416	Foot gaiters with forefoot restraint made of neoprene	<input type="checkbox"/> size 1 <input type="checkbox"/> size 2 <input type="checkbox"/> size 3	PA	197.40 €	<input type="text"/>

Total

plus statutory VAT

Final total

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh

Notes

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh