

MAGICLIGHT[®] REDESIGN

HMV no.: 26.11.04.0011 · Customs tariff no.: 90211010

Order Cost estimate

Name of the orderer	Company
Order/Patient name	Your order no.
Date	Signature/Company stamp

MAGICLIGHT[®] REDESIGN – multifunctional, dynamic seating system

- ▶ with integrated, adjustable abduction and windswept
- ▶ dynamic backrest continuous adjustable and lockable
- ▶ continuous adjustable back extension
- ▶ continuous adjustable seat depth
- ▶ continuous adjustable seat width
- ▶ anatomical molded seat
- ▶ adjustable side supports (adjustment in height and width)
- ▶ adapter system for interco undercarriages
- ▶ fabric cover for all cushions in different colours and with lettering Magiclight



Technical Data

	Size 1	Size 2	Size 3
Max. load	20 kg	25 kg	35 kg
Max. total width (without armrests)	28 - 35 cm	35 - 42 cm	42 - 49 cm
Abduction angle	-10° to + 20°	-10° to + 20°	-10° to + 20°
Total height including adapter	41 - 48 cm	48 - 55 cm	55 - 62 cm
Backrest adjustment	-10° to + 30°	-10° to + 30°	-10° to + 30°
Seat width	20 - 27 cm	27 - 34 cm	34 - 41 cm
Seat depth	20 - 27 cm	27 - 34 cm	34 - 41 cm
Back height	30 - 37 cm	37 - 44 cm	44 - 51 cm
Height of side supports from seat upholstery	23 - 28 cm	30 - 35 cm	37 - 42 cm
Width between side supports	12 - 19 cm	19 - 26 cm	26 - 33 cm

Art. no. 37101_NE 37102_NE 37103_NE

Suggested retail price * **2,174.50 €** **2,196.80 €** **2,219.10 €**

Unit

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh

Fabric cover

Outside the fabric cover of the MAGICLIGHT is standard in UNI Black, for inside please choose from the following options.

Stoffbezug Innen

<input type="checkbox"/> DOTS Red HW10907-12	<input type="checkbox"/> UNI Red HW10907-22	<input type="checkbox"/> UNI Grey HW10907-27
<input type="checkbox"/> DOTS Blue HW10907-13	<input type="checkbox"/> UNI Blue HW10907-24	<input type="checkbox"/> UNI Black HW10907-21
<input type="checkbox"/> DOTS Green HW10907-15	<input type="checkbox"/> UNI Green HW10907-26	<input type="checkbox"/> Mesh Material Black HW10907-98S
<input type="checkbox"/> DOTS Black HW10907-18	<input type="checkbox"/> UNI Purple HW10907-23	<input type="checkbox"/> Oxidon Noir HW10907-94

Art. no.	Article name Article description Options	UNIT	Price *	Quantity
----------	--	------	---------	----------

Covers



37205-1_NE 37205-2_NE 37205-3_NE	Replaceable cover, complete for MAGICLIGHT NE Seat and backcover in different types and side cushions in different colours	<input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 <input type="checkbox"/> Size 3	SA	286.60 €	<input type="text"/>
47011	Water-repellent coating for complete seat shell		ST	304.40 €	<input type="text"/>
47012	Water-repellent coating for seat unit		ST	152.80 €	<input type="text"/>

Head



37250_NE	ALULINE head rest bracket for MAGICLIGHT NE Size 1-3 fixed model, including 1 x Aluline universal bracket with one fastening element		ST	268.70 €	<input type="text"/>
37251_NE	ALULINE head rest bracket for MAGICLIGHT NE Size 1-3 with depth adjustment option		ST	328.00 €	<input type="text"/>
37260_NE 37262_NE 37264_NE	Head rest cushion, standard with removable fabric cover like MAGICLIGHT NE side elements or with black synthetic leather	<input type="checkbox"/> Size 0 <input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2	ST	141.60 €	<input type="text"/>
37261_NE 37263_NE 37265_NE	Head rest cushion, shell-shaped with removable fabric cover like MAGICLIGHT NE side elements or with black synthetic leather	<input type="checkbox"/> Size 0 <input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2	ST	141.60 €	<input type="text"/>



Upper body

37211_NE 37212_NE 37213_NE	Chest belt in neoprene for boys or girls including mounting material and strap guide plates	<input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 <input type="checkbox"/> Size 3	ST	320.00 €	<input type="text"/>
----------------------------------	---	---	----	----------	----------------------



Upper extremities

37245-1_NE 37245-2_NE 37245-3_NE	Upper arm support including mounting material as refitting and cover	<input type="checkbox"/> Size 1 for MAGICLIGHT NE Size 1 <input type="checkbox"/> Size 2 for MAGICLIGHT NE Size 2 <input type="checkbox"/> Size 3 for MAGICLIGHT NE Size 3	PA	267.60 €	<input type="text"/>
--	---	--	----	----------	----------------------

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh

Art. no.	Article name Article description Options	UNIT	Price *	Quantity
37270-1_NE 37270-2_NE 37270-3_NE	ALULINE arm rests height- and angle-adjustable including universal holder and cushion	<input type="checkbox"/> Size 1 for MAGICLIGHT NE Size 1 <input type="checkbox"/> Size 2 for MAGICLIGHT NE Size 2 <input type="checkbox"/> Size 3 for MAGICLIGHT NE Size 3	PA 312.10 €	<input type="text"/>
37380-1_NE 37380-2_NE 37380-3_NE	Therapy tray table made of acrylic with adjustable holder	<input type="checkbox"/> Size 1 for MAGICLIGHT NE Size 1 <input type="checkbox"/> Size 2 for MAGICLIGHT NE Size 2 <input type="checkbox"/> Size 3 for MAGICLIGHT NE Size 3	ST 338.60 €	<input type="text"/>
37390-1_NE 37390-2_NE 37390-3_NE	Tabletop cushions including cover	<input type="checkbox"/> Size 1 for MAGICLIGHT NE Size 1 <input type="checkbox"/> Size 2 for MAGICLIGHT NE Size 2 <input type="checkbox"/> Size 3 for MAGICLIGHT NE Size 3	ST 179.90 €	<input type="text"/>
37240-1_NE 37240-2_NE 37240-3_NE	Lateral side support swing-away including upholstery and cover	<input type="checkbox"/> Size 1 for MAGICLIGHT NE Size 1 <input type="checkbox"/> Size 2 for MAGICLIGHT NE Size 2 <input type="checkbox"/> Size 3 for MAGICLIGHT NE Size 3	PA 312.10 €	<input type="text"/>
37241-1_NE 37241-2_NE 37241-3_NE	Lateral side support pluggable including upholstery and cover	<input type="checkbox"/> Size 1 for MAGICLIGHT NE Size 1 <input type="checkbox"/> Size 2 for MAGICLIGHT NE Size 2 <input type="checkbox"/> Size 3 for MAGICLIGHT NE Size 3	PA 312.10 €	<input type="text"/>



Pelvis

62540	Abductionpommel removable and plug-in version including cover	ST	238.10 €	<input type="text"/>
62541	Cover for abductionpommel removable or plug-in		49.70 €	<input type="text"/>
37220-1_NE 37220-2_NE 37220-3_NE	2-point hip belt with two pulls with neoprene cushion and 25 mm straps	<input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 <input type="checkbox"/> Size 3	ST 190.40 €	<input type="text"/>
37221-1_NE 37221-2_NE 37221-3_NE	4-point hip belt with two pulls with neoprene cushion and 25 mm straps	<input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 <input type="checkbox"/> Size 3	ST 190.40 €	<input type="text"/>
37222-1_NE 37222-2_NE 37222-3_NE	2-point hip belt with two pulls with neoprene cushion and and magnet buckle	<input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 <input type="checkbox"/> Size 3	ST 211.60 €	<input type="text"/>
37223-1_NE 37223-2_NE 37223-3_NE	4-point hip belt with two pulls with neoprene cushion and and magnet buckle	<input type="checkbox"/> Size 1 <input type="checkbox"/> Size 2 <input type="checkbox"/> Size 3	ST 211.60 €	<input type="text"/>



Lower extremities

37310-1_NE 37310-2_NE 37310-3_NE	Footrest system with separated footrest for mounting on the adduction guides, foot plates in black, including corner guards and overlay	<input type="checkbox"/> Size 1 for MAGICLIGHT NE Size 1 <input type="checkbox"/> Size 2 for MAGICLIGHT NE Size 2 <input type="checkbox"/> Size 3 for MAGICLIGHT NE Size 3	ST 396.80 €	<input type="text"/>
37150	Foot rest cushion including black synthetic leather cover		PA 132.30 €	<input type="text"/>
46552	Heel stop for leg rest with individual foot supports including upholstery and cover		PA 116.40 €	<input type="text"/>
37120	Calf support, single version including black synthetic leather cover		PA 208.40 €	<input type="text"/>

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh

Art. no.	Article name Article description Options	UNIT	Price *	Quantity	
46414-1	Foot gaiters made of neoprene	<input type="checkbox"/> Size 1	PA	165.00 €	<input type="text"/>
46414-2		<input type="checkbox"/> Size 2			
46414-3		<input type="checkbox"/> Size 3			
46416-1	Foot gaiters with forefoot restraint made of neoprene	<input type="checkbox"/> Size 1	PA	197.40 €	<input type="text"/>
46416-2		<input type="checkbox"/> Size 2			
46416-3		<input type="checkbox"/> Size 3			

Sonderausstattung



39050	KUSCHLY snuggle bag custom-made for different seat shells	ST	556.50 €	<input type="text"/>
39050i	KUSCHLY inner bag made of real lamb skin for KUSCHLY snuggle bag	ST	390.30 €	<input type="text"/>

Total	<input type="text"/>
plus statutory VAT	<input type="text"/>
Final total	<input type="text"/>

Notes

By fax to: +49 (0) 2243/8807-29

By email to: info@interco.gmbh